
INFOREX - 4.2 (2001)

Melting-Experiment Database for Igneous Petrologists and Geochemists



VERNADSKY INSTITUTE
RUSSIAN ACADEMY OF SCIENCES
Moscow, Russia

CONTENTS

Preface	3
Chapter 1. Introduction	5
1.1. Significance of the INFOREX database	5
1.2. Hardware requirements	6
1.3. Availability	6
Chapter 2. General Description	7
2.1. Database files	8
2.1.1. Structure and linkage of basic data files	8
2.1.2. Data management system files	11
2.1.3 Auxiliary files	14
2.2. Directory structure	14
Chapter 3. Setting Up INFOREX	15
Chapter 4. Description of User Interface	16
4.1. System Configuration (SC)	17
4.2. General System Information (GSI)	21
4.3. View / Edit / Input Data Files (<i>V/E/I</i>)	23
4.3.1. The use of View option	23
4.3.2. The use of Edit option	25
4.3.3. The use of Input option	26
4.4. Selecting Run Conditions (SRC)	29
4.5. Setting Phase Assemblages (SPA)	32
4.6. Exchange INFOREX Data (<i>EID</i>)	39
4.7. Petrological Calculations (<i>PC</i>)	41
4.7.1. Mineral-Melt Geothermometers	42
4.7.2. Fe ³⁺ /Fe ²⁺ ratios	48
4.7.3. Water Solubility	51
4.7.4. Special Tools	55
4.7.5. Ternary Diagrams	58
Chapter 5. Troubleshooting	60
5.1. Memory problems	60
5.2. Correction of mismatch problems	62
References	63
Appendix 1	65
On-line help (HELP.TXT file)	65
INFOREX abbreviations (NOTES.TXT file)	71
Appendix 2	73
List of INFOREX-4.0 (1997) experimental studies	73


Preface

<u>Welcome to the INFOREX melting-experiment database!</u> This is a sophisticated database manager for igneous petrologists and geochemists designed as an aid in searching the published experimental database, and extracting specific subsets of that data. In addition, the program allows the user to derive numerical mineral-melt expressions based on any data set, as well as to verify phase equilibria models in the field of igneous petrology.

Its major function is as a computerized melting-experiment reference manual that presents data in a fixed format, with key word search functions and flags as an integral part of the program. At present, the system accesses information on 290 experimental studies carried out from 1962 to 2000. The database includes 10,591 individual runs with 14,933 coexisting phase compositions for 44 minerals plus melt. Most of the data are based on "dry" experiments, but 4,638 include volatile components, such as H₂O and CO₂. INFOREX is hot-key driven, allowing users to sort easily and quickly through the data using any set of experimental conditions or rock type parameters.

<u>The main options include:</u> (1) configuration of INFOREX; (2) extraction of general information on the current state of the system; (3) updating and editing the database files; (4) selection of run conditions of interest (setting rock types plus a range of pressures, temperatures, oxygen fugacities, run duration as well as types of containers and resultant phase assemblages); (5) setting phase assemblages for the selected set of experiments; (6) export/import operations; and (7) petrological calculations, including calibrations of mineral-melt geothermometers, computations of Fe³⁺/Fe²⁺ ratios in experimental glasses, estimates of water solubility for H₂O-saturated experiments, and the calculations of projection parameters to plot glass compositions onto ternary diagrams.

The INFOREX database is a powerful tool for the use in the development of petrogenetic models by providing easy access and a means for manipulating mineral-melt equilibria data. As far as the fashion of its compositional files is same as used in COMAGMAT (Ariskin et al., 1993; Ariskin, 1999) and METEOMOD (Ariskin et al., 1997b) programs, the INFOREX database might be especially useful for petrologists who have been working with these phase equilibria models.






Figure 1. Statistics of melting-experiment information included into the main archive of the INFOREX-4.0 (1997) petrological database

Now the system includes as much as 10,600 experiments from 290 studies

Chapter 1. Introduction

1.1. Significance of the INFOREX Database

The experimental investigation of phase equilibria in natural and synthetic silicate systems is one of the basic fields of igneous petrology. Almost all petrologic models used for interpretation of diverse natural volcanic or intrusion rocks use experimental phase equilibria information. The development of microprobe analysis and advanced experimental techniques has produced literally thousands of major and trace element analyses of experimental runs.

These data have been presented in the literature and are preserved in personal or institutional databases around the world. This information, obtained at considerable expense, is a scientific resource of great importance for the further development of igneous petrology, especially in the field of geothermometry and computer simulation of magma differentiation processes (*Nielsen and Dungan, 1983; Ariskin et al., 1987-1999; Nielsen, 1990; Weaver and Langmuir, 1990; Longhi, 1991; Ghiorso, 1995; Yang et al., 1996*).

At present, to use the available experimental data one has to laboriously navigate through the literature and evaluate the quality of the information case by case regardless of the type of application. In fact the ability to test two-phase geothermometers and geobarometers available in the literature, or to develop them independently may require different criteria on data selection. There is a need, therefore, for a quick and easy means to access and manipulate the experimental data, select data by petrochemical criteria (igneous rock types), or run conditions (e.g., temperature, pressure, or oxygen fugacity).

Our approach to this problem has been to develop an experimental database, and to make that database available to all interested investigators (*Ariskin et al.*, 1992-1997a, *Meshalkin et al.*, 1996; *Meshalkin and Ariskin*, 1996).

1.2. Hardware Requirements

<u>INFOREX - 4.2 USER</u> MANUAL

The main program manager, called *INFOMAN*, and other codes of the INFOREX User Interface (*UI*) were written in Microsoft C 2.5 for an IBM PC or compatible, whereas all of the *Petrological Calculation* programs were written in FORTRAN-77 with Microsoft Fortran 5 compiler. Thus, the INFOREX-4.2 programs operate under DOS 3.0 or higher on any IBM PC compatible computer with at least 1 MB RAM and a hard disk with 7 MB available disk space (data base files of 3 MB, programs of 2 MB, and a space for writing operation of 2 MB). The INFOREX system can be successfully used also under MS-DOS emulated from WINDOWS 95.

The system is distributed on a 3.5" diskette and provided with an installation program; all of the the INFOREX datafiles are compressed. A hard disk drive with fast access is more critical to efficient operation. To speed up the searches and preserve your hard disk from too many writing operations one can use a virtual disk with ≈7 MB space. The software uses SVGA monitors.

1.3. Availability

Input of the melting-experiment information as well as the progamming work has been conducted by a group of scientists working at the Vernadsky Institute (Moscow, Russia): Alexei A.Ariskin, Sergei S. Meshalkin, Galina S. Barmina, Renat R.Al'meev, and Georgy S. Nikolaev were engaged in the INFOREX development process. A complete version of the INFOREX database (including the User's Manual) is distributed for a fee and can be requested from the leader of the group using the address:

Dr. Alexei A. Ariskin
Vernadsky Institute, Kosygin Str. 19
Moscow 117975, RUSSIA

Fax: 7-095-938-2054, E-mail: ariskin@geokhi.ru, ariskin@rambler.ru

Chapter 2. General Description

6

In this chapter we describe the directory structure and organization of the INFOREX database files. We developed the database so as to be useful for the full range of experimental studies in igneous petrology. Towards that end, the database was designed to have a structured, rather than flat-file architecture (Figs 2, 5).




Figure 2. Structure of the INFOREX database files. The letters j, k, u, v, w denote different run numbers

2.1. Database Files

In general, the INFOREX database files can be divided into three main categories:

basic data files, data management system files, and auxiliary files.

<u>The basic data files</u> contain information on the reference source (BIBL.TXT file), experimental conditions and phase assemblages (CONDIT.EXP), phase compositions (LIQ.EXP, PLAG.EXP, SPIN.EXP, etc.), volatile components (VOLAT.EXP) and others. Examples of these files are given in Fig. 3. <u>The data management system files</u> are designed to store information on topics users can access using the help functions, plus the general parameters of the INFOREX environment and configuration: these include working subsets of the names of phase composition files, rock types, and container materials. <u>The auxiliary files</u> facilitate the manipulation and data swapping among the basic data files using the User Interface.

<u>Another group of INFOREX files</u> includes those formed after searching through the database or created during petrological calculations. These output INFOREX files will be discussed in a detail in *chapters 4.5-4.7*.

2.1.1. Structure and linkage of basic data files

All INFOREX database files are ASCII alphanumeric rows with fixed length. The basic unit upon which the database is built is an individual experimental study. Each study has a reference record in the BIBL.TXT file and a set of master records in the CONDIT.EXP file (Fig. 3). Records in BIBL.EXP consist of up to four strings and include a working number and reference description (author, source, year, etc.). Note that the user can utilize any sequence of the bibliographic information if necessary. The same working number is given in the first column of the master records in CONDIT.EXP. This is used to link the master records to the BIBL.TXT file records. The second column in the CONDIT.EXP file represents the run numbers of the experimental study considered.

1. Reference source records in BIBL.TXT file:

N	Reference: Athors. Title //Source, Year, V. ?, N ?, P. ??-??.
	Stolper E. Experimental petrology of eucritic meteorites // Geochim. Cosmochim. Acta, 1977, V. 41, P. 587-611.

2. Master records of the CONDIT.EXP file:

N	n	Bulk Syst	P,kb	Temp	lgf02	Dur,hr	Con		Pha	ase a	asser	nbla	ges	
32	1	MEL - 0 V	30.0	1200	MH*	0.3	AGP	LQ	OL1					
34	1	LHR - 1 N	60.0	1850	ND	0.2	GRA	LQ1	OL2		İ	•	İ	
104	1	BAS + 0 N	0.00	1200	-10.98	48.0	COR	LQ1	PL	OL	SF		İ	
111	1	AND + 1 F	2.0	800	-14.18	102.0	AUC	LQ1	PL1	SF1	CP	HB	İ	
124	1	KIM - 1 V	10.0	1400	ND	0.2	PTC	LQ	OL2		1		1	
147	1	ALB - 1 N	0.0	1200	-8.00	36.0	AGP	LQ1	RE		1		1	
149	8	BAS + 0 N	0.0	1202	-8.01	17.5	PTL	LQ1	OL	RE	I	ĺ	I	
211	7	CHAR- 2 V	6.9	925	ND	168.0	AUC	LQ1	PL2	AF	OP2	IL	MT	AP
226	1	RHY - 1 F	0.8	850	NNO+1.0	144.0	PTC	LQ1	WS		1		1	
231	43			1075							I	l	I	
232	2	DIO - 1 V	6.9	950	ND	96.0	AUC	LQ	PL2	AF	OP2	GR2	SP	*

3. Records of a regular phase composition file (LIQ.EXP):

SiO2	Ti02 Al203	FeO	MnO	MgO	CaO	Na20	К20	P205	Cr203	Н2О	N	n
46.30	0.18 4.80	8.90	0.17	32.90	4.40	0.46	0.00	0.00	0.00	0.00	34	1
56.81	2.47 18.50	.88	.00	7.61	7.97	3.10	1.16	0.00	0.00	0.00	104	1
64.80	0.20 15.11	1.90	0.11	0.15	1.59	4.52	5.15	0.00	0.00	6.02	111	1
47.08	3.02 15.90	11.71	0.16	6.14	10.71	2.86	0.86	0.48	0.00	0.00	147	1
47.99	1.96 15.09	11.07	0.28	5.84	8.28	3.07	4.75	0.00	0.00	0.00	149	8
72.92	0.31 14.56	1.82	0.08	0.37	0.89	2.58	6.17	0.33	0.00	0.00	211	7
77.70	0.07 13.00	0.38	0.04	0.05	0.52	4.08	4.19	0.00	0.00	1.88	226	1
77.30	0.13 12.42	1.03	0.03	0.05	0.35	4.12	4.53	0.00	0.00	6.80	231	1

4. Examples of records of irregular phase compositions files: \boldsymbol{a} - VOLAT.EXP

CHAI	RGE	BULK (+	-sil),	wt.%	BULK (+	sil)	,mol%	BULK(-	-sil),	mol%	X(vap	oor),m	01%
N	n	Н2О	C02	Н2	Н2О	C02	Н2	Н2О	CO2	Н2	Н2О	C02	Н2
32	1	10.0	0.0	0.0	ND	0.0	0.0	100.0	0.0	0.0	100.0	0.0	0.0
124 226	1 1	6.2 ND	4.8 ND	0.0	ND 3.0	ND ND	0.0	76.0 ND			ND 32.0		0.0

b - MNLS.EXP

N	n	An (PL)	Fo(OL)	En (AU) W	o(AU)	En (PG)W	o (PG)	En (OP) W	o(OP)	(GR)_	_(GR)
34	1	.0	94.4	.0	.0	.0	.0	.0	.0	.0	.0
211	7	32.0	0.0	0.0	0.0	0.0	0.0	53.0	0.6	0.0	0.0
232	2	36.0	0.0	0.0	0.0	0.0	0.0	38.0	1.0	72.0	24.0

c - SULF.EXP

N	n	lgfS2	Melt S,wt.%	Imr S	niscible Fe	Sulfic Mn	de Liqui Ni	id, wt.9 Co	add A	REMARKS
104 111	1 1	-1.00 -3.32	0.060 0.020	00.000 37.810	00.000 61.070	00.000	00.000	00.000	00.000	

d - REDOX.EXP

N	n	WT% FeO(II)					Mol.ratio Cr3+/Cr2+	
147 149 231	1 8 43	9.55 8.80 0.88	2.40 2.51 2.61	.257	.00	.00	.000	

Figure 3. Examples of the INFOREX-4.2 (2001) basic data files Usually this run number corresponds to the sequence of run conditions and phase assemblages given in tables in the source, but this sequence has been changed

in cases where retention of the original numbers would have caused difficulty in the manipulation of the data.

<u>The CONDIT.EXP master file records</u> (Figs 2 and 3.2) include essential information about each run, such as: (1) the type of igneous rock or synthetic system, (2) a reference number for the starting material (which is linked to a composition in the START.EXP file, as well as (3) flag of the presence ($\mathbf{V/F}$) or absence (\mathbf{N}) of volatile components (where \mathbf{F} denotes water/fluid saturated experiments and \mathbf{V} - only undersaturated with volatiles runs). These switches refer the user to the VOLAT.EXP file (Fig. 3.4.a). Other important information fields are total pressure, temperature, oxygen fugacity, run duration, and up to 7 codenames for the observed phases in the run products. In addition, the master record can include flags for special "water solubility" experiments (\mathbf{WS}), where a value of water content in melt is available, or for estimated quantities of Fe³⁺ and Fe²⁺ species in the silicate glasses (\mathbf{RE}).

A number of other records are contained in additional subfiles linked to the main master record files. The individual subrecords contain information on the experimental product compositions written in the form of regular and irregular formatted strings. Names of the "normally" formatted files are the abbreviations of mineral names with the ".EXP" extensions (see Appendix 1). These files (Fig. 3.3) contain records formatted as 12 floating point number fields plus two integer number fields, which correspond to the basic experimental study and the current run number in BIBL.TXT and CONDIT.EXP files. The first twelve fields are to store contents of major oxides (in weight percents) in the sequence: SiO₂, TiO₂, Al₂O₃, FeO_{tot}, MnO, MgO, CaO, Na₂O, K₂O, P₂O₅, Cr₂O₃, and H₂O (same format and sequence as in the START.EXP file). If the chemical composition of a phase is available, the flag in the main master record is set at "1" instead of as a space as the third character in the codename of a specific particular phase present in the assemblage.

<u>The irregularly formatted files</u> were developed to store information that is beyond the scope of the normally formatted files. One of the most important files is VOLAT.EXP. Subrecords of this file include weight and molar percents of H₂O, CO₂,

and H₂ referenced to the bulk system (silicates+volatiles) or bulk volatile composition (minus silicates) as well as the volatile component contents in a vapor phase for fluid-saturated experiments (if present) - Fig. 3.4.a.

The second irregular file shown in Fig. 3.4.*b* is MNLS.EXP. It was designed to accumulate phase compositions recalculated by authors to end-member component contents, such as forsterite for olivine, anorthite for plagioclase, etc. Our preference in creating this database was to store data that was as complete as possible. Therefore, presentation of mineral compositions as endmembers only results in significant loss of information. Nevertheless, we have created these files so that we will not lose even this rudimentary experimental information, if other data on phase compositions are absent. The MNLS.EXP file is linked to the master record by the symbol "2" just after the indicator of a particular mineral phase (see Fig. 3.2).

Two more irregular files are SULF.EXP and REDOX.EXP. The first one (Fig. 3.4.*c*) can contain weight percent S, Fe, Mn, Ni, and Co - the components of majority of sulfide phases (if present). The second file (Fig. 3.4.*d*) contains the estimated quantities of Fe³⁺ and Fe²⁺ species in the silicate glasses. They are linked to the main master record in CONDIT.EXP by using the codenames *SF* and *RE*, respectively.

2.1.2. Data management system files

The core of the INFOREX database is composed of three kinds of system files:

- 1) files containing information on the INFOREX environment, including the general directory structure, run products configuration and searching control system,
 - 2) HELP.TXT file responsible for the on-line INFOREX help, and
- 3) a file called NOTES.TXT that clears the main INFOREX abbreviations, being simultaneously linked with the INFOREX *System Configuration* option (see *chapter 4.1*).

<u>The main environment file</u> **CONFIG.INF** is located in the INFOREX root directory and consists of one record describing the current operating session ...\DATA subdirectory, e.g.:

D:\INFOREX\DATA.97\

In the D:\INFOREX\DATA.97\ subdirectory we have placed the second system file **CONFIG.CNF**, which contains four strings linked to four INFOREX directories developed to work with the individual system options:

D:\INFOREX\
D:\INFOREX\SAVE\
D:\INFOREX\DATA.97\
D:\INFOREX\EXCHANGE\

CONFIG.SYS file (STANDARD) Peridotite PER PXT Pyroxenite LHR Lherzolite HAR Harzburgite Kimberlite KIM PTC Picrite KOM Komatiite MAR MiAmRu-rock AMP Amphibolite ECL Eclogite Melilitite ME.L. NEP Nephelinite NOR Norite BON Boninite BAS Basalt ALB AlkalineBAS HAB High-Al BAS ANB Andesit.BAS AND Andesite DAC Dacite RHY Rhyolite SYE Syenite TRA Trachyte DIO Diorite GRD Granodior. GRN Granite GNE Gneiss SED Sedimentary EUC Eucrite Howardite HOW DGN Diogenite SNC Shergottite ${\tt WIA}$ WIncAllende CHN Chondrite ANT ANT series Green glass GGL Pyrolite PYR SYN Synthetic MNT Aug Minette TON Tonalite CHA Charnockite

These working directory pointers provide the user with the flexibility necessary to be able to simultaneously work with different subsets of experimental data, for example with "volatile", "alkali basalt" or "high pressure" systems.

The third group of data management system files were developed to store current working sets of igneous rock names (CONFIG.SYS), mineral phases (CONFIG.PHA) and container materials (CONFIG.CON). The number of names given in these files is limited, while the current not configuration given in the NOTES.TXT file supports exactly 14 different containers, 36 mineral phases, and 42 rock types. Thus, by using this specific configuration program option (see chapter 4.1) one can automatically customize the INFOREX system to work with any dataset.

The D:\INFOREX\DATA.97\ directory contains also a binary **SEARCH.CNF** file to store a current information on the searching configuration.

	ONFIG.	PHA	file	(STANDARD)
1	LIQ		LQ	Glass(melt)
1	PLAG		PL	Plagioclase
	OLIV		OL	Olivine
1	AUG		AU	Augite+Di
1	PIG		PG	Pigeonite
1	CPX		CP	Clino-Px(*)
1	OPX		OP	Ortho-Px
1	PRPX		PR	ProthoPx(*)
1	ILM		IL	Ilm-Hem_sol
1	MAGN		MT	MagnUlv_sol
1	SPIN		SP	Spinel
1	GARN		GR	Garnet
٧,	NEPH		NP	Nepheline
<pre>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>></pre>	LEUC		LC	Leucite
٧,	MEL		ML	Melilite
٧,	AMPH		AM	Ortho-Amph
1	HORN		HB	Clino-Amph
٧,	MICA		MI	Mica_(*)
٧,	ALFS		AF	AlkFeldspar
1	APAT		AP	Apatite
٧,	ARM		AR	Armalcolite
<pre>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>>></pre>	CARB		CB	Carbonate
ν,	EPID		EP	Epidote
٧,	ALSI		AS	Alumina-sil
7	QUAR		QU	Quartz
-V	PERV		PV	Perovskite
N	PSBR		PB	Ps-brookite
N N	SULF		SF	Sulfide
N	CAPV		Pc	Ca-perovsk
2/	MGPV TOUR		Pm	Mg-perovsk
N N	CRND		TU CO	Tourmaline Corundum
1	ZRCN		ZR	Zircon
1	MONT		MO	Monticellit
1	RUT		RU	Rutile
1	SPHE		SH	Sphene
_	CORD		CD	Cordierite
_	STRL		ST	Staurolite
_	MNZT		MZ	Monazite
	QUEN		QP	Ouenched
_	FEBS		QР FB	FeBustamite
	MET		ME	Metal
_	LLIQ		LL	Second_LIQ
				2200110_112

<u>File HELP.TXT</u> was designed to support online help in the data management system (see Appendix 1). It is located in the main INFOREX directory and is written in the rigorously defined format so we do not recommend changing anything in the file.

The user can get help anywhere in the UI using [F1] key. "Help" informs you how to work with the data management system, explains assignements of all active keys in option windows. You can scan the information by [PgUp] and [PgDn] keys; beep signals mean you are at the top or bottom of the HELP.TXT file.

To return to the user interface you should press [*Esc*] key (see prompt at the left bottom angle).

```
CONFIG.CON file (STANDARD)
 PTL Pt loops
 GRA
      Graphite caps.
      Fe_(±Pt_sat)c.
Mo_capsule
 FEC
 MOC
 PTA
      Pt-Ag capsule
       Pt_capsule
  PTC
 AGP
       Ag-Pd capsule
 AUC
       Gold capsule
  IRL
       Iridium loops
       Peridotite cap
  PER
       Al203 crucible
 COR
 REC
       Rh capsule
 PRH
       Pt-Rh_capsule
\sqrt{}
 AUP
       Au-Pd capsule
 CUC
       Cu capsule
```

<u>File NOTES.TXT</u> explains all abbreviations of INFOREX database files (see *Appendix 1*). It contains codenames for igneous rocks and systems, names of phase composition files and their codenames in CONDIT.EXP file, abbreviations for containers and oxygen buffers, parameters of oxygen buffers equations describing temperature and pressure dependence of oxygen fugacity, as well some additional comments.

You can navigate through the on-line "Notes" exactly as through the "Help" except the calling key is [**F2**] instead of [**F1**]. One should not change the NOTES.TXT file.

2.1.3. Auxiliary files

The auxiliary files facilitate the manipulation of data and data swapping among the basic data files using the UI. Usually, the auxiliary files are deleted after each working session, but some of them are stored and overwrited during other sessions.

Few examples of the most important auxiliary files are listed in Table 1 and will be discussed in detail in the following chapters.

Table 1. The most important auxiliary files

File Name (s) Initially present Creating Option

File Name (s)	Initially present	Creating Option	After Work
DATA0INF DATA1INF DATA2INF	NO	SRC	Overwrited
DATA3INF ASS_REF.INF			
#	NO	SRC	Overwrited
BASE.ERR	NO	SPA	Overwrited
DATA2)\$(DATA3)\$(NO	SPA	Deleted
PROCESS.CNF PROCESS.DAT COEXIST.DAT REGRESS.OUT	NO	Mineral-Melt Geothermometers	Overwrited
FE2FE3	NO	Fe3+/Fe2+ ratios	Deleted
H2OMODEL	NO	Water Solubilty	Overwrited
TERNARY	NO	Ternary Diagrams	Overwrited

2.3. Directory structure

The root INFOREX directory is initially named *D*:\INFOREX\, where "*D*" represents a disk drive. It contains a set of executable files to be run by the main management program called INFOMAN.EXE as well as some management system files. Basic experimental data files, BIBL.TXT file and configuration files are located in the *D*:\INFOREX\DATA\.. subdirectory.

Two other subdirectories called *D*:\INFOREX\EXCHANGE\ and *D*:\INFOREX\SAVE\ begin as empty ones after installation. The first one is created for the *IMPORT/EXPORT* operations (*Exchange INFOREX Data* option), which provide an easy exchange of the experimental data files between different users.

The second subdirectory serves mainly to store datafiles resulted from the searching through the INFOREX database (e.g., files in the PROCESS.DAT fashion), as well as output files incuding results of thermodynamic processing and some other petrological calculations. The GENERAL.INF file (see comments to Fig. 6) and files created to store configurations of searching throughout the INFOREX database (*.SEL) will be also placed into the *D*:\INFOREX\SAVE\.. directory.

<u>Users can change the names of these subdirectories, but must keep its general structure intact.</u>

Chapter 3. Setting Up INFOREX

The Setting Up procedures include installation, start and configuration of the INFOREX system. During the installation one should create a directory structure and transfer properly all of the working files to a hard disk (for this purpose we provide an installation program). After installation and running INFOREX, if necessary, one can change a standard configuration of the database, icluding modifications of igneous rock indicators, set of minerals, etc. (see *chapter 4.1*).

<u>Note:</u> At least 1 MB free conventional memory should be available on your computer to work properly with the INFOREX-4.2 database (see *Troubleshooting*).

Before installation of the INFOREX system, make sure that the disk you are going to install INFOREX on does not contain the directory ..\INFOREX. If this subdirectory exists on the disk, rename or remove it.

Installation

To install INFOREX on your hard disk, insert your floppy disk *INFOREX* in drive A or B, make it active and type in the command line:

INSTALL D:

where *D* is drive letter (C, D, E...); note the space between INSTALL and the hard disk name, and colon sign after the name.

The installation program begins to create the INFOREX directory structure and to transfer the system files into the above mentioned subdirectories. After the installation is completed, pull out the floppy disk and start INFOREX.

Starting INFOREX

To start the INFOREX database you should run the main program manager called INFOMAN.EXE located in the main *D*:\INFOREX\... directory. Make the *D*:\INFOREX\ directory the current one, type in the command line **INFOMAN** and press "Enter".

Chapter 4. Description of User Interface

If you first started INFOREX you were find yourself in the **Main Menu** window, including 7 main INFOREX options (Fig. 4). All of the options can be called by pressing the [*F1*] to [*F7*] keys, provided the INFOREX User Interface is hot-key driven, with function keyboard keys serving as the hot keys. Detailed descriptions of working sessions with the functions are followed below.




Figure 4. Structure of the INFOREX-4.2 (2001) data management system

<u>Note:</u> To scroll through a page (*frame or window*) corresponding to a selected *Main Menu* option one should use the arrow keys, whereas for going to another page within a function use the [PgUp] and [PgDn] keys.

Press the [Shift]+[F3] keys, if you want to come back to the initial screen or to a previous window.

4.1 System Configuration (SC)

Feedback on use of the earlier versions of INFOREX showed us that working with smaller subbases is sometimes more useful than with the entire database. Moreover, it was necessary to develop a facility which allows users to change or to input new mineral abbreviations, igneous rock or synthetic system indicators, indecies of other oxygen buffers and containers. To address all of the problems *the System Configuration option* has been developed.

<u>Note:</u> This option has been developed for users experienced in the practical applications of the INFOREX system. Therefore, at an initial stage of using the INFOREX we strongly recommend starting with a standard INFOREX configuration provided during the installation process (see *chapter 3.1*).

So, if you are already familiar with the INFOREX system and really need to reconfigure the database, run INFOREX and press [*F1*] key to choose the *SC* option from the *Main Menu*.




Figure 5. Screen image of the System Configuration option

A screen image will appear at your monitor (Fig. 5), including five windows "Paths", "Phases", "Systems", "Containers", and "Settings" which could be edited during a reconfiguration, as well as the "Buffers" window just demonstrating the oxygen buffers included. One can switch between the windows using the [F3] (Window - Go Forward) or [Shift]+[F3] (Go Backward) keys.

<u>Note:</u> For all but EXCHANGE INFOREX DATA functions the use of the [*F1*], [*F2*] and [*F9*] keys is similar:

[F1] (HELP) provides the on-line help which displays and comments all specific options of INFOREX-4.2 (2001); [F2] (NOTES) involves 5 pages with tables deciphering indicies, search words and flags used in the database; [F9] (QUIT) quits the current operation.

Paths window

This is the first active window of the SC option designed to define paths to the four working INFOREX subdirectories (see *Directory structure*). To make some changes, just type the new ones over the default paths and press "Enter". Go up and down (*only if you are no typing a new path !*) within the window by [\uparrow] or [\downarrow] keyboard keys.

Phases window

To get this window, just press the [F3] key. Now you can add a new mineral phase name in the default set or replace an existing one with another. To do this you should first disable an old mineral phase name (in which you probably do not need during a forthcoming work with INFOREX). It is necessary to do this because the current version of INFOREX-4.2 requires that the total number of *phases* should not exceed 36 for each working session.

So, choose the name to disable and press the [F5] (On/Off) key. As a result, the enable flag ($\sqrt{}$) disappear and the number in the lower right corner of the window is decreased by one. After that press the [F4] (New) key to call input fields as follows:

Enter file name, short name, and commentary for the new phase:

TEST TS Testing-PHA

Type sequently the name of a new file where microprobe data for the new phase will be stored (3-4 symbols, *e.g. TEST*) and a short name of this phase which will be referred to (exactly 2 symbols, *e.g. TS*), followed with a commentary (usually a complete mineral name, up to **11** symbols).

<u>Note:</u> the commentary both here and within the other windows <u>must be</u> an inseparatable single word (use "_" or "-" symbols to connect separate words, e.g. Testing-PHA).

After completing the operations the new phase name will appear in the list of phases just after the old name is disabled. You can change any number of default names in a similar fashion: *monitor the total number in the lower right corner which should be equal to 36.*

Systems window

To get the window press the [F3] (Window) key again. Here, one can create a new rock type or synthetic system indicators with a commentary (usually, a complete name of igneous rock or a commonly used silicate system abbreviation, up to 11 symbols). All of the input operations are absolutely the same as for the *Phases* window - the input fields look like the following:

Enter short name and commentary for the new system:

NEWS NEW-SYSTEM

<u>Note:</u> the length of the new system name must be 3 or 4 symbols, whereas the total number of system indicators should not exceed 42.

Containers window

This window is designed to display or to change the working set of containers included into INFOREX-4.2 database. Container materials are known to be an important feature of experimental conditions, because they can control oxygen fugacity or change the bulk chemical composition of the initial charge. To allow users to filter experimental data based on the materials, we have included in the database a list of commonly used containers (see NOTES.TXT in *Appendix 1*). If there is a need to input a new container name, you should follow the rules as described above. Two input fields are shown in the screen message:

Enter short name and commentary for the new container:

NEW NEW-CONTAINER

Keep in mind that the input field for the container name is of exactly 3 symbols, whereas the commentary field can contain up to **14** symbols. The total number of containers must also be equal to **14**.

Settings window

This window is designed to save your custom settings to a separate file or load into INFOREX previously saved settings.

<u>To save</u> your new settings press the [F7] (Save) key. If you are going to save the settings as a new file, press [F4] (New) and type the new name (without extension!) in the frame that will appear in the screen:

Enter the new name: PLAG-LIQ

and press "Enter". If you wish to store the settings in a file that already exists on the disk, select this file (turn the enable flag of the selected file name on and that of other file names - off by the [*F5*] key) and press the [*F7*] key again.

<u>Note:</u> all of the changes described above are always written directly into the family of INFOREX configuration files (CONFIG.CNF, CONFIG.PHA, CONFIG.SYS, CONFIG.CON) and during other working sessions will be called automatically as default settings.

In order to come back to a standard configuration provided with INFOREX just reload the STANDARD.CNF file.

If you wish to work with other previous settings, load a corresponding ***.CNF file, which has been created earlier. <u>To load</u> previously arranged settings, get the Settings window by the [F6] (Load) key, disable a current configuration file name (CONFIG) and enable the name of settings you are going to load from the file list. Then press the [F6] key to see if the loaded settings in the Paths, Phases, Systems, and Containers windows are those you need.

Buffers window is now inaccessible: the ability to change the list of oxygen buffers currently available in INFOREX-4.2 will be present in the forthcoming versions of the database.

4.2. General System Information (GSI)

Press the [F2] key in the *Main Menu* and wait for a moment until two windows appear in the frame called *General System Information*.

<u>The left GSI window</u> is a passive informational one serving to display the number of master records currently available in the database (written in the CONDIT.EXP file - Fig. 3.2). These records are subdivided in accordance with the *Natural/Synthetic* or *Dry/Wet* system indicators, as well the *Presence/Absence* of volatiles, and run durations.

<u>The right GSI window</u> is an active one, incuding a list of filenames for experimental glasses and minerals currently available in the database and enabled during a working session (the numbers to right of the filenames inform how many records are present in the file considered - Fig. 6). These files can be browsed or printed out using the [**F5**] functional key.




Figure 6. Screen image of the **General System Information** option This image can be saved as GENERAL.INF file or printed out using the F3 or F4 keys

This is an important informational option which allows one to understand how the INFOREX database is organized and to see at the monitor different types of records used to develop the informational system. After pressing the [*F5*] (*Browse*) key one can see a small menu asking for 3 kinds of information you can browse now:

Browse what ?
Phase compositions
Run Conditions
References

Use the $[\uparrow]$ or $[\downarrow]$ keys to select an informational object and press "Enter" to activate a special internal viewer which allows one to scroll a datafile (e.g., CONDIT.EXP), to print out a part of the viewed information or to conduct a search for a text symbol(s) in the file considered (some other options are shown at the bottom of the viewer window).

If you selected *Phase compositions,* the first name in the list of phase composition files will be highlighted, so that one can choose any file to browse navigating using the arrow keys and pressing "Enter".

If you selected *Run Conditions* or *References*, you can browse the CONDIT.EXP and BIBL.TXT files. During scrolling the BIBL.TXT file one can search for the names of authors whose publications or unpublished experimental studies were included into INFOREX database, as well to search for any kind of specific petrological (geochemical, geographic, etc.) terms based on given keys words.

To initiate the searching option, press the [F7] key in the Browse mode and a prompt will appear:

```
Search for: Grove
After a string is found press PgDn
Then F7 to search for the next string
```

To conduct a search, type a word of interest (e.g., *Grove*) and press "Enter". If a string including the word will be found in the BIBL.TXT file, the in-built viewer will highlight the string(s) and stop the operation. To continue the searchig, one should press the [*PgDn*] key, followed by the [*Shift*]+[*F7*] combination, if you are going to search for the same symbols. Press the [*F7*] key to search for a new entry.

After completing the operations with the *Browse* option press "Esc" (see prompt at the bottom) to return to the *GSI windows*. Then press the [**F9**] key to quit the *GSI* option.

4.3. View / Edit / Input Data Files (V/E/I)

As mentioned above, the basic unit of the database is an individual experimental study recorded in the BIBL.TXT file and linked with master records in the CONDIT.EXP file (Fig. 2). The *V/E/I* option of the Main Menu offers the user an ability to work with the main database units. Using this option one can scan the bibliography records page by page or select them at once.

4.3.1. The use of View option

After you pressed the [F3] key of the Main Menu you will see a frame including two windows (Fig. 7). The upper one is active and contains the first reference available in the INFOREX database, whereas the second one, called "Experimental conditions and resultant phase assemblages", is passive inculding only empty fields below the title. One can use the [PgUp] and [PgDn] keys to scan the bibliography records page by page, followed with pressing "Enter" to select an experimental study.




Figure 7. Screen image of the View option at an initial stage

If you need to quickly find out a reference based on its INFOREX number (which might be known, if you have browsed the BIBL.TXT file using the *General Information* option), just press the [*F2*] (*Go to*) key and type the number in a small window, e.g.

Reference Numb: 31

After few seconds of searching through the INFOREX database, the UI will display master records of the CONDIT.EXP file corresponding to a selected experimental study (Fig. 8).




Figure 8. Screen image of the **View** option after selecting a reference

Simultaneously, 2 functions from the initial frame will be changed: **F2:Go to** is transformed to **F2:Notes**, and **F3:Input** is transformed to **F3:Next**. Use the [**F3**] key to get compositional files concerned with the selected study: starting materials or phase composition records linked with the main master records will subsequently appeare on the screen. Calling for the last mineral composition file will automatically result in a sound and appearence of the string:

No more phase compositions! Press any key to continue.

One should then press any alphanumeric key either the [**F9**] key to return to the *Main Menu*, or [**Shift**]+[**F3**] to return to the upper "Reference" window.

Using the [F5] (Write) key, one can specify a file name to store different master and compositional records at any stage of the operations: this option was developed to allow the user to extract specific subsets of experimental information available in the study of interest. If you press the [F5] key a small window will appear:

where "N" is a system prompt: we suppose you are going to save these data as a file named "N" (where "N" corresponds to the INFOREX reference). Sure, one can overtype any appropriate name. To printout the information requested just type *PRN* and press "Enter".

4.3.2. The use of Edit option

Sometimes, the user may encounter an error in the database, e.g. during performing a search for mineral-melt compositional pairs (see *chapter 4.5*, BASE.ERR file), or errors due to misprints in the source papers. To correct the problems or to make changes in the selected study one should use the *Edit* option. To initiate the *Edit* option, press the [*F4*] (*Edit*) key at any stage of the use of the *View* option. A highlighted field will appear in the upper left part of the lower window: move the highlighted field to the field of interest with arrow and *PgUp/PgDn* keys and make all of the necessary changes, followed with pressing "Enter" after each correction (Fig. 9).

Numb		Refere	ence: A						Vol.,	Numb.	Reference: Authors. Title //Source, Year, Vol., Numb., Pages. Umb Grove T.L., Juster T.C. Experimental investigations of low-Ca pyroxene														
67	stabi: basalt	lity an	nd oliv	vine-py sitic l	roxene iquids	e-liqui	d equi	libria	at 1-	atm i	Ca pyron natur	al													
Contents of major components in LIQ - phase, wt.%																									
Num	SiO2	TiO2	A1203	FeO	MnO	MgO	CaO	Na20	K20	P205	Cr203	Н2О													
1 56.20 1.27 15.20 8.76 0.19 5.43 7.80 3.42 1.16 0.19 0.05 0. 2 57.20 1.31 14.50 8.74 0.17 4.87 7.58 3.52 1.34 0.17 0.04 0. 3 56.90 1.59 13.90 9.15 0.21 4.49 7.53 3.55 1.43 0.22 0.03 0. 4 57.50 1.61 14.00 8.99 0.21 4.10 7.09 3.62 1.63 0.25 0.05 0. 5 57.40 1.74 1 9.38 0.21 3.99 6.80 3.62 1.69 0.31 0.06 0. 6 57.70 1.75 13.90 9.41 0.17 3.69 6.79 3.75 1.77 0.23 0.03 0. 7 57.80 1.87 13.70 9.35 0.20 3.41 6.33 3.82 1.94 0.30 0.04 0. 8 57.30 1.97 13.60 9.63 0.19 2.99 6.21 3.96 2.23 0.33 0.10 0. 9 58.00 2.25 13.50 9.22 0.21 2.77 5.91 3.86 2.42 0.29 0.10 0.												0.00 0.00 0.00 0.00 0.00 0.00 0.00													
F1:F	Help	F2:Not	tes	F3:Nex	t.	F4:Edi	.t F	5:Writ	e F9	:Quit	PgUp/E	gDn .													

Figure 9. Screen image of the **Edit** option after selecting a reference and scanning the LIQ.EXP compositional datafile

When you press the [F3] (Next) key in order to go to the next dataset available for the selected experimental study the system will ask if you want to save the corrected data:

Data have been changed Save these data ? Yes No

Your choice! Press the [**F9**] key to quit the operation.

<u>Note:</u> as a rule, our team provides a ternary personal control of the experimental information to be input into the INFOREX database.

4.3.3. The use of Input option

This option works as an in-built editor developed to introduce into INFOREX a new experimental information from the literature. In fact, its use is convinient only if you have a small batch of experimental data (no more than 10-15 experiments with 10-20 phase compositions), because the Input option does not allow one to save uncompleted datasets or to continue the loading at another time. This is the main reason why we recommend users to call for the EID-option to input new experimental data (see the chapter 4.6). Nevertheless, a detailed description of the Input option is presented below.

To initiate the *Input* option one should press the [*F3*] (*Input*) key in the initial *View* option frame (Fig. 7). After pressing the [*F3*] key, the UI gives the new number (*Last+1*) to the dataset to be entered and clears the input fields of the "*Reference*" window. A full blinking cursor appears therein: this means a new bibliography information is requested (at this stage the system is in overwrite mode, therefore the cursor responds only to the arrow keys). All other *Input* options are disabled now but the *F3* one called "*End of reference Input/Edit*": type the reference in any format (*no more 4 strings!*) and press the [*F3*] key. After that there will be a question if you want to save the new bibliography information:

n rows have been input Save these data ? Yes No Choose "Yes" by pressing "Enter", if you are sure all is correct with your new reference. Same requests will appear in the screen when you are working with the master or compositional records.

<u>Note:</u> the input fields are usually highlighted on the screen, therefore if you are going to input a number or a string that is shorter than the input field you should press the "Enter" key to finish the operation. If the string to be input is equal to the length of the highlighted field the input control system will automatically close the current operation followed with the movement to the next input field.

After finishing the "Reference" operation the input highlighted field will be moved to a set of empty fields corresponding to those of master records in the CONDIT.EXP file. All of the fields are tabulated allowing one to type the run numbers, system indecies, pressures, and other experimental conditions (Fig. 10). The system automatically controls the number of units ("1") presented after each phase indicator (Lq1, PL1, etc.), sending in the UI memory the number of empty strings which should be filled later in the phase composition datafiles.

Nielsen R.L. Experimental study of high temperature plagioclase-melt equilibria in high magnesium tholeiitic system (unpublished dat 1992). Experimental conditions and resulted phase assemblages	Numb		Re	efer	en	ce: Au		of expe				Vol	., N	umb.	, Pa	ges.	
N n SystemN/V P,kb Temp lgf02 Dur,hr Con Phase assemblages 102 1 BAS - N 0.0 1250 QFM-1.0 72.0 PTL LQ1 PL1 OL 2 BAS - N 0.0 1260 QFM-1.0 60.0 PTL LQ1 PL1 OL 3 BAS - N 0.0 1270 QFM-1.0 60.0 PTL LQ1 PL1 OL 4 BAS + N 0.0 1270 QFM-1.0 60.0 PTL LQ1 PL1 OL 5 BAS + N 0.0 1280 QFM-1.0 60.0 PTL LQ1 PL1 CQ1 PL1 6 BAS + N 0.0 1290 QFM-1.0 48.0 PTL LQ1 PL1 CQ		Nie mel	lt eq			-		_	_		-		-	_			ì,
102			Exp	peri	me	ntal c	conditi	ions and	resulte	ed pl	nase	ass	embl	ages			
2 BAS - N 0.0 1260 QFM-1.0 60.0 PTL LQ1 PL1 OL 3 BAS - N 0.0 1270 QFM-1.0 60.0 PTL LQ1 PL1 OL 4 BAS + N 0.0 1270 QFM-1.0 60.0 PTL LQ1 PL1 S BAS + N 0.0 1280 QFM-1.0 48.0 PTL LQ1 PL1 S BAS + N 0.0 1290 QF 0.0 7 - 0.0 0 0 0.0 0 0.0 9 - 0.0 0 0.0 9 - 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0 0.0 0.0 0 0.	N																

Figure 10. Screen image of the **Input** option during work with main master records This is an example of 1992 where only 101 studies were included into the INFOREX database

The *Input* control system identifies strings that have been ready for the automatic input by the first nonempty field in the phase assemblage codenames area for the run conditions master records, and by first 3 nonzero positions for the phase composition records (it means at least 3 first positions in each highlighted field should be filled in).

Note: most of the *Input* fields are valid only for numbers or alphanumeric strings, and if you press an alphabetic key in the numeric field you will hear a beep and nothing appears on the screen. So, be sure you press permissable keys.

But some fields allow both alphanumeric string and numbers to be input, for example the "IgfO2" field of CONDIT.EXP file permits both a number like "-12.58" or a string like "QFM-1.2". If you mistakenly type the "number" -12.58 instead of -12.58, this error will not be detected here by the *Input* control system and later the program will interprete the "number" -12.58 as -12.00.

In the case of large batches of experimental data (of about 50 experiments with hundreds of microprobe analyses) there is an ability to load these data in two stages. To do this one should first complete the fields of "Phase Assemblages" for all master records to be entered. Due to such a preliminary operation the Input control system "knows" how many strings need to be reserved for the input of chemical compositions in each of the phase composition files. So, now one can conduct an input of "uncompleted" datasets including major element contents, e.g. by filling in only the first three columns (SiO₂ contents) within the compositional strings, followed with saving these uncompleted datasets. Such a "template" can be completed later during a second stage of the input using the Edit option.

Press any alphanumeric key followed with either [**F9**] to return to the *Main Menu*, or [**Shift**]+[**F3**] to return to the "Reference" window.

<u>Note:</u> Beginninig from the INFOREX-2.0 release the data input is rather easier due to a new ability to *Mark* and Copy the whole input string(s), if there are a few similar strings in your source data. This option is now available for the "*New reference input*" mode only.

To copy string(s) in the "New reference input" mode you should place the cursor at the first string to be copied (using the arrows and PgUp/PgDn keys) and press [F4] (Mark) key (we assume that you have already pressed it once to initialize the Edit routine and you are seeing now the full-block cursor in your screen). A blinking "Mark" label and bell indicate that this Mark-option is activated. Use arrows and PgUp/PgDn keys to highlight strings to copy and press [F4] key again to fix the highlighted block. Place the cursor before the string where you are going to insert the block and press [F4]: the block will be inserted before the string. Press [Esc] to cancel copying at any stage of the operation.

Again we emphasize, that our own experience in the practical work with the INFOREX database has shown that the use of Export/Import operations is much more preferable to input new experimental information.

4.4. Selecting Run Conditions (SRC)

This is one of the most important INFOREX functions oriented to the following petrologic applications. The *SRC* option is implemented to perform search-match operations throughout the database according to some "filters". The groups of the filters are:

- 5 general types of silicate systems investigated, such as Terrestrial, Lunar, Meteoritic, Planetary (Venus, Martian, etc.), and Synthetic
- names of igneous rocks (BAS basalts, KOM komatiites, SYN synthetic, etc. see NOTES.TXT in Appendix 1)
- combinations of volatile components (CO₂, H₂, H₂O), including an ability to select fluid saturated and/or undersaturated experiments
- flags of the presence or absence of a direct chemical information on water contents and Fe³⁺/Fe²⁺ ratios in experimental glasses
- types of containers used in the experiments (*PTL* Pt loops, *AGP* Ag-Pd capsules, etc. see NOTES.TXT in *Appendix 1*)
- experimental conditions, including the range of temperatures, pressures, run durations, and oxygen fugacities)
- the range of major component contents or some routinely used petrochemical ratios for available experimental glasses (here is also a new ability to specify a set of glasses corresponding to igneous rock series based on the SiO₂ - Na₂O+K₂O diagram).




Figure 11. Screen image of the Selecting Run Conditions option

The screen display for presenting the *SRC* option is given in Figure 11. It contains a set of countered boxes: transition from one box to another is by means of pressing the **[F3]** (*Next_box*) key (the **[Shift]**+**[F3]** combination allows one to come

back to a previous box or window). To setup a parameter to filter data one should move the highlighted field to the parameter field and press "Enter" or type a number.

Note: All items within the first 4 groups of the filters mentioned above are related to each another by logical "OR", whereas logical "AND" is the relation between the main filter groups. Items of "Conditions" and "Series/Liquid Parameters" are related by logical "AND", just like the groups itself. It means if you selected the runs with Lunar and Meteoritic materials you will find the CONDIT.EXP master records including both eucrites (EUC), SNC meteorites (SNC), lunar basalts (BAS), green glasses (GGL), and others. But if you have used an additional restriction, such as SNC in the "Igneous Rocks" window, you will find only experiments on SNC meteorites.

An attractive feature of the SRC option is the ability to specify the range of redox conditions of interest. To perform this operation, one should specify a type of oxygen buffer and the range of oxygen fugacities with respect to the given buffer ($log\ units$). To give the buffer name, move the highlighted field to the " $fO2\pm\delta$ " field in the "Conditions" window and press the [F4] ($fO2\ Buffer$) key. As a result a window will appear in the screen:

F=[Esc=Exi	t]====[Enter=Select]===
IW	Myers, Eugster, 1983
MM	Myers,Eugster,1983
MI M	Huebner J.S.,1971
MH	Myers,Eugster,1983
QFM	Myers, Eugster, 1983
IQF	Myers, Eugster, 1983
NNO	Huebner J.S.,1971
cco	Myers, Gunter, 1979
MMO	Mn304 -Mn0
COC>5kb	Woermann et al.,1977
COC<5kb	French B.H., 1966
GRS	Graphite-COH
GCH	Graphite-methane
	RESERVED —
[Select	: buffer by ↓/↑ keys] —

Use the arrow keys to select a needed buffer by a highlighted field and press "Enter": selected buffer name appears at the place of **fO2** label in the "Conditions" window.

<u>Searching throughout the database</u> can be started from any field in the screen: press the [**F5**] (**Search**) key and wait for a few seconds until the **Search** system informs you that 3 files, characterizing the experiments selected, are created in the D:\INFOREX\SAVE\.. directory:

These auxiliary files involve short references, a list of phase assemblages and master records for all individual runs satisfying the previously defined criteria. Those files (as well some others - see *Table 1*) serve as temporary data files for the following *SPA* procedure. If no data match your restrictions, another window will appear in the sceen with the mismatch information.

Beginning from the INFOREX-4.0 (1997) release, the SRC option includes the ability to save configurations of the searching to a separate file or load some previously saved configurations from the screen. For the purpose a special function "File" was developed. We recommend using this function first before calling for the Search procedure. During use of Setting Phase Assemblages option the INFOREX UI creates output files (e.g., including olivine-melt or spinel-garnet equilibria information see below) which in the first string contain a file name of the "Search" configuration, as well a date of conduction the search. So, the function "File" allows one to find out a previously used "Search" configuration, if you have forgotten what kinds of filters and restrictions were used previously.

To save or to load a "Search" configuration, press the [F6] (File) key from any field in the SRC screen - a small menu will appear on the monitor:

Read Write Delete

Select any of the functions by the *PgDn/PgUp* keys and press "Enter". A place for writing the "*Search*" configuration, as well a list of previously saved setting will appear in a frame:

Type in file name or
Select one from a list

TEST
SNC-PLAG
OL-GARN
AUG-MELT
WATERSAT

Type a new file name (e.g., *TEST*) if you are going to save (*Write*) a "*Search*" configuration, or select another one if you want to load (*Read*) previously saved settings. Each name written using the "*File*" function will be saved as *.SEL file (e.g., *TEST.SEL*) in the D:\INFOREX\SAVE\.. directory. Your choice!

<u>To perform a new search</u>, one should press the [F9] (Quit) key to return to the Main Menu, followed by pressing the [F4] key to call for the SRC option again.

4.5. Setting Phase Assemblages (SPA)

This option is designed for the purpose of creating phase equilibria compositional files which could be used in further thermodynamic processing and analysis of experimental information. *It works only if the SRC option has created a set of auxiliary files satisfying some previously defined criteria (Table 1)*.

Three windows appear on the screen when one calls the *SPA* option from the *Main Menu* by pressing the [*F5*] functional key. In the first, all phase assemblages are listed that have appeared in the data subset search-matched in the *SRC* procedure, whereas the second contains a list of shortcut references to experimental investigations where these data were described (Fig. 12).




Figure 12. Screen image of the Setting Phase Assemblages option (after selecting Olivine-Melt equilibrium)

Using the [F3] (Window) key one can switch between these upper two windows. Pressing the [F4] (On/Off) key followed with the [F5] (Exclude) key allows one to exclude from these lists any phase association or an entire experimental investigation. This is convenient especially if only a single multiphase assemblage is required, such

as locating the compositions of melts equilibrated with olivine, plagioclase, pyroxenes, and spinel.

<u>Note:</u> The *On/Off* function sets the disabling flag " $\sqrt{}$ " in front of a selected phase assemblage or a shortcut reference: all strings marked to disable can be excluded from further processing by the *Exclude* function. <u>Important:</u> none of the erased strings can be enabled with the *SPA* option again! One should repeate a previous *SRC* session to correct a disabled mistake if necessary.

When the desired arrangement is completed, the user can proceed to the third window containing the full list of phases (experimental products). It can be done by pressing the [F9] (Phases) key - a small menu appears on the screen:

Mineral-Melt Equilibria Two mineral Equilibria Only melt compositions

Significance of these functions with some practical examples will be given in the following paragraph.

If you have selected "Mineral-Melt Equilibria" or "Only melt compositions" (followed by the "Melts Equilibrated with a Given Mineral" function) the SPA option allows the user to select one of the mineral phases in the third window (assumed that liquid is to be the second coexisting phase). A highlighted field appearing at the first name field should be moved by the arrow keys to the name of interest. Press "Enter" and in a few seconds the number of mineral-melt coexisting pairs extracted will be displayed on the screen (Fig.12); simultaneously a prompt will appear:

Press any key to continue...

Press "Enter" and a message will inform you that a corresponding file has been created and stored in the D:\INFOREX\SAVE.. directory. For example, if at this stage you selected OPX this message will look like:

File OP_LQ.DAT containing phase compositions has been created. Press any key to return to main menu

<u>If the "Two mineral equilibria" function has been selected</u>, move the highlighted field to the name of first phase requested and press "Enter", then do the same for the second phase. The system prompts will look like those described above.

If you have selected "Only melt compositions" (followed by the "Melts Equilibrated with All Phases" or "Melts To test Fe3+/Fe2+ Ratio" functions) the SPA procedure will also show how many liquid compositions were found and a corresponding file (see below) created.

```
OP LQ.SEL (03/21/96 16:50:12)
      2 | SNCM- 2 | N | 1 | CHNM- 1 | N |
                                             49.0 PTL LQ1 OL1 OP1 SP 72.0 PRH LQ1 OL1 OP1
                      0.0
                                    -7.70
                            1270
166
                      0.0
                            1350
                                   -11.20
      2 CHNM- 1 N
166
                      0.0
                            1328
                                   -11.43
                                             92.0 PRH LQ1
                                                            OL1
                                                                OP1
        CHNM- 1 N
166
                      0.0
                            1300
                                   -11.73
                                            122.0 PRH LQ1
      5
                                                            OL1
                                                                OP1
               1
166
         CHNM-
                 N
                      0.0
                            1300
                                   -11.73
                                             67.0
                                                  PRH LQ1
                                                            OL1
                                                                OP1
        CHNM- 1 N
CHNM- 1 N
                                             97.0 PRH LQ1
                            1275
                                   -12.01
166
                      0.0
                                                            OL1
                                                                OP1
166
                      0.0
                            1220
                                   -12.66
                                            120.0 PRH LQ1
                                                            OL1
                                                                OP1
                                            120.0 PRH LQ1
         CHNM- 1 N
166
      9
                      0.0
                            1200
                                   -12.90
                                                            OL1
                                                                OP1
               1
                                   -13.16
                                                  PRH LQ1
166
     10
        CHNM-
                 N
                      0.0
                            1180
                                            120.0
                                                            OL1
                                                                OP1
     11 CHNM-
               1 N
                      0.0
                            1170
                                   -13.32
                                            168.0 PRH LQ1
166
                                                            OL1
                                                                OP1
166 12 CHNM- 1 N
166 17 CHNM- 2 N
166 18 CHNM- 2 N
166 20 CHNM- 2 N
                      0.0
                                   -13.32
                                            114.0 PRH LQ1
48.0 PRH LQ1
                            1170
                                                                OP1
                                                            OL1
                      0.0
                            1220
                                   -12.66
                                                            OL1
                                                                OP1
                                   -12.90
                                            120.0 PRH LQ1
                      0.0
                            1200
                                                                OP1
                                                            OL1
                      0.0
                            1170
                                  -13.32 168.0 PRH LQ1 OL1 OP1
OPX.EXP
 55.74
                                                  .02
                                                         .00
          .08
                 .50 12.62
                              .37 29.88
                                          1.24
                                                                .00
                                                                       .40
                                                                              .00 77
 55.70
         0.00
               0.29 14.20
                             0.28 28.90
                                          0.40
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.61
                                                                             0.00 166
 55.00
                                                 0.00
                                                        0.00
                                                               0.00
                                                                             0.00 166
               0.37 15.80
                             0.30 27.60
         0.00
                                           0.54
                                                                      0.71
 55.30
         0.00
               0.41
                     16.40
                             0.29
                                   27.60
                                           0.61
                                                 0.00
                                                        0.00
                                                               0.00
                                                                             0.00
                             0.31 27.30
                                                               0.00
 55.10
         0.00
               0.43 16.80
                                           0.61
                                                 0.00
                                                        0.00
                                                                      0.70
                                                                             0.00 166
                             0.33 27.40
 55.10
         0.00
               0.48 17.40
                                           0.78
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.79
                                                                             0.00 166
 54.30
         0.00
               0.66 18.20
                             0.35 25.50
                                           1.38
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.69
                                                                             0.00
                                                                                  166
 54.40
               0.47 17.40
                             0.36 25.60
                                           1.82
                                                 0.00
         0.00
                                                        0.00
                                                               0.00
                                                                      0.43
                                                                             0.00 166
 54.20
         0.00
               0.44 17.80
                             0.37 24.60
                                           2.44
                                                  0.00
                                                        0.00
                                                               0.00
                                                                      0.49
                                                                             0.00 166 10
                                                        0.00
 54.10
        0.00
                                                 0.00
                                                               0.00
                                                                             0.00 166 11
               0.42 18.00
                             0.39 24.10
                                           2.86
                                                                      0.50
 54.50
         0.00
               0.60 17.00
                             0.40 25.70
                                           1.84
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.45
                                                                             0.00
        0.00
                                                        0.00
                                                               0.00
 53.40
               0.70 20.60
                             0.27 24.00
                                           0.90
                                                 0.00
                                                                      0.17
                                                                             0.00 166 17
 53.30
         0.00
               1.17 20.80
                             0.29 22.80
                                          1.42
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.83
                                                                            0.00 166 18
 52.40
         0.00
               1.09 20.90
                            0.31 22.10
                                          2.12
                                                 0.00
                                                        0.00
                                                               0.00
                                                                      0.72
                                                                            0.00 166 20
LIQ.EXP
          .98
 50.47
               7.98 17.38
                              .56 10.98
                                           9.52
                                                  .70
                                                         .00
                                                                .00
                                                                              .00 77
                                           4.52
 52.10
         0.23
               5.64 23.30
                             0.44 14.60
                                                 0.01
                                                        0.00
                                                               0.00
                                                                      0.74
                                                                            0.00 166
 50.90
         0.26
                6.48 23.20
                             0.42 12.90
                                           5.12
                                                 0.02
                                                        0.00
                                                               0.00
                                                                      0.77
                                                                             0.00
                                                                                  166
 50.60
         0.29
               7.21 24.30
                             0.40 12.10
                                           5.69
                                                 0.01
                                                        0.00
                                                               0.00
                                                                      0.77
                                                                             0.00 166
                                                               0.06
                                                                      0.70
 50.80
         0.28
               6.95 24.50
                             0.41 11.80
                                           5.50
                                                 0.03
                                                        0.00
                                                                             0.00 166
 50.60
               8.50 23.40
                             0.41 10.40
                                           6.85
                                                 0.03
                                                        0.00
         0.38
                                                               0.07
                                                                      0.72
                                                                             0.00
 48.90
         0.44 11.90 21.20
                             0.35
                                    8.54
                                           8.59
                                                 0.03
                                                        0.00
                                                               0.40
                                                                      0.44
                                                                             0.00 166
 50.70
         0.47 13.00 17.80
                             0.32
                                    7.31
                                           8.39
                                                 0.69
                                                        0.00
                                                               0.52
                                                                      0.38
                                                                             0.00 166
                                                        0.00
 53.80
         0.44 13.50 16.00
                             0.27
                                    6.13
                                           8.16
                                                 1.29
                                                               0.16
                                                                      0.27
                                                                             0.00 166 10
 49.60
                                    6.70 10.10
         0.50
              14.80 16.80
                             0.30
                                                 0.12
                                                        0.00
                                                               0.14
                                                                      0.29
                                                                             0.00 166
                                    5.80
                                                 1.58
 54.20
         0.45 14.00 15.00
                             0.27
                                           8.00
                                                        0.00
                                                               0.36
                                                                      0.22
                                                                             0.00 166 12
                                    7.97
                                           7.84
                                                        0.00
 49.10
         0.35
               9.47 24.60
                             0.30
                                                 0.05
                                                               0.20
                                                                      0.43
                                                                             0.00 166 17
 50.80
              11.10 22.00
         0.40
                             0.29
                                    6.85
                                           8.11
                                                 0.33
                                                        0.00
                                                               0.14
                                                                      0.43
                                                                             0.00
 52.90
         0.42 12.10 18.60
                             0.24
                                    5.79
                                          8.64
                                                 1.01
                                                        0.00
                                                               0.19
                                                                      0.27
                                                                             0.00 166 20
```

Figure 13. An example of OP_LQ.DAT file used to develop *Opx*-melt geothermometers for meteoritic igneous systems (Ariskin et al., 1997*b*) created using the "*Mineral-Melt Equilibria*" function

Significance of the main "Phases" functions

Use the arrow keys followed by pressing "Enter" to select one of the 3 functions:

- (1) <u>"Mineral-Melt Equilibria"</u> provides an ability to find out an experimental information on the equilibrium of a given mineral with melts (glasses), e.g. orthopyroxene-melt or spinel-melt. Files like **_LQ.DAT (OP_LQ.DAT or SP_LQ.DAT) will be created by the function and stored at the D:\INFOREX\SAVE.. directory (Fig.13).
- (2) selecting <u>"Two mineral equilibria"</u> means you are going to request information on a two solid phases equilibrium, e.g. *olivine-garnet* ot *augite-orthopyroxene*. Files like **_**.**DAT** (*OL_GR.DAT* or *AU_OP.DAT*) will be created by the function and stored at the D:\INFOREX\SAVE.. directory:

```
OL SP.SEL (05/14/97 12:50:37)
158
      3 PERS- 1 N
                    10.0
                                           58.0 GRA LQ1 OL1 OP1 SP1
                                     ND
                                                GRA LQ1
      4 PERS- 1 N
158
                    10.0
                          1350
                                     ND
                                           53.7
                                                         OL1
                                                             OP1
                                                                 SP1
158
      8 PERS- 1 N
                    10.0
                          1330
                                     ND
                                           71.9 GRA LQ1
                                                         OL1
                                                             OP1
                                                                 AU1
                                                                     SP1
                                           72.0 GRA LQ1
51.6 GRA LQ1
158
     14
        PERS-
              1 N
                    10.0
                          1300
                                     ND
                                                         OL1
                                                             OP1
                                                                 AU1
                                                                     SP1
     16 PERS- 1 N
158
                    10.0
                          1280
                                                             OP1 AU1 SP1
                                     ND
                                                         OL1
158
    18 PERS- 1 N
                    10.0
                          1270
                                           97.9 GRA LQ1 OL1 OP1 AU1 SP1
                                     ND
     19 PERS- 1 N
                          1270
158
                    10.0
                                     ND
                                          151.3 GRA LQ1 OL1 OP1 AU1 SP1
175
     48
        PERS+ 0 N
                    15.0
                          1320
                                     ND
                                           48.0 GRA LQ1
                                                        OL1
                                                             OP1
                                                                 AU1
175 | 53 | PERS+ 0 | N | 15.0 | 1335 |
                                     ND
                                          48.0 GRA LQ1 OL1 OP1 SP1
OLIV.EXP
 40.70
        0.00
              0.06
                     8.44
                           0.09 51.10
                                        0.30
                                               0.00
                                                     0.00
                                                            0.00
                                                                  0.29
                                                                         0.00 158
                     8.50
                           0.11 50.80
                                        0.32
                                                     0.00
                                                                         0.00 158
 40.40
              0.05
                                               0.00
                                                            0.00
                                                                  0.25
        0.00
 40.50
        0.00
              0.07
                           0.10 50.50
                                        0.34
                                               0.00
                                                     0.00
                                                            0.00
                                                                  0.24
                     8.60
                                                                         0.00 158
              0.06
                                               0.00
                                                            0.00
 40.80
                     9.10
                                        0.28
                                                     0.00
                                                                         0.00 158 14
        0.00
                           0.14 51.00
                                                                  0.14
 40.80
        0.00
              0.03
                     8.98
                           0.13 49.80
                                        0.22
                                               0.00
                                                     0.00
                                                            0.00
                                                                  0.13
                                                                         0.00 158
 40.20
        0.00
              0.04
                     9.20
                           0.16 50.40
                                        0.25
                                               0.00
                                                     0.00
                                                            0.00
                                                                         0.00 158 18
                                                                  0.12
 40.70
        0.00
              0.00
                     9.14
                           0.11 50.50
                                        0.22
                                               0.00
                                                     0.00
                                                            0.00
                                                                  0.13
                                                                         0.00 158 19
 38.60
        0.00
              0.08 20.60
                           0.34 40.30
                                        0.32
                                               0.00
                                                     0.00
                                                            0.00
                                                                  0.31
                                                                         0.00 175
              0.11 21.90
                                                                         0.00 175 53
 39.00
                           0.37 38.70
                                        0.33
                                               0.00
        0.00
                                                     0.00
                                                            0.00
                                                                  0.15
SPIN.EXP
        0.10 24.00 10.22
                           0.19 18.00
                                        0.11
                                               0.00
                                                     0.00
                                                            0.00 47.60
                                                                         0.00 158
  0.24
  0.23
        0.14 24.50 10.90
                           0.14 18.00
                                        0.15
                                               0.00
                                                     0.00
                                                            0.00 46.80
                                                                         0.00 158
  0.32
        0.17 28.50
                    9.90
                           0.17 18.16
                                        0.07
                                               0.00
                                                     0.00
                                                            0.00 42.10
                                                                         0.00 158
  0.05
        0.09 37.60 10.00
                           0.12 19.78
                                        0.09
                                               0.00
                                                     0.00
                                                            0.00 32.90
                                                                         0.00 158 14
                                               0.00
  0.18
        0.17 44.20
                     9.30
                           0.11 20.00
                                         0.08
                                                     0.00
                                                            0.00 25.00
                                                                         0.00
        0.09 44.20
                                        0.05
                     9.50
                                                     0.00
                                                                         0.00 158 18
  0.19
                           0.11 20.40
                                               0.00
                                                            0.00 25.60
  0.27
        0.14 43.40
                     9.80
                           0.09 20.20
                                        0.09
                                               0.00
                                                     0.00
                                                            0.00 27.00
                                                                         0.00 158 19
                                 11.80
  0.48
        0.66 27.60 22.40
                           0.31
                                        0.00
                                               0.00
                                                     0.00
                                                            0.00 36.70
                                                                         0.00 175 48
  0.49
        0.57 26.60 23.10
                           0.36 11.40
                                        0.11
                                               0.00
                                                     0.00
                                                            0.00 36.60
                                                                         0.00 175 53
```

Figure 14. Example of OL_SP.DAT file characterizing *Olivine-Spinel* equilibrium created using the "*Two mineral equilibria*" function for the experiments conducted at 10-15 kbar on *Dry Synthetic* low-alkaline systems (run duration ≥ 48 hours)

(3) selecting <u>"Only melt compositions"</u> permits one to extract only experimental glass compositions, which can correspond to one of 3 additional restrictions:

```
Melts Equilibrated with All Phases
Melts Equilibrated with a Given Mineral
Melts To test Fe3+/Fe2+ Ratio
```

(3a) "Melts Equilibrated with All Phases" allows one to retrieve all experimental glass compositions, satisfying previously defined criteria: LQ_LQ.DAT file will be created here. This function is especially useful if one is interested in melt compositions representing a given set of phase assemblages, specified earlier in the upper left window (Fig. 12) using the [F4] (On/Off) and [F5] (Exclude) keys.

One can readily create a file including all glass compositions corresponding to a multiphase equilibrium, e.g. *Olivine-Orthopyroxene-Augite-Spinel-Melt*. Such an information might be used in further projectings of the liquid compositions onto ternary diagrams (see *chapter 4.7.5*), or for calibrations of equations describing liquidus temperatures of the mineral assemblage in terms of the liquid composition (*Longhi, 1991; Ariskin et al., 1997a*). An example of such a file is given in Fig. 15.

```
LQ LQ.SEL (05/14/97 14:01:35)
    43 LHR - 1 N 1 LHR - 1 N
                                           17.0 GRA LQ1 OL1 AU1 OP1 GR1 6.0 GRA LQ1 OL1 OP1 AU1 GR1
                    30.0
                                     ND
156
                   45.0
                          1650
                                     ND
    27 LHR - 1 N
                   30.0
                          1470
                                     ND
                                           11.1 GRA LQ1 OL
                                                                  AU
183 28 LHR - 1 N
                    30.0
                          1460
                                     ND
                                           10.0 GRA LQ1 OL
                                                             OP
                                                                 AU
                                                                     GR
LIQ.EXP
        0.90 11.00 7.80
                           0.20 19.20 12.20
                                               1.20
                                                     0.00
                                                            0.00
                                                                  0.40
 46.90
                                                            0.00
        0.00 9.70 10.40
                           0.00 19.00 10.40
                                               0.70
                                                     0.00
                                                                   0.50
 45.90
                                                                         0.00 156
        2.00 12.54 11.20
                            .19 15.24
                                        8.98
                                               2.06
                                                                          .00 183 27
 45.63
                                                     1.88
                                                             .00
                                                                    .28
        2.23 13.09 10.14
                                        7.74
 46.73
                             .13 12.35
                                              2.55
                                                     4.91
                                                             .00
                                                                    .13
                                                                           .00 183 28
```

Figure 15. Example of LQ_LQ.DAT file characterizing a high-pressure Olivine-Orthopyroxene-Augite-Garnet-Melt cotectics created using "Only melt compositions" function coupled with "Melts Equilibrated with All Phases" for the experiments conducted on Dry Terrestrial Harzburgites at 30-45 kbar (run duration ≥ 4 hours)

(3b) "Melts Equilibrated with a Given Mineral" allows one to get experimental glass compositions being in equilibrium with a given mineral, independent of the type of phase assemblages or the presence/absence of microprobe information on the solid phase compositions. Files like **_LQ.DAT will be created by the function and stored at the D:\INFOREX\SAVE.. directory.

For example, one can readily create a file including all glass compositions in equilibrium with hornblende (HB_LQ.DAT), even if for part of the data the hornblende

compositions are unknown. Such an information might be used in further projectings of the liquid compositions onto ternary diagrams (see *chapter 4.7.5*).

Fig. 16 provides an example of such a file characterizing water-saturated melts equlibrated with hornblende (*Ariskin et al., 1997a*).

HB I	Q.SI	EL (05	5/14	1/97	14:4	1:19)											
$11\overline{1}$	36	BAS -	+ 2	F	2.0	800	-14.18	102.0	AUC	LQ1	PL1	SF1	НВ	OP			
111	66	BAS -	+ 2	F	2.0	1000	-6.42	144.0	AUC	LQ1	PL	AH	МТ	CP	HB	OL	
112	11	BAS -	- 1	F	3.0	850	ND	115.0			PL2	HB1					
112	12	BAS -	- 1	F	3.0	875	ND	137.0	AUC	LQ1	PL2	HB1					
112	27	BAS -	- 2	F	3.0	850	ND	115.0	AUC	LQ1	PL2	HB1					
112	28	BAS -	- 2	F	3.0	875	ND	137.0	AUC	LQ1	PL2	HB1			I		
112	46	BAS -	- 3	F	3.0	850	ND	115.0	AUC	LQ1	PL2	HB1					
112	47	BAS -	- 3	F	3.0	875	ND	137.0	AUC	LQ1	PL2	HB1			I		
112	79	BAS -	- 5	F	3.0	850	ND	137.0	AUC	LQ1	PL2	HB1					
133	3	BAS -	- 1	F	5.0	750	QFM	120.0	AGP	LQ1	PL1	MT	IL	нв1	*		
133	28	BAS -	- 2	F	4.9	725	QFM	216.0	AGP	LQ1	OL	PL1	IL1	HB1	SH1	PB1	
133	42	BAS -	- 2	F	4.9	725	MH	216.0	AGP	LQ1	PL1	IL	нв1	PB	SH1	MI	
229	12	BAS -	+ 0	F	2.0	940	QFM	168.0	AGP	LQ1	AU?	нв1	*		ĺ		
LIQ.	EXP	•			•	•											
67.	04	0.11	14.	. 40	1.32	0.03	0.42	1.35	5.14	3.0	9 (0.00	0.0	00	5.95	111	36
53.	54	0.56	19.	. 31	2.50	0.20	0.81	6.90	3.89	0.9	96 (0.00	0.0	00	5.96	111	66
74.	37	0.36	16.	.15	1.76	0.09	0.25	3.25	2.70	1.0)1 (0.09	0.0	00	00.0	112	11
71.	85	0.38	16.	. 95	1.95	0.15	0.56	3.99	3.10	0.8	39 (0.18	0.0	0.0	00.0	112	12
74.	13	0.43	17.	.06	1.68	0.10	0.24	2.38	3.77	0.1	.2 (0.13	0.0	00	00.0	112	27
70.	16	0.42	18.	.77	1.72	0.11	0.28	3.55	4.71	0.1	.0 (0.21	0.0	00	00.0	112	28
74.	29	0.35	16.	.20	1.89	0.11	0.43	3.52	2.57	0.5	54 (0.12	0.0	00	00.0	112	46
71.	63	0.37	17.	.26	1.87	0.15	0.49	4.07	3.54	0.4	13 (0.21	0.0	00	00.0	112	47
72.	68	0.31	17.	.05	1.52	0.10	0.12	2.91	3.59	1.5	6 (0.19	0.0	0.0	00.0	112	79
74.	40	0.10	16.	.00	1.43	0.00	0.00	2.54	3.70	1.8	30 (0.00	0.0	0.0	00.0	133	3
75.	80	0.10	15.	. 80	1.06	0.00	0.00	2.44	2.40	2.4	10 (0.01	0.0	0.0	0.00	133	28
76.	00	0.05	15.	. 30	1.07	0.00	0.04	2.58	2.60	2.4	10 (0.01	0.0	0.0	0.00	133	42
60.	29	1.46	21.	. 44	2.60	0.17	0.29	6.71	5.02	2.0)2 (0.00	0.0	0.0	0.0	229	12

Figure 16. Example of HB_LQ.DAT file characterizing *Hornblende-Melt* equilibrium created using "*Only melt compositions*" function coupled with "*Melts Equilibrated with a Given Mineral*" for the experiments conducted on *Water-saturated Terrestrial Basalts* at 2-5 kbar (run duration ≥ 100 hours)

(3c) selecting "Melts To test Fe3+/Fe2+ Ratio" function makes it possible to retrieve information, including experimental conditions, glass compositions and Fe^{3+}/Fe^{2+} ratios in the glasses, if available.

File called **RE_LQ.DAT** will be created by the function and stored at the D:\INFOREX\SAVE.. directory.

						MANUAL
	BAS - 0 N 0.0	1459 -9.21	5.0 PTL 1			
152 11	BAS - 0 N 0.0	1459 -6.15	5.0 PTL 1	LQ1 RE		
- 1	BAS - 0 N 0.0		5.0 PTL 1			
152 20	BAS - 0 N 0.0	1362 -7.93	5.0 PTL 1			
152 35	1 1 1	1366 -8.91	5.0 PTL 1	LQ1 RE		
152 38			1 11	LQ1 RE		
	BAS - 0 N 0.0	1363 -6.81		LQ1 RE		
152 51	BAS - 0 N 0.0	1363 -6.81	5.0 PTL 1	LQ1 RE		
REDOX.E						
152 1		.095 .00	.00		'e203=1.111	(FeOt-FeO)
152 3		.078 .00	.00	.000		
152 4		.070 .00	.00	.000		
152 6	8.81 .58	.059 .00	.00	.000∥		
152 11		.194 .00	.00	.000		
152 12		.156 .00	.00	.000		
152 20		.102 .00	.00	.000		
152 35		.062 .00	.00	.000∥		
152 38		.049 .00	.00	.000		
152 50		.189 .00	.00	.000		
152 51	9.54 1.96	.185 .00	.00	.000		
LIQ.EXP						
51.36	1.90 14.02 11.3				20 0.02	0.00 152 1
52.57	1.94 14.23 9.4				13 0.02	0.00 152 3
51.76	1.95 14.16 10.7				15 0.02	0.00 152 4
52.03	1.95 14.23 9.3				14 0.02	0.00 152 6
51.28	1.91 13.99 11.3				16 0.02	0.00 152 11
49.97	2.86 14.04 10.5				25 0.04	0.00 152 12
51.11	1.95 14.09 10.6				35 0.04	0.00 152 20
51.93	1.97 14.41 10.7		1.72		56 0.03	0.00 152 35
51.41	3.29 13.64 11.0		2.09		43 0.10	0.00 152 38
52.39	1.96 13.96 11.7				17 0.03	0.00 152 50
52.35	1.98 14.25 11.3	30 0.21 7.11 11	31 2.33	0.17 0.	21 0.02	0.00 152 51

Figure 17. Example of RE_LQ.DAT file including information on Fe³+/Fe²+ ratios created using "Only melt compositions" function coupled with "Melts To test Fe3+/Fe2+ Ratio" for the experiments conducted on Dry Terrestrial Tholeiite-Basaltic systems at 1 atm, 1350-1500°C, and oxygen fugacities ≤ WM

Correction of mismatch problems

Occasionally data mismatches may occur in the INFOREX database (usually, if a phase composition is absent whereas the flag of its presence "1" is present in the initial CONDIT.EXP file master record). If such a mismatch is recognized by the *SPA* option, the warning message will alert you of the problem:

```
ATTENTION! A mismatch in INFOREX database encountered.

See 'base.err' file to correct the problem.
```

There are two ways to solve the problem:

- (1) contact Alexei Ariskin or Sergei Meshalkin at the address given above we will correct all of the necessary INFOREX files and send the corrected files to the user;
- (2) correct the database files independently, using the *Main Menu* options and information available in the BASE.ERR file (see *chapter 5.2*).

4.6. Exchange INFOREX Data (EID)

The *EID* option was implemented to facilitate data exchange between laboratories. We provide two functions in the data exchange option: "*Export*" and "*Import*". The first allows the user to save all the data from an experimental investigation to a separate file, whereas the second permits a complete dataset to be included in the INFOREX database. The exported file will appear in the D:\INFOREX\EXCHANGE\.. subdirectory and the file prepared for import must be in the same place and have ".IMP" extension.

So, if you need to conduct the *Exchange INFOREX Data* operations, press the [**F6**] key to choose the *EID* option from the *Main Menu* - a message will appear on the screen:

```
WARNING! Data import causes all previously created *.BAK files to be deleted. To save the files, leave the window now!

Press any key to continue...
```

This is a warning against entering the *.BAK files created during a previous "Import" operation. Press "ESC " or "Enter" and you will find yourself in the window as shown in Fig. 18.




Figure 18. Screen image of the Exchange INFOREX Data option <u>Export operations</u>

When you press the [F2] (EXPORT from INFOREX) key a window appears on the screen like that provided with the View option (Fig. 7). Here one can choose an experimental work to save to a file: do this same way as in the View of Experimental Data window, using the PgUp/PgDn keys or "Go to (Reference)" function. Press "Enter" to select a study and type a file name to be saved (up to 8 symbols) in the D:\INFOREX\EXCHANGE\... subdirectory.

Import operations

When you press the [F3] (IMPORT to INFOREX) key, a highlighted field will appear in the list of files to be imported, if there is at least one file having the ".IMP" extension (the number of such *.IMP files available in the D:\INFOREX\EXCHANGE\.. subdirectory is displayed - Fig. 18). Using the arrow or PgUp/PgDn keys move the highlighted field to a file to import and press "Enter". The EID control system will check out if all is well with your imported dataset and in a few seconds a message might appear on the screen:

O'K with your import data.
Press F5 to ACCEPT them
or ESC to CANSEL

In such a case press the [**F5**] key and your data will be properly inserted into the INFOREX database. If The *EID* control system encounters a problem in the initial datafile, a message will appear on the screen indicating what kind of errors should be corrected by the user in the file to be imported.

<u>Note:</u> any file prepared for import should have the same format as files created during *Export* operations.

The following fields must be present in the records of the input file: number of the reference (the last number in your databaseincreased by 1) in the first string, up to four strings no longer than 71 positions - for bibliography, then word "CONDIT" followed by master records in the CONDIT.EXP file format (see Fig. 3). If this study has data on starting compositions, type word "START" in the next string followed by compositional master records. The next strings represent nonstandard file data (if any) in the formats of VOLAT.EXP, MNLS.EXP, SULF.EXP, and REDOX.EXP files preceded by corresponding headings - VOLAT, MNLS, SULF, REDOX. Finally, strings of glasses and mineral compositions appear with their headings. The format of all these strings must be exactly the same as in the corresponding files of the database.

4.7. Perological Calculations (PC)

Although we did not design the INFOREX system to be used as a comprehensive compilation of techniques for the thermodynamic processing of experimental data, a set of additional programs has been included into the system for use in petrological studies. Most of the programs were written in FORTRAN and include the development of mineral-melt equilibria equations (geothermometers), calculations of Fe³⁺/Fe²⁺ ratio and water solubility in experimental glasses, different major-element reformulations as well as calculations of normalization parameters for projecting glass compositions onto ternary diagrams (Fig. 19).




Figure 19. Main functions of the *Petrological Calculations* option

To perform *Petrological Calculations*, press the [*F7*] key from the Main Menu and a submenu will appear on the screen:

```
Mineral-Melt Geothermometers
Fe3+/Fe2+ Ratios
Water Solubility
Special Tools
Ternary Diagrams
```

Select one of the items using $[\uparrow]$ and $[\downarrow]$ keys and press [Enter].

4.7.1. Mineral-Melt Geothermometers

The INFOREX database has been primarily designed to provide access to published mineral-melt compositional information with the goal of developing empirical geothermometers, which in turn can be applied to the modeling of natural igneous systems (*Ariskin et al., 1987-1997*). To facilitate this function, the "*Mineral-Melt Geothermometers*" (*MMG*) option was incorporated into INFOREX which makes it possible to create a system of equations describing low pressure mineral-melt equilibria for rock-forming minerals such as olivine, plagioclase, high-Ca (augites from diopside to sub-calcium species) and low-Ca (pigeonite and orthopyroxene) pyroxenes as well as spinels (*Ariskin and Nikolaev, 1996*).

Following the extraction of a set of experimental data using INFOREX, including tables of coexisting mineral and glass compositions (e.g., see Fig. 13), the calibration of the equations is accomplished by multiple regression in the form of a modified Arrhenius relation:

$$\ln K_i = A/T + B \lg f_{O2} + C \ln R_{L1} + D \ln R_{L2} + E \ln R_{L3} + F, \tag{1}$$

where K_i is an equilibrium constant of a formation reaction, distribution coefficient or exchange constant for component-i, T is the absolute temperature, f_{O2} is oxygen fugacity, R_{L1} , R_{L2} , and R_{L3} are the melt structure-chemical parameters, such as Si/O, NBO/T, Al/Si, (Na+K)/Al, Fe³⁺/Fe²⁺, etc., whereas the parameters A, B, C, D, E, and F are the regression constants.

To specify an equilibrium constant for the formation reaction (crystallization) of an end-member component, the degree to which K_i is compositionally dependent could be determined as a function of the melt component activities. The INFOREX MMG-function enables the user to apply single oxide, two-lattice (*Nielsen and Dungan, 1983*) or a modified two-lattice model (*Ariskin and Barmina, 1990*) for the melt component activities calculations.

The use of melt structure chemical parameters as independent variables in the Equations (1) is a means to fine-tune the model. <u>This approach assumes that the effects of nonideality on mineral-melt equilibria can be attributed primarily to the liquid phase</u>.

The multiple linear regression calculations are build into the *MMG*-function, allowing the user to numerically study different combinations of equilibrium constants

versus reciprocal temperature, oxygen fugacity or melt structure-chemical parameters. After the regression parameters are calculated, the *MMG*-function calculates the statistics on goodness of fit by comparing the calculated and experimental values. This permits the user to independently estimate the accuracy of the geothermometers. If they do not satisfy the requirements for modeling, the calculations may be repeated for another combination of variables, until a good fit with the experimental data is obtained.

Working with the MMG-function

The *MMG* window includes 7 countered boxes with function keys abbreviations given at the bottom (Fig. 20). These boxes can be switched by the [F3] (Next_Box) key, with the exception of the 5th box when you are working with the LIQUID PARAMETERS mode (in such a case one should use the [F5] key to go to the following REJECTING-mode). In other cases the [F5] (Browse) key serves to view or print out initial datafiles and output files including the calculated parameters of mineral-melt geothermometers.

The [F4] (Geotherm) key is not active at the initial stage of working with MMG-function: it can be used to browse the GEOTHERM.DAT file (see below) only after completing the whole cycle of operations including specification of output files.




Figure 20. Screen image of the **Mineral-Melt Geothermometers** function After completing operations for all 5 orthopyroxene end-members

Note: After completing respective calculations (i.e., the creation of an output file) the [F3] (Next_Box) key is transformed to the [F3] (New Calc.) key which serves to come back to the first box in order to start calculations for another mineral-melt pair (if

desired). If you are going to continue work with same mineral, press the [Shift+F3] combination to come back to the COMPONENT OR RATIO box!

OL 4 FeSi½02:a(OL)/a(LQ) MnSi½02:a(OL)/a(LO) $MgSi\frac{1}{2}O2:a(OL)/a(LQ)$ CaSi½02:a(OL)/a(LQ) PT. 4 FeAl2Si2O8:aPL/aLQ CaAl2Si2O8:aPL/aLO NaAlSi308 :aPL/aLQ KAlSi308 :aPL/aLQ AU 5 FeSiO3 :a(AU)/a(LQ) MnSiO3 :a(AU)/a(LQ)MgSiO3 : a(AU)/a(LQ)CaSiO3 : a(AU)/a(LQ)Al(AU)/Al(LQ) FeSiO3 :a(PG)/a(LQ) MnSiO3 : a(PG)/a(LQ)MgSiO3 : a(PG)/a(LQ)CaSiO3 : a(PG)/a(LQ)Al (PG) /Al (LQ) OP 5 FeSiO3 :a(OP)/a(LQ) MnSiO3 :a(OP)/a(LQ)MgSiO3 : a(OP)/a(LQ)CaSiO3 :a(OP)/a(LO) Al(OP)/Al(LQ) SP 5 FeCr2O4:a(SP)/a(LQ)Cr/Al :r(SP)/r(LQ) Mg/Fe2+:r(SP)/r(LQ) Fe3+/Fe2+(SP)Ti(SP)/Ti(LQ) MT 5 FeFe2O4:a(MT)/a(LQ)Fe3/Ti :r(MT)/r(LQ) Fe2TiO4:a(MT)/a(LQ) Fe3+/Fe2+(MT) Ti(MT)/Ti(LQ)

The left box in the *MMG* window always contains names and dates for files written in the "mineral-melt" format using the *SPA*-option (OP_LQ.DAT, PL_LQ.DAT, etc). One can select any file of interest by moving a highlighted field to the name desired using the arrow keys and pressing "Enter". In response to this action a FORTRAN program will create an auxiliary file called COEXIST.DAT (Table 1) which will be used in further thermodynamic processing.

A list of end-member components for the designated mineral will appear in the second *COMPONENT OR RATIO* box (see the complete list to left). After selecting an end-member component, Al distribution coefficient, or a partition coefficient ratio one should specify *MELT ACTIVITY MODEL*: at this stage of the work with INFOREX we prefer to use the two-lattice model (*Nielsen and Dungan, 1983*) which was successfully used in the development of COMAGMAT (*Ariskin et al., 1993*) and METEOMOD (*Ariskin et al., 1997b*) phase equilibria models.

The third and fourth boxes serve to specify a set of parameters which can be used as independent variables in the least square fitting (see Eq. 1). The *MAIN VARIABLES* box allows one to include logf_{O2} and total pressure (in bars) in the expressions being calibrated (reciprocal temperature 1/T,K is always present), whereas the *LIQUID PARAMETERS* box includes logarithms of 9 melt structure-chemical ratios (note, that one can add to the fitting simultaneously no more than 3 liquid parameters).

After specification of the independent variables, press the [**F5**] key and indicate if you wish to activate a special *REJECTING* procedure (YES/NO). If YES, a message will appear on the screen asking for the maximum temperature deviation (precision of

calculations, T_{exp} - T_{cal}) if one is to use the calculated mineral-melt equilibria expressions as geothermometers:

```
Max.Temp. Deviation: \pm 2.00~\sigma
```

For example, typing 2.0 means that all experimental points for which calculated temperatures more than 2 times exceed in absolute value an average temperature deviations T_{exp} - T_{cal} will be excluded from consideration at the next stage of thermodynamic processing.

Structure of output files

After calling for the REJECTING-procedure wait for a few seconds and results of the least square fitting will be displayed on the screen. Simultaneously an auxiliary file called REGRESS.OUT (Table 1) containing all of the independent variables will be created in the D:\INFOREX\SAVE directory:

```
Numb
          K3.dat
                   1/Temp,K
                                   Si/O P(bar)/T,K
                                                         lg f02
                              .3172E+00
                                         .6481E-03 -.7700E+01
       .2465E+01
                  .6481E-03
  1
                              .3219E+00
       .1906E+01 .6161E-03
                                         .6161E-03 -.1120E+02
                  .6246E-03
                                         .6246E-03
       .2035E+01
  3
                              .3200E+00
                                                    -.1143E+02
                                          .6357E-03 -.1173E+02
                   .6357E-03
                              .3161E+00
       .2195E+01
       .2197E+01
                  .6357E-03
                              .3183E+00
                                         .6357E-03
                                                    -.1173E+02
       .2479E+01
                  .6460E-03
                                         .6460E-03
                              .3151E+00
  6
                                                    -.1201E+02
       .2945E+01
                  .6698E-03
                              .3040E+00
                                          .6698E-03
                                                    -.1266E+02
                  .6789E-03
                                          .6789E-03 -.1290E+02
  8
       .3245E+01
                              .3123E+00
                              .3250E+00
  9
       .3480E+01 .6882E-03
                                         .6882E-03 -.1316E+02
                  .6930E-03
                                          .6930E-03
 10
       .3494E+01
                              .3062E+00
                                                    -.1332E+02
                                          .6930E-03 -.1332E+02
       .3729E+01
                   .6930E-03
                              .3258E+00
 11
       .2904E+01
                  .6698E-03
                              .3116E+00
                                         .6698E-03
                                                    -.1266E+02
                  .6789E-03
                                          .6789E-03
       .3137E+01
                              .3165E+00
                                                    -.1290E+02
 1.3
                                          .6930E-03 -.1332E+02
       .3403E+01
                   .6930E-03
                              .3242E+00
```

Figure 21. An example of REGRESS.OUT file used to develop *Opx*-melt expressions created from OP_LQ.DAT shown in Fig. 13

At this point, the program will ask if you wish to save the calculated parameters to GEOTHERM.DAT file (YES to save/"Enter" to skip) followed by a message asking to indicate the extension of output file to be created:

```
Enter file extension: TST
```

Examples of combined names of basic output files addressed to *Opx*-melt equilibrium are shown in Fig. 20 (note that digits "1","2",..."5" correspond to the number of an end-member component) .

```
2. Two Lattice Model by Nielsen, Dungan (1983)

Model fitting results for: MgSiO3

Independent variables Coefficient Stand. error

Constant -4.18669 .38848
```

```
245.97290
                                                    8211.36300
            1/Temp,K
    14 observations fitted

        Numb
        lnK-Comp
        l0**4/T
        lnK(exp)
        lnK(cal)
        D(lnK)
        T,C(exp)
        T,C(cal)

        1
        1.1353
        6.4810
        .9022
        .9020
        .0002
        1270.0
        1269.9

        2
        .8751
        6.1610
        .6450
        .6422
        .0028
        1350.1
        1349.2

        3
        .9418
        6.2460
        .7105
        .7108
        -.0003
        1328.0
        1328.1

        4
        1.0200
        6.3570
        .7862
        .7995
        -.0133
        1300.1
        1304.1

        5
        1.0195
        6.3570
        .7871
        .8009
        -.0138
        1300.1
        1304.2

        6
        1.1423
        6.4600
        .9079
        .8834
        .0244
        1275.0
        1267.9

        7
        1.3218
        6.6980
        1.0801
        1.0716
        .0085
        1220.0
        1217.7

        8
        1.4134
        6.7890
        1.1771
        1.1518
        .0254
        1200.0
        1193.3

        9
        1.4752
        6.8820
        1.2470
        1.2362
        .0108
        1180.1
        1177.3

                                                                                                                                  D(T,C)
                                                                                                                     1269.9
                                                                                                                      1349.2
                                                                                                                      1328.1
                                                                                                                                          -4.0
                                                                                                                      1304.2
                                                                                                                                          -4.2
                             6.6980 1.0801 1.0716
6.7890 1.1771 1.1518
6.8820 1.2470 1.2362
6.9300 1.2510 1.2635
                                                                                                                                            2.3
                                                                                                                                            6.7
                             6.8820 1.2470
6.9300 1.2510
6.9300 1.3161
6.6980 1.0661
            1.4752
1.4913
                                                                                    .0108
                                                                                                    1180.1
                                                                                                                      1177.3
                                                                                                                                            2.8
                                                                                -.0125
                                                                                                   1170.0
    10
                                                                                                                      1173.2
                                                               1.2761 .0400
1.0766 -.0105
                                                                                                                      1159.9
    11
             1.5438
                                                                                                   1170.0
                                                                                                                                         10.1
    12
             1.3028
                                                                                                   1220.0
                                                                                                                      1222.8
                                                                                                                                         -2.9
                                                               1.1545 -.0112
1.2751 -.0505
iations: .0160
                             6.7890 1.1433
6.9300 1.2247
                                                                                                                                         -3.0
    13
             1.3768
                                                                                                   1200.0
                                                                                                                     1202.9
             1.4533
                                                                                                   1170.0
                                                                                                                     1182.9
                                                                                                                                        -12.9
                                       Average Deviations:
                                                                                                                                           4.3
ABS[lnK(exp)-lnK(cal)] < STATISTICS > ABS[T(exp)-T(cal)]
                                                                                 ± 5.0° - 71.4%
±10.0° - 85.7%
          ± .05 - 92.9%
          ± .10 -100.0%
                                                                                ±10.0° - 85.7%

±15.0° -100.0%

±20.0° -100.0%

±25.0° -100.0%

±30.0° -100.0%
          ± .15 -100.0%
          ± .20 -100.0%
          ± .25 -100.0%
          ± .30 -100.0%
Rejected: 11 D(Temp) = 10.1
Rejected: 14 D(Temp) = 12.9
Model fitting results for: MgSiO3
                            -----
    Independent variables Coefficient Stand. error
 _____
                                                    -4.15677 .27784
8284.78900 201.46160
            Constant
                                                                                  201.46160
            1/Temp, K
                                                           .26996
          ln Si/O
    12 observations fitted
                                                                               D(lnK) T,C(exp) T,C(cal) D(T,C)
-.0004 1270.0 1270.1 -.1
.0035 1350.1 1349.0 1 1
Numb lnK-Comp 10**4/T lnK(exp) lnK(cal)
                             6.4810 .9022 .9026
           1.2122
                                                  .6450
                .9510
                                                                   .6415
                                                                                                                                            1.1
                             6.1610
6.2460
                                                                                                                   1328.0
                                                                                                  1328.0
            1.0181
                                                 .7105
                                                                  .7103
                                                                                    .0002
                             6.3570 .7862
6.3570 .7871
                                                               .7103 .0002
.7990 -.0128
.8008 -.0137
.8834 .0244
1.0709 .0092
1.1536 .0235
                                                                                                                                         -3.8
            1.0971
                                                                                                   1300.1
                                                                                                                     1303.9
                                               .7871
             1.0961
                                                                                                   1300.1
                                                                                                                      1304.2
                                                                                                                                         -4.1
                             6.4600 .9079 .8834
6.6980 1.0801 1.0709
            1.2196
                                                                                                   1275.0
                                                                                                                      1268.0
             1.4016
                                                                                                   1220.0
                                                                                                                      1217.5
                             6.7890 1.0601 1.0709 .0092
6.7890 1.1771 1.1536 .0235
6.8820 1.2470 1.2414 .0056
6.9300 1.2510 1.2651 -.0140
6.6980 1.0661 1.0776 -.0115
            1.4913
                                                                                                  1200.0 1193.8
              1.5504
                                                                                                   1180.1
                                                                                                                      1178.6
                                                                                                  1170.0 1173.5
1220.0 1223.1
           1.5705
                                                                                                                                          -3.5
             1.3809
                                                                                                                                         -3.1
          1.4538 6.7890 1.1433 1.1572
                                                                                -.0139
                                        Average Deviations:
                                                                                     .0111
ABS[lnK(exp)-lnK(cal)] < STATISTICS > ABS[T(exp)-T(cal)]
                                                                                ± 5.0° - 83.3%
±10.0° -100.0%
          ± .05 -100.0%
          ± .10 -100.0%
                                                                                ±15.0° -100.0%
          ± .15 -100.0%
          ± .20 -100.0%
                                                                                 ±20.0° -100.0%
                                                                                 ±25.0° -100.0%
          ± .25 -100.0%
                                                                                 ±30.0° -100.0%
          ± .30 -100.0%
```

Figure 22. An example of **OP_GTM_3.***** file developed with *MMG*-option two points were rejected after 1st stage of the thermodynamic processing

After completing the main *MMG* operations one can browse or print out the output files using the [*F5*] (*Browse*) key. The [*F4*] (*Geotherm*) key is now active allowing the user also to look at the GEOTHERM.DAT file which was changed if the

calculated regression parameters for a specified end-member component were written to the file (Fig. 23).

LAM=2 Regr B A/T,K					, ,	ln(LCP3)
PLAG: Kor=0.814 ======						
Fe -10.235 14202.680	.000	.000		-1.041(1)	.000(6)	.000(8)
An -5.515 10121.080	.000	.000	.000	.000(1)	, ,	.000(8)
Ab -10.684 17056.690	.000	.000	.000	.000(1)	.000(6)	.000(8)
or -7.906 10561.180	.000	.000	.000	.000(1)	.000(6)	.000(8)
OLIV: Kca=0.033						
Fa -4.343 5951.826	.000	.000	.000	206(1)	559(2)	.000(8)
Mn -6.161 6592.316	.000	.000		-1.265(1)	.000(2)	.000(8)
Fo -3.992 5623.459	.000	.000	.000	946(1)	798(2)	.000(8)
Ca -6.931 5047.081	.000	.000	.000	.000(1)	.000(2)	.000(8)
AUG : Kal=0.240						
Fs -2.654 8638.068	.000	.000	.000	3.464(1)	832(2)	.000(1)
Mn -5.367 11083.250	.000	.000	.000	2.528(1)	-2.279(2)	.000(1)
En -6.601 10247.640	.000	.000	.000	386(1)	.551(2)	.000(1)
Wo .054 589.070	.000	.000	.000	.000(1)	.572(2)	.000(1)
PIG : Kal=0.115						
Fs -2.011 4628.454	.000	.000	.000	1.313(1)	-1.219(2)	.000(1)
Mn584 5009.566	.000	.000	.000	2.301(1)	.000(2)	.000(1)
En -4.378 7820.873	.000	.000	.000	164(1)	.236(2)	.000(1)
Wo -15.910 15087.660	.000	.000	.000	-4.165(1)	2.674(2)	.000(1)
OPX: Kal=0.057						
Fs -3.296 4978.494	.000	.000	.000	.371(1)	-1.014(2)	.000(8)
Mn .202 2182.190	.000	.000	.000	2.153(1)	-3.442(2)	.000(8)
En -2.481 6256.046	.000	.000	.000	.743(1)	-1.002(2)	.000(8)
Wo -11.478 13457.620	.000	.000	.000	587(1)	.000(2)	.000(8)
SPIN: ???=0.000						
FeCr2-14.598 43656.360	0.457	.000	.000	-6.514(2)	.000(1)	0.000(5)
Cr/Al -7.478 22277.570	0.267	.000		-1.115(3)		0.000(5)
Mg/Fe-10.832 17996.040	0.581	.000		-3.203(2)	-1.192(8)	0.000(5)
FE3/2-30.623 58107.130	2.007	.000		-3.648(5)	-3.026(8)	0.000(5)
K(Ti)-16.919 29209.050	0.840	.000		-2.442(5)	, ,	0.000(5)
=======================================	=======	======	=======	========	========	=======

Figure 23. An example of **GEOTHERM.DAT** file developed using the *MMG*-option for calibration of a numeric model describing phase equilibria in meteoritic igneous systems (*Ariskin et al., 1997b*)

The GEOTHERM.DAT file is a basic datafile used to calculate phase equilibria in mafic igneous systems by means of programs addressed to the COMAGMAT series (*Ariskin et al., 1993*). This ability to calculate mineral-melt geothermometers for a given range of conditions or compositions is one of the main advantages of the INFOREX database, in that the user can independently calibrate phase equilibria equations, test them and use in calculations of melting-crystallization relationships in the systems of interest.

4.7.2. Fe³⁺/Fe²⁺ ratios

An important parameter for magma differentiation models is the Fe³⁺/Fe²⁺ ratio in the magma. Several equations for these calculations have been proposed, with the dependence of Sack et al. (1980) being one of the most popular among igneous petrologists for the last 15 years. We have collected the initial experimental information on the Fe³⁺-Fe²⁺ redox equilibrium which is available in INFOREX including over 440 experiments. This makes it possible to test the proposed equations for a given range of oxygen fugacities or melt compositions, or even to independently develop numerical constraints on the empirical dependencies (*Nikolaev et al., 1996*).

The Fe3+/Fe2+ ratios function has been designed to calculate Fe³⁺/Fe²⁺ ratios in experimental glasses at a given temperature, oxygen fugacity (lgf_{O2}), and melt composition. Calculations can be conducted using one of nine empirical equations given below. Also, this function allows one to test the precision of these equations using data on direct Fe³⁺/Fe²⁺ ratio determinations in the melts (see example of RE_LQ.DAT file in Fig. 17).




Figure 24. Screen image of the **Fe3+/Fe2+ RATIOS** calculation function After work with OP_LQ.DAT and RE_LQ.DAT files

Working with the Fe3+/Fe2+ ratios function

The main Fe3+/Fe2+ ratios function window is configured in a way similar to the Mineral-Melt Geothermometers window and includes 5 boxes with the names of data files at the left side, and the names of output files at the right side (Fig. 24). Input files for the calculations are **_LQ.DAT, where double asterisk denotes RE, OP, LQ, OL, etc. Output files receive names such as **_OXI_*.OUT, where double asterisk denotes the same abbreviation as in the data file, whereas the character replacing single asterisk ranges from 1 to 9 depending on the calculation model selected.

There is no switching between the countered boxes: going from one box to another is done automatically after selection of a datafile or choice of an equation for the Fe³⁺/Fe²⁺ calculations. The [*F5*] (*Browse*) key serves to view or print out initial datafiles or output files including results of the calculations. The [*F4*] (*New Calc.*) key is not active at initial stage of working with the function. It can be used for repeating or conducting another set of calculations after specification of names of output files.

Basic	Equat	ion:	Sack e	t al.,	1980						
Numb	INFO	REX	P,kb		emp,C	lg	f02		Fe3	+/Fe2+	
	=Work=	=Run=		==EXP=	CAL	====EXP=	CAL=	$=\delta(i)=$	EXP=	CAL=	 δ(i)
1	77	2	.0	1270	ND	-7.70	ND	ND	ND	.109	ND
2	166	1	.0	1350	ND	-11.20	ND	ND	ND	.009	ND
3	166	2	.0	1328	ND	-11.43	ND	ND	ND	.009	ND
4	166	5	.0	1300	ND	-11.73	ND	ND	ND	.009	ND
5	166	6	.0	1300	ND	-11.73	ND	ND	ND	.009	ND
6	166	7	.0	1275	ND	-12.01	ND	ND	ND	.010	ND
7	166	8	.0	1220	ND	-12.66	ND	ND	ND	.010	ND
8	166	9	.0	1200	ND	-12.90	ND	ND	ND	.011	ND
9	166	10	.0	1180	ND	-13.16	ND	ND	ND	.012	ND
10	166	11	.0	1170	ND	-13.32	ND	ND	ND	.011	ND
11	166	12	.0	1170	ND	-13.32	ND	ND	ND	.013	ND
12	166	17	.0	1220	ND	-12.66	ND	ND	ND	.010	ND
13	166	18	.0	1200	ND	-12.90	ND	ND	ND	.011	ND
14	166	20	.0	1170	ND	-13.32	ND	ND	ND	.012	ND

Figure 25. An example of OP_OXI_1.OUT file

The Fe3+/Fe2+ ratios function provides an access to nine basic equations proposed since 1980 (note, that Nikolaev et al. proposed 3 different equations for the melts/magmas which could be addressed to tholeiitic, subalkaline, and alkaline series):



Examples of output files which are created in the D:\INFOREX\SAVE directory are given in Figs 25-27: the column **Numb** includes successive numbers of the calculations; two columns addressed to the **INFOREX** section contain reference numbers (**Work**) and corresponding run numbers (**Run**); run pressure is given in the column **P,kb**, whereas experimental and calculated temperatures are in **EXP** and **CAL** columns of the **Temp,C** section.

Basic Equation: Kilinc et al., 1983

Numb	INFO	REX	P,kb	T	emp,C	lg	f02		Fe3+,	/Fe2+	
	Work=	=Run=		==EXP=	CAL=	==EXP=	CAL=	$=\delta(i)=$	====EXP===	CAL-	= δ(i)
1	152	1	.0	1461	1405.7	-7.17	-6.69	48	.095	.075	.020
2	152	3	.0	1459	1337.6	-8.21	-7.11	-1.10	.078	.045	.033
3	152	4	.0	1459	1357.4	-8.21	-7.30	91	.070	.044	.026
4	152	6	.0	1459	1295.3	-9.21	-7.69	-1.52	.059	.027	.032
5	152	11	.0	1459	1368.9	-6.15	-5.35	80	.194	.130	.064
6	152	12	.0	1459	1425.4	-6.15	-5.86	29	.156	.135	.021
7	152	20	.0	1362	1325.4	-7.93	-7.58	35	.102	.085	.017
8	152	35	.0	1366	1301.0	-8.91	-8.28	63	.062	.045	.017
9	152	38	.0	1366	1354.2	-8.91	-8.80	11	.049	.046	.003
10	152	50	.0	1363	1298.4	-6.81	-6.18	63	.189	.137	.052
11	152	51	.0	1363	1300.5	-6.81	-6.20	61	.185	.136	.049
$\sigma = \sqrt{(\Sigma)}$	δ(i))	² /n)			84.2			.77			.035
δ(av)=	$\delta(av) = \Sigma \delta(i) / n$				73.3			.68			.030

Figure 26. An example of **RE_OXI_1.OUT** file created from RE_LQ.DAT file given in Fig. 17

One can find experimental and calculated oxygen fugacities in the columns **EXP** and **CAL** of the **IgfO2** section (if the user intends to use the equations as oxybarometers). Differences between experimental and calculated values are placed in the field $\delta(i)$ =EXP-CAL. Experimental and calculated FeO_{1.5}/FeO ratios are given in the columns **EXP** and **CAL** followed by the $\delta(i)$ = EXP-CAL value. If a parameter is unknown (there were no data in the reference source) the ND (*No Data*) index will appear in the corresponding field.

General statistics on the comparisons between experimental and calculated values for temperature, oxygen fugacity, and Fe^{3+}/Fe^{2+} ratio is always given in the two lower rows: $\sigma = (\Sigma(\delta(i))2/n)1/2$ and $\delta(av) = \Sigma|\delta(i)|/n$ (av - average).

Basic Equation: Nikolaev et al., 1996 (Tholeiitic)

Numb	INFOREX	P,kb	Temp,C	lgf02	Fe3+/Fe2+
	Work==Run=		EXP—CAL—	==EXP===CAL=	$=\delta$ (i) $===EXP==CAL==\delta$ (i)

									<u>INF</u> (<u>OREX</u>	- 4.2	USER
											MA	NUAL
1	152	1	.0	1461	1518.0	-7.17	-7.60	.43	.095	.109	014	
2	152	3	.0	1459	1461.6	-8.21	-8.23	.02	.078	.079	001	
3	152	4	.0	1459	1497.6	-8.21	-8.51	.30	.070	.077	007	
4	152	6	.0	1459	1451.5	-9.21	-9.15	06	.059	.058	.001	
5	152	11	.0	1459	1380.6	-6.15	-5.51	64	.194	.158	.036	
6	152	12	.0	1459	1461.9	-6.15	-6.17	.02	.156	.157	001	
7	152	20	.0	1362	1422.2	-7.93	-8.44	.51	.102	.120	018	
8	152	35	.0	1366	1427.5	-8.91	-9.43	.52	.062	.073	011	
9	152	38	.0	1366	1533.4	-8.91	-10.24	1.33	.049	.075	026	
10	152	50	.0	1363	1309.4	-6.81	-6.32	49	.189	.162	.027	
11	152	51	.0	1363	1302.9	-6.81	-6.26	55	.185	.155	.030	
$\sigma = \sqrt{(\Sigma)}$	(δ(i))	² /n)			69.3			.57			.020	
δ(av)=	-Σ δ(i	.) /n			53.6			.44			.016	

Figure 27. An example of RE_OXI_7.OUT file

If we examine the data of Fig. 26, one can see that the model of Nikolaev et al. (1996)

works better than that of Kilinc et al. (1983)

<u>Note:</u> As follows from the Figs 26 and 27 the Fe3+/Fe2+ ratios function can not only be used to test existing models calculating Fe^{3+}/Fe^{2+} ratios in the melt, but also to calculate $\lg f_{02}$ values for melts with known Fe^{3+}/Fe^{2+} ratios. To perform this operation one should create first an ASCII file like RE_LQ.DAT (see the chapter 4.5).

4.7.3. Water Solubility

Another useful application of the INFOREX system is its ability to calculate water solubility in magmas. In the last decade, many investigations have been undertaken to develop thermodynamic models of water solubility. Each, to different extents, do so by accounting for the changes in the partial molar volume of H₂O and for the relations between H₂O species in silicate melts, see the review by *Holloway and Blank (1994)*. One problem that we encounter in applying this information to natural systems is the fact that little experimental data on water solubility exists for natural systems, particularly for basalts and andesites. Therefore, a classic thermodynamic approach based on the use of independent constraints, such as the effect of temperature on the partial molar volume of H₂O, may generate additional errors to calculations of water solubility. Towards the end of producing useful, if empirical constraints, we have extracted the water solubility experiments from the INFOREX database to develop a purely empirical equation which can be applied to a wide range of systems and conditions (Almeev and Ariskin, 1996).

Processing of the data was carried out by multiple linear regression calculations with the melt structure-chemical parameters used as independent variables. A set of

17 parameter combinations were investigated. The best internal precision was obtained using the expression:

$$ln(C_{H2O})_{L} = [-14710.2 + 38438.3(Si/O)_{L}](1/T) + 0.5868lnP -$$

$$21.455(Si/O)_{L} + 3.894(Al/Si)_{L} + 4.393,$$
(2)

where $(C_{H2O})_L$ is water content in silicate liquid (wt.%), P is total pressure (bars), and Si/O and Al/Si are atomic ratios in the melt. Comparison of the experimental values with those calculated using the Equation (2) indicates an average accuracy of 0.3 wt.% H_2O ; solving the inverse problem for pressure (if water content in a saturated melt is given) also indicates an average accuracy of 250 bars (Almeev and Ariskin, 1996). This is consistent with our current understanding of the strong effect of pressure on water solubility.

Derivation of the equation for the calculation of water solubility in magmatic melts shows how the information in the INFOREX database can be developed further. As mentioned above, there is a large existing dataset of water saturated experiments, but direct analytical information on the water contents in the experimental glasses is absent for almost all of them. To accomplish the goal we have implemented the *Water Solubility* function providing an access to 3 basic equations proposed and improved since 1979, including the above mentioned model of Almeev and Ariskin (1996):

- 1. Burnham's model (1979-1994)
- 2. Stolper's model (1982-1995)
- 3. Almeev & Ariskin model (1996)

The first function is based on the most recent Burnham's model parameters (*Burnham*, 1994), whereas the second one uses parameters of water solubility equation given by *Dixon et al.* (1995) with an algorithm provided by *Holloway and Blank* (1994). Working with the *Water Solubility* function

The Water Solubility function window includes 5 boxes (Fig. 28). In the ENTER DATAFILE box one should select an initial datafile, containing information on the experimental conditions as well as quenched glasses compositions. These files are

implied to have been created earlier using the *SRC*- and *SPA*-options: they have names like **_LQ.DAT, where double asterisk denotes LQ, HB, OL, etc. The *CHOOSE A MODEL* box includes an ability to select a *Water Solubility* model (note the pressure restrictions given below at the box bottom - see Fig. 28). References for each model appear automatically in the *FULL REFERENCE* box accompanied by some additional comments. There is no switching between the countered boxes: going from one box to another is done automatically after pressing "Enter" to select a datafile or choose an equation for the water solubility calculations.




Figure 28. Screen image of the **Water Solubility** function After work with LQ_LQ.DAT and HB_LQ.DAT files

The [F5] (Browse) key serves to display or print out initial datafiles or output files including results of the calculations. The [F4] (New Calc.) key is not active during the initial stage of working with the function: it can be used for repeating or conducting another set of calculations after specification of names of output files.

Note: In contrast to other equations the Stolper's model requires few parameters you can change in the STOLPER'S MODEL PARAMETERS box - the default values for these parameters are taken from Holloway and Blank (1994).

Structure of output files

Output files receive names such as **_H2O_*.OUT, where double asterisk denotes the same abbreviation as in the data file, whereas the character replacing single asterisk ranges from 1 to 3 depending on the calculation model selected (Fig.

28). Examples of output files which are created in the D:\INFOREX\SAVE directory are given in Figs 29-30: the column **Numb** includes successive numbers of the calculations; two columns addressed to the **INFOREX** section contain reference numbers (**Work**) and corresponding run numbers (**Run**); **FLUID** column includes flags of the presence (*Yes*) or absence (*No*) of a vapor phase; run pressure and experimental temperature are given in the column **P,kb** and **Temp,C**, respectively. Experimental (if present) and calculated water contents are given the next 3 columns addressed to the **H2O solubility** title (Fig. 29).

Model for H2O solubility calculations: Almeev & Ariskin (1996)

Numb	INFO	REX	FLUID	P,kb	Temp,C	H2	20 solubi	lity
	Work=	=Run==				EXP	CAL	δ (i)
1	111	36	Yes	2.0	800	5.95	6.05	_
2	111	66	Yes	2.0	1000	5.96	6.87	_
3	112	11	Yes	3.0	850	ND	7.98	ND
4	112	12	Yes	3.0	875	ND	7.61	ND
5	112	27	Yes	3.0	850	ND	8.16	ND
6	112	28	Yes	3.0	875	ND	7.91	ND
7	112	46	Yes	3.0	850	ND	7.92	ND
8	112	47	Yes	3.0	875	ND	7.64	ND
9	112	79	Yes	3.0	850	ND	7.95	ND
10	133	3	Yes	5.0	750	ND	11.59	ND
11	133	28	Yes	4.9	725	ND	12.29	ND
12	133	42	Yes	4.9	725	ND	12.18	ND
13	229	12	Yes	2.0	940	ND	6.35	ND

Figure 29. An example of **HB_H2O_2.OUT** file Obtained from the HB LQ.DAT file given in Fig. 16

In some cases water contents in experimental glasses were measured directly (in opposite to the water contents for the first two glasses of Work 111 which represent calculated values - Fig. 29). For such experiments, the *Water Solubility* function allows one to estimate the accuracy of the water contents calculations using 3 models considered above. To do this you should first choose the *H2O Content* - function of the *SRC*-option (see Fig. 11), followed by the *Only melt compositions* + *Melts Equilibrated* with *All Phases* functions of the *SPA*-option. These operations will result in the creation of a LQ LQ.DAT file, including only glass compositions with measured H₂O contents.

In such cases, the *Water Solubility* function calculates the difference between experimental and computed water contents placed in the field $\delta(i)$. General statistics on the comparisons between experimental and calculated values is always present in the two lower rows: $\sigma=(\Sigma(\delta(i))2/n)1/2$ and $\delta(av)=\Sigma|\delta(i)|/n$ (av - average) - see Fig. 30.

Numb	INFO	O solubi	llity ===δ(i)====					
1	work = 142	=Run== 7	Yes	1.0	1100	4.49	3.36	1.13
2	142	8	Yes	2.0	1100	6.04		
3	142	9	Yes	3.0	1100	7.40		
4	142	10	Yes	4.0	1100	8.70		
5	142	11	Yes	5.3				
6	201	20	Yes	1.0	1170	4.40		
7	201	22	Yes	1.0	1100	4.50		1.12
8	201	24	Yes	1.0	1050			
9	201	40	Yes	3.0	1030	6.50		.57
10	201	45	Yes	3.0	1000	6.80		.99
11	201	72	Yes	5.0	1100	8.80	8.52	.28
12	201	79	Yes	5.0	1000	8.60		
13	201	89	Yes	3.0	1200	5.65		
14	201	90	Yes	3.0	1200	5.62		
15	201	91	Yes	3.0	1100	6.53	6.18	.35
16	201	92	Yes	3.0	1100	6.27	6.18	.09
			rated rur	ns with	directly	determin	ed H2O	16
	(δ(i))							.96
δ(av)	$=\Sigma \mid \delta$ (i) / N						.89

Model for H2O solubility calculations: Stolper (1982-1994)

Figure 30. An example of **LQ_H2O_2.OUT** file

Obtained from a LQ_LQ.DAT file including data on water solubility experiments in andesitic systems at 1 to 5 kbar pressures

4.7.4. Special Tools

This is an additional option which permits one to perform some simple recalculations and formatting of the INFOREX data files created using the *SPA*-option. These 3 functions (Fig. 31) include:

1. Recast to 100% including H₂O contents as follows:

 Σ =SiO₂+TiO₂..+Cr₂O₃+H₂O, SiO₂'=100*(SiO₂/ Σ), TiO₂'=100*(TiO₂/ Σ), etc. Name of output file: **RECAST_1.DAT** (Fig. 32)

2. Recast to 100% excluding H₂O (dry residue) as follows:

 Σ =SiO₂+TiO₂..+Cr₂O₃, SiO₂'=100*(SiO₂/ Σ), TiO₂'=100*(TiO₂/ Σ), etc. Name of output file: **RECAST_2.DAT** (Fig. 33)

3. Creation of data files which permit the user to test mineral-melt geothermometers used to develop phase equilibria models, such as COMAGMAT (*Ariskin et al., 1993*) and METEOMOD (*Ariskin et al., 1997b*). This option replaces an INFOREX database number of an experiment in the melt composition row with the

value of experimental temperature. Such a file can be used without further reformatting in a special *Thermometry* mode of COMAGMAT and METEOMOD programs to compare calculated and experimental temperatures for each experimental point.

Name of output file: **_THMEQ.DAT (Fig. 34)




Figure 31. Screen image of the **Special Tools** function After work with LQ LQ.DAT and OP LQ.DAT files

Si02	Ti02 Al203	FeO	MnO	MgO	CaO	Na20	K20	P205	Cr203	Н2О	N	n
56.04	1.10 17.51	6.05	0.10	3.25	6.43	4.17	0.79	0.25	0.00	4.31	142	7
55.22	1.09 17.25	5.97	0.09	3.20	6.33	4.11	0.78	0.25	0.00	5.71	142	8
54.52	1.07 17.03	5.89	0.09	3.16	6.25	4.06	0.77	0.24	0.00	6.91	142	9
53.86	1.06 16.83	5.82	0.09	3.13	6.18	4.01	0.76	0.24	0.00	8.02	142	10
53.19	1.05 16.62	5.75	0.09	3.09	6.10	3.96	0.75	0.24	0.00	9.18	142	11
56.80	0.49 18.58	4.34	0.08	2.48	7.31	3.96	1.53	0.17	0.00	4.25	201	20
56.75	0.49 18.56	4.34	0.08	2.48	7.30	3.96	1.53	0.17	0.00	4.35	201	22
60.56	0.53 18.10	3.49	0.11	2.77	6.86	3.65	1.66	0.00	0.02	2.25	201	24
55.67	0.48 18.21	4.25	0.08	2.43	7.16	3.88	1.50	0.17	0.00	6.16	201	40
55.52	0.48 18.16	4.24	0.08	2.43	7.14	3.87	1.49	0.17	0.00	6.42	201	45
54.49	0.47 17.82	4.16	0.07	2.38	7.01	3.80	1.46	0.17	0.00	8.16	201	72
54.59	0.47 17.85	4.17	0.07	2.39	7.02	3.81	1.47	0.17	0.00	7.99	201	79
56.13	0.49 18.36	4.29	0.08	2.45	7.22	3.92	1.51	0.17	0.00	5.40	201	89
56.14	0.49 18.36	4.29	0.08	2.46	7.22	3.92	1.51	0.17	0.00	5.37	201	90
55.66	0.48 18.20	4.25	0.08	2.43	7.16	3.88	1.50	0.17	0.00	6.18	201	91
55.80	0.48 18.25	4.26	0.08	2.44	7.18	3.89	1.50	0.17	0.00	5.95	201	92

Figure 32. An example of RECAST_1.DAT file

Obtained from the LQ LQ.DAT used to demonstrate calculations in Figs 28 and 30

SiO2	TiO2 Al2O3	FeO	MnO	MgO	CaO	Na20	K20	P205	Cr203	Н2О	N	n
58.56	1.15 18.30	6.33	0.10	3.40	6.72	4.36	0.82	0.26	0.00	0.00	142	7
58.56	1.15 18.30	6.33	0.10	3.40	6.72	4.36	0.82	0.26	0.00	0.00	142	8
58.56	1.15 18.30	6.33	0.10	3.40	6.72	4.36	0.82	0.26	0.00	0.00	142	9
58.56	1.15 18.30	6.33	0.10	3.40	6.72	4.36	0.82	0.26	0.00	0.00	142	10
58.56	1.15 18.30	6.33	0.10	3.40	6.72	4.36	0.82	0.26	0.00	0.00	142	11
59.33	0.51 19.40	4.53	0.08	2.59	7.63	4.14	1.59	0.18	0.00	0.00	201	20

4.14

4.14

1.59

1.59

0.18

0.18

0.00

0.00

0.00 201 91

0.00 201 92

Figure 33. An example of **RECAST_2.DAT** file Obtained from the LQ_LQ.DAT used to demonstrate calculations in Figs. 28 and 30

7.63

7.63

2.59

2.59

Both RECAST_1.DAT and RECAST_2.DAT can be further used in graphics plotting programs, whereas contents of major elements in the **_THMEQ.DAT still needed to be recast to 100%.

Figure 34. An example of **OP_THMEQ.DAT** file Obtained from the OP LQ.DAT shown in Fig. 13

<u>Note:</u> In contrast to the Fe3+/Fe2+ ratios and Water Solubility functions, the [F4] key in the Special Tools function window is labeled as New Exec: in fact this is same as New Calc.

4.7.5. Ternary Diagrams

59.33 0.51 19.40

59.33 0.51 19.40

4.53

4.53

0.08

0.08

One of traditional techniques of the use of melting-experiment information is to project experimental glass compositions onto ternary diagrams, such as AFM (alkalis-

FeO-MgO), *OLIV-CPX-QTZ*, or *AB-AN-QTZ*. Despite the fact that such a projection results in the lost of a part of potentially important chemical information, many scientists still use that approach for making comparisons of natural glasses (rocks, inclusions) compositions with those obtained in experiments. To make this work easier and faster we have developed the *Ternary Diagrams* function which allows one to calculate the projection parameters based on the techniques of Grove (*Tormey et al., 1987; Grove, 1993*) and *Walker et al. (1979*).




Figure 35. Screen image of the **Ternary Diagrams** function After work with HB LQ.DAT file

Working with the *Ternary Diagrams* function

The *Ternary Diagrams* function window includes 5 boxes (Fig. 35), similar to the *Special Tools* function window (Fig. 31). In the *ENTER DATAFILE* box one can select an initial datafile, whereas output files are displayed in the right *OUTPUT FILES* box. There is no switching between the countered boxes: going from one box to another is done automatically after pressing "Enter" to select a datafile or choose a type of ternary diagram for which the projection parameters will be calculated.

Structure of output files

Output files receive names such as **_TRN**.OUT, where the first double asterisk denotes the same abbreviation as in the data file, whereas the second ranges from 01 to 10 depending on the ternary diagram selected (Fig. 35). Examples of output files which are created in the D:\INFOREX\SAVE directory are given in Figs. 36-38.

For the files named **_TRN01.OUT to **_TRN09.OUT the column Numb includes successive numbers of the calculations, two columns addressed to the INFOREX section contain reference numbers (Work) and corresponding run numbers (Run), Temp,C column contains the value of experimental temperature; the next 3 columns include the calculated projection parameters followed with norms in oxygen units.

1 01#	,CPx	‡ , Qt	z# -	projection	parame				_	-+ -1	1007	
	INFOR Work		TempC	01#	CPx#	Qtz#	QTZ	•	om Tormey OLIV CPX		,1987 ORT	AP
1	111 111	36	800 1000	.068	011 .083	.944 .828	25.4 14.8	54.0 72.9	1.83 1.6 1.5		18.5 6.1	.0
•	112	11	850	.074	213	1.139	49.2	48.5	3.2 -9.2		5.4	.2
-	112 112	12 27	875 850	.081 .108	196 319	1.115 1.212	43.5 47.3	52.8 57.3	3.1 -7.6 4.2-12.5		4.8	.3
•	112 112	28	875 850	.108	312 219	1.204	36.1	66.3 49.2	3.2 -9.3 3.7 -9.8		.5	. 4
	112	46 47	875	.082 .081	219 203	1.137	50.9 42.3	56.7	3.7 -9.8		2.9	.2
_	112 133	79 3	850 750	.080	245 164	1.165 1.083	41.3	53.4	2.9 -8.7 3.1 -6.4		8.5 9.9	.3
11	133	28	725	.087	210	1.123	48.8	43.2	3.8 -9.1	3	13.1	.0
	133 229	42 12	725 940	.061 .042	163 .188	1.102 .771	47.8 9.9	42.5 73.2	2.6 -7.1 .5 2.4		13.0 11.6	.0

Figure 36. An example of **HB_TRN01.OUT** file Obtained from the HB_LQ.DAT shown in Fig. 16

2 01	#,Pl#	, Qt:	z# -	projection	paramet			ternary from Walk	er et al	.,1979===
Numb	INFO		Temp(01#	Pl#	Qtz#	QTZ	PLAG	OLIV	DIOPS
1	111 111	36 66	800	.024	.408	.568	56.9 38.2	40.9 55.1	2.4	2 1.4
3	112	11	850	.039	.211	.750	78.6	22.1	4.1	-4.8
-	112 112	12 27	875 850	.043	.245 .228	.713 .725	74.1 77.4	25.4 24.3	4.4 5.0	-4.0 -6.7
-	112 112	28 46	875 850	.045	.311 .199	.644 .758	67.9 79.6	32.8 20.9	4.8 4.4	-5.5 -5.0
•	112 112	47 79	875 850	.041	.254	.704 .695	73.2 72.8	26.5 28.1	4.3	-4.0 -4.8
10	133	3	750	.029	.259	.712	73.9	26.9	3.0	-3.7
12	133 133 229	28 42 12	725 725 940	.031 .027 .037	.217 .221 .676	.752 .753 .287	78.9 78.2 28.0	22.8 22.9 66.0	3.2 2.8 3.6	-4.9 -3.9 2.5

Figure 37. An example of **HB_TRN07.OUT** file Obtained from the HB LQ.DAT shown in Fig. 16

First 4 columns in the file **_TRN10.OUT are the same as in other output files, whereas next 9 columns contain 3 projection parameters for the *AFM* diagram followed by 6 additional routinely used petrochemical coefficients and ratios.

³ A#, F#, M# - projection parameters for AFM ternary diagram ALK - Na2O+K2O; MG# - Mg/(Mg+Fe*); Ca# - Ca/(Ca+Al); AL# - (Na+K)/Al

Numb	INFO	REX '	TempC	A#	F#	M#	SiO2	FeO/MgO	ALK	MG#	CA#	AL#
	Work	Run					 wt% 	wt	 wt%=	 mol=	=mol=	<u></u> mol
1	111	36	800	.825	.132	.042	67.04	3.14	8.23	.362	.079	.819
2	111	66	1000	.594	.306	.099	53.54	3.09	4.85	.366	.245	.385
3	112	11	850	.649	.308	.044	74.37	7.04	3.71	.202	.155	.343
4	112	12	875	.614	.300	.086	71.85	3.48	3.99	.339	.176	.358
5	112	27	850	.670	.289	.041	74.13	7.00	3.89	.203	.113	.371
6	112	28	875	.706	.253	.041	70.16	6.14	4.81	.225	.147	.419
7	112	46	850	.573	.348	.079	74.29	4.40	3.11	.289	.165	.297
8	112	47	875	.627	.295	.077	71.63	3.82	3.97	.318	.177	.364
9	112	79	850	.758	.224	.018	72.68	12.67	5.15	.123	.134	.445
10	133	3	750	.000	.000	.000	74.40	.00	5.50	.000	.126	.502
11	133	28	725	.000	.000	.000	75.80	.00	4.80	.000	.123	.414
12	133	42	725	.818	.175	.007	76.00	26.75	5.00	.062	.133	.449
13	229	12	940	.709	.262	.029	60.29	8.97	7.04	.166	.221	.487

Figure 38. An example of **HB_TRN10.OUT** file Obtained from the HB_LQ.DAT shown in Fig. 16

Chapter 5. Troubleshooting

5.1. Memory problems

The most common problem you can encounter working with the INFOREX database is memory problems. These problems look like INFOREX "does not want to work" with the SCR-, SPA-, or Petrological Calculations options. This is due to the fact that the INFOREX database needs at least 590K conventional memory (see the chapter 1.2). However, often the users have too many resident programs stored in the memory to run the system properly. To avoid this problem, we recommend that you make sure there is enough conventional memory on your computer before starting INFOREX.

To do this type in the DOS command line **mem /c /p** and a message on the screen will provide you a complete information on the distribution of conventional memory between resident programs and free memory available.

All is well if there is more than 590K free conventional memory. If not, one should make some arrangements in the computer configuration.

1. If you are working under **DOS 5.0** or higher, just use the **MEMMAKER.EXE** program in ..\DOS directory to improve your computer memory allocation (follow all of the instructions appeared on the screen). Use the command **mem /c /p** to make sure you have accomplished the goal.

If there is still not enough memory after the *memmaker.exe operations*, one should make changes in your **AUTOEXEC.BAT** file to free some additional memory.

Contact a technician to do that (in such a case it is convenient to have 2 copies of AUTOEXEC.BAT files with one of them addressed to the use of the INFOREX database).

2. If you are working with DOS emulated from **WINDOWS 95** or higher, it is also necessary to make changes in the **AUTOEXEC.BAT** and **CONFIG.SYS** files. A simple example of these files is given below (the main idea is to load as many programs into the upper memory as possible). If you use another programs you can load them like given in these examples.

AUTOEXEC.BAT:

```
@echo off
PATH=C:\WINDOWS;C:\WINDOWS\COMMAND
SET TMP=D:\TMP
mode con codepage prepare=((866) C:\WINDOWS\COMMAND\ega3.cpi)
mode con codepage select=866
```

CONFIG.SYS:

```
DEVICE=C:\WINDOWS\HIMEM.SYS
DEVICEHIGH=C:\WINDOWS\EMM386.EXE RAM
DEVICEHIGH=C:\WINDOWS\COMMAND\display.sys con=(ega,,1)
INSTALLHIGH=C:\WINDOWS\COMMAND\NLSFUNC.EXE C:\WINDOWS\COMMAND\COUNTRY.SYS
Country=007,866,C:\WINDOWS\COMMAND\country.sys
DOS=HIGH
DOS=UMB
BUFFERS=15,0
FILES=30
```

3. There are should not be any memory problems if you are working under **OS/2 Merlin**, in that it provides 610K free memory automatically without a special configuration.

5.2. Correction of mismatch problems

In the case of mismatch problems encountered during *SPA* operations (working with the *Mineral-Melt Equilibria* or *Two mineral equilibria* functions), a message on the creation of a **BASE.ERR** file appears on the screen. This auxiliary file is located in the D:\INFOREX\SAVE\... directory and includes a set of strings addressed to the INFOREX numbers of experiments where these errors were detected. These strings are formatted

in the CONDIT.EXP fashion (see Fig. 3) with the detected errors being usually caused by mismatches in the phase assemblages flags. For example, if index OL1 is given there while the olivine composition for same experiment is absent in the OLIV.EXP file.

To correct the problem the user must write the INFOREX database numbers of experiments where these errors were detected and run INFOREX again. Using the *Browse* function of the *GSI* option one should examine the mineral composition file which was involved in the *Mineral-Melt Equilibria* or *Two mineral equilibria* processing in detail.

After determining what compositions are absent, the user must return to the *Main Menu*, select the *V/E/I* option and using the [*F2*] (**Go to**) key to call for the study considered. Press the [*F4*] (*Edit*) function key to initialize the *Edit* function, and correct the corresponding strings of run conditions by deleting the mistaken symbols <1>. Press the [*F3*] (*Next*) key to call for a function to save the changes (see *the chapter 4.3.2*). After that, one should repeat the search process beginning from the *SRC*-option followed by SPA-functions of interest. These operations need to be done until the BASE.ERR file is empty.

Acknowledgements

We are very grateful to Chris Chalokwu and Roger Nielsen for critically reading and correcting the manuscript. Financial support of Russian Foundation of Basic Research (grants NN 94-05-16098, 96-05-64231, and 96-07-89054) and International Science Foundation (grant N M1N000) is gratefully acknowledged.

6. References

- Al'meev R.R. and Ariskin A.A. (1996) Mineral-melt equilibria in a hydrous basaltic system: computer modeling. Geochem. Intern., 34 (7): 563-573.
- Ariskin A.A. (1999) Phase equilibria modeling in igneous petrology: use of COMAGMAT model for simulating fractionation of ferro-basaltic magmas and the genesis of high-alumina basalt. J. Volcanol. Geotherm. Res., 90: 115-162.
- Ariskin A.A., Barmina G.S., and Frenkel M.Ya. (1987) Computer simulation of basalt magma crystallization at a fixed oxygen fugacity. Geochem. Intern., 24 (6): 85-98.
- Ariskin A.A. and Barmina G.S. (1990) Equilibria thermometry between plagioclases and basalt or andesite magmas. Geochem. Intern., 27 (10): 129-134.

- Ariskin A.A., Bouadze K.V., Meshalkin S.S., and Tsekhonya T.I. (1992) INFOREX: A database on experimental studies of phase relations in silicate systems. Amer. Miner., 77: 668-669.
- Ariskin A.A., Frenkel M.Ya., Barmina G.S., and Nielsen R.L. (1993) COMAGMAT: a Fortran program to model magma differentiation processes. Computers and Geosciences, 19: 1155-1170.
- Ariskin A.A. and Nielsen R.L. (1993) Application of computer simulation of magmatic processes to the teaching of petrology. J. Geol. Education, 41: 438-441.
- Ariskin A.A., Barmina G.S., Meshalkin S.S., Nikolaev G.S., and Almeev R.R. (1996) INFOREX-3.0: A database on experimental phase equilibria in igneous rocks and synthetic systems. II. Data description and petrological applications. Computers and Geosciences, 22: 1073-1082.
- Ariskin A.A. and Nikolaev G.S. (1996) An empirical model for the calculation of spinel-melt equilibrium in mafic igneous systems at atmospheric pressure: I. Chromian spinels. Contrib. Mineral. and Petrol., 123: 282-292.
- Ariskin A.A., Meshalkin S.S., Almeev R.R., Barmina G.S., and Nikolaev G.S. (1997a) INFOREX information retrieval system: analysis and processing of experimental data on phase equilibria in igneous rocks. Petrology, 5 (1): 28-36.
- Ariskin A.A., Petaev M.I., Borisov A.A., and Barmina (1997b) METEOMOD: a numerical model for the calculation of melting-crystallization relationships in meteoritic igneous systems. Meteoritics and Planetary Sciences, 32: 123-133.
- Borisov A.A. and Shapkin A.I. (1990) A new empirical equation rating Fe3+/Fe2+ in magmas to their composition, oxygen fugacity, and temperature. Geochem. Intern. 27(1): 111-116.
- Burnham C.W. (1994) Development of the Burnham model for prediction of the H₂O solubility in magmas. Reviews in Mineralogy (*Volatiles in Magmas*), Eds. Carroll M.R. and Holloway J.R. Mineralogical Society of America, 30: 123-129.
- Dixon J.E., Stolper E.M., and Holloway J.R. (1995) An experimental study of water and carbon dioxide solubilities in mid-ocean ridge basaltic liquids. Part 1: Calibration and solubility models. J. Petrol., 36: 1607-1631.
- Ghiorso M.S. and Sack R.O. (1995) Chemical mass transfer in magmatic processes IV. A revised and internally consistent thermodynamic model for the interpolation and extrapolation of liquid-solid equilibria in magmatic systems at elevated temperatures and pressures. Contrib. Mineral. Petrol., 119: 197-212.
- Grove T.L. (1993) Corrections to expressions for calculating mineral-components in "Origin of calc-alkaline series lavas at Medicine Lake volcano by fractionation, assimilation and mixing" and "Experimental petrology of normal MORB near the Kane Fracture Zone: 22°-25°N, mid-Atlantic ridge". Contrib. Mineral. Petrol., 114: 422-424.
- Holloway J.R. and, Blank J.G. (1994) Application of experimental results to C-O-H species in natural melts. Reviews in Mineralogy (*Volatiles in Magmas*), Eds. Carroll M.R. and Holloway J.R. Mineralogical Society of America, 30: 187-230.
- Longhi J. (1991) Comparative liquidus equilibria of hypersthene-normative basalts at low pressure. Amer. Mineral., 76: 785-800.
- Kilinc A., Carmichael I.S.E., Rivers M., and Sack R.O. (1983) The Ferric-ferrous ratio of natural silicate liquids equilibrated in air. Contrib. Mineral. Petrol., 83: 136-140.

- *Kress V.C. and Carmichael I.S.E.* (1988) Stoichiometry of the iron oxidation reaction in silicate melt. Amer. Miner., 73: 1267-1274.
- Kress V.C. and Carmichael I.S.E. (1991) The compressibility of silicate liquids containing Fe2O3 and the effect of composition, temperature, oxygen fugacity and pressure on their redox states. Contrib. Mineral. Petrol., 108: 82-92.
- Meshalkin S.S., Ariskin A.A., Barmina G.S., Nikolaev G.S., and Almeev R.R. (1996) Development of a database on experimental crystal-melt equilibria for igneous rocks: INFOREX system, version 3.0. Geochem. Intern., 34 (2): 85-90.
- Meshalkin S.S. and, Ariskin A.A. (1996) INFOREX-3.0: A database on experimental phase equilibria in igneous rocks and synthetic systems. I. Datafile and management system structure. Computers and Geosciences, 22: 1061-1071.
- Mysen B.O. (1991) Relations between structure, redox equilibria of iron, and properties of magmatic liquids. Advances in Physical Geochemistry (*Physical Chemistry of Magmas*). New-York (Springer-Verlag), 9: 41-98.
- *Nielsen R.L. and Dungan M.A.* (1983) Low-pressure mineral-melt equilibria in natural anhydrous mafic systems. Contrib. Mineral. Petrol., 84: 310-326.
- Nielsen R.L. (1990). The theory and application of a model of open magma system processes. Reviews in Mineralogy (*Modern Methods of Petrology*), Eds. Nicholls J. and Russel J.K. Mineralogical Society of America, 24: 65-106.
- Nikolaev G.S., Borisov A.A., and Ariskin A.A. (1996) Calculation of the ferric-ferrous ratio in magmatic melts: testing and additional calibration of empirical equations for various magmatic series. Geochem. Intern., 34 (8): 641-649.
- Sack R.O., Carmichael I.S.E., Rivers M., and Ghiorso M.S. (1980) Ferric-ferrous equilibria in natural silicate liquids at 1 Bar. Contrib. Mineral. Petrol., 75: 369-376.
- Tormey D.R., Grove T.L. and Bryan W.B. (1987) Experimental petrology of normal MORB near the Kane Fracture Zone: 22°-25°N, mid-Atlantic ridge. Contrib. Mineral. Petrol., 96: 121-139.
- Walker D., Shibata T., and DeLong S.E. (1979) Abyssal tholeiites from the Oceanographer Fracture zone: phase equilibria and mixing. Contrib. Mineral. Petrol., 70: 111-125.
- Weaver J.S. and Langmuir C.H. (1990) Calculation of phase equilibrium in mineral-melt systems. Computers and Geosciences, 16: 1-19.
- Yang H.-J., Kinzler R.J., and Grove T.L. (1996) Experiments and models of anhydrous basaltic olivine-plagioclase-augite saturated melts from 0.001 to 10 kbar. Contrib. Mineral. Petrol., 124: 1-18.

Appendix 1

On-line help (HELP.TXT file)

On-Line Help Facility

INFOREX EXPERIMENTAL DATABASE (Ver. 4.0, 1997)
Copyright (C) 1991-97. All rights reserved.
A.A.Ariskin and S.S.Meshalkin (Moscow, Russia)

INFOREX is a sophisticated manager for igneous petrologists and geochemists interested in the use of experimental information on phase equilibria in molten silicate systems. It is designed as an aid in searching the published experimental data, and extracting specific subsets of that data.

Its major function is as a computerized melting-experiment

2

3

reference manual that presents data in a fixed format, with key word search functions and flags as an integral part of the program. In addition, the program allows one to derive numerical mineral-melt equilibria expressions (geothermometers) based on any data set.

INFOREX is hot-key driven, allowing users to easily and quickly sort through the data using any set of experimental conditions or rock type parameters. There are $7\ \text{main}$ INFOREX options including:

Esc - exit PgDn 26 25

= On-Line Help Facility =

<f2> <f3> <f4> <f5> <f6></f6></f5></f4></f3></f2>	System Configuration General System Information View/Edit/Input data files Selecting Run Conditions Selecting Phase Assemblage Exchange INFOREX Data Petrological Calculations Quit	[SC] [GSI] [V/E/I] [SRC] [SPA] [EID] [PC]

- SC $\,$ includes an utility to change INFOREX settings, such as PATHS, PHASES, SYSTEMS, CONTAINERS, etc.
- GSI allows the user to receive General System Information on the current state of the INFOREX database as well as to browse basic data files on the screen.
- $\begin{tabular}{lll} $V/E/I$ provides operations for VIEWing, EDITing, and INPUT (updating) \\ & of the database files. It includes also an ability to write \\ & master records to a file or print them out. \\ \end{tabular}$
- SRC permits one to Specify Run Conditions, such as rock types plus a range of pressures, temperatures, fO2, and melt compositions.

Esc - exit - PgUp/PgDn - 26 25

On-Line Help Facility =

- SPA allows the user to Select Phase Assemblages and to create phase equilibria compositional datafiles corresponding to the previously selected set of experiments (see the SRC-option).

 These files can be used in further thermodynamic processing.
- EID provides Export/Import operations including an ability to write results of new experimental studies directly to INFOREX database or save information under a given reference to a separate file.
- PC includes 5 kinds of Petrological Calculations, such as:
 - 1. Creation of Mineral-Melt equilibria geothermometers;
 - 2. Computation Fe3+/Fe2+ ratios for the experimental glasses;
 - Estimates of Water Solubility for the experimental glasses produced in water-saturated experiments;
 - Special Tools, including creation of files for testing of the COMAGMAT and METEOMOD phase equilibria models;
 - Computations of projection parameters for Ternary Diagrams routinely used in petrology.

QUIT - same as ESC, back 1 level.

26 25

On-Line Help Facility =

4

GENERAL SYSTEM INFORMATION

If you used the <F2> function key in the MAIN MENU you should wait for the program to read database information.

The available data base parameters will appear on the screen and further actions in GSI routine can be initiated using function keys keys or cursor:

Using the Function Keys in GSI-Routine

- F1: Help displays all specific options of the INFOREX system and always is the same as in MAIN MENU.
- F2: Notes involves 5 pages with tables deciphering indices, search words and flags used in the database.
- F3: Save saves a copy of the GSI screen in the 'general.inf' file on '\INFOREX\SAVE\' disk directory.

F4: Print - prints the current text screen on the printer. F5: Browse - permits browsing of any INFOREX database files. - same as ESC, quits the current operation. Esc - exit⊨ PqUp/PqDn 26 25 On-Line Help Facility = To browse the INFOREX system database files press <F5> key and select "Phase compositions"/"Run conditions"/"References" highlighted field using the cursor and the "Enter" key. If you selected the "Phase compositions" mode you must select a desired compositional file using the cursor keys located on the righthand keypad. Move the highlighted field to file name and and press "Enter". VIEW/EDIT/INPUT DATA FILES The V/E/I-routine screen is always divided on two windows with a reference in the upper part. One can scan the 'bibl.txt' reference file here using "PageDn" or "PageUp" keys. If you need to quickly find a reference based on its number, use the <F2> function key. To change a mode of the routine, press <F3> or <F4> function keys. Esc - exit d PgUp/PgDn 🗗 26 25 On-Line Help Facility = 6 Using the Function Keys in VIEW/EDIT-modes - displays all specific options of the INFOREX system and F1: Help always is the same as in MAIN MENU. F2: Goto - allows you to leave or come back to a reference by means of scanning the reference list. F3: Input - calls INPUT mode. - calls EDIT mode. F4: Edit writes to a file on '\INFOREX\SAVE\' disk directory or prints on printer if type 'PRN' instead of a file name. F5: Write - same as ESC, quits the current operation. F9: Onit To browse or edit information concerning a given experimental study in V/E modes, $\,$ press "Enter" and wait for a few seconds until 'condit__.exp' file data appear in the lower window. To scan this file use "PageDn" or "PageUp" keys (see next page). To make changes in the selected study information press <F4> and move a highlighted field to the field of interest with arrow keys. Make changes and press "Enter". desc - exit⊨ d PgUp/PgDn 📙

26 25

On-Line Help Facility =

Note, that after pressing "Enter" at a selected reference the current meaning of the <F2> key will be transformed to "Notes" and the <F3> key will be transformed to "Next".

To set to a new file representing the starting or phase compositions press the $\ensuremath{<F3>}$ [INPUT] function key.

Using the Function Keys in INPUT-mode

F1: Help $\,$ - displays all specific options of the INFOREX system and

F2: Notes - involves 5 pages with tables deciphering indices, search words and flags used in the database.


F3: Next - saves the data input and goes to a new dataset.

F4: Edit - permits marking and copying blocks of the strings selected.

(see User's Manual).

F5: Write - writes to a file on '\INFOREX\SAVE\' disk directory or prints on printer if type 'PRN' instead of a file name.

```
F9: Quit - same as ESC, quits the current operation.
   desc - exit⊨
                                                                         = PgUp/PgDn =
26 25
                            On-Line Help Facility =
                                                                                    8
      \mathtt{NOTE:} <\mathtt{Shift-F3}\mathtt{>} - enables you to come back to the initial screen
      from any stage of execution of the V/E/I-mode procedures.
                            SELECTING RUN CONDITIONS
      The SRC-option is implemented to perform search-match operations through the database according to some "filters".
                   Using the Function Keys in SRC-mode
      F1: Help
                   - same as in other MAIN MENU options.
                  - same as in other MAIN MENU options.
      F2: Notes
      F3: Next Box- go to another active box (system type, name of igneous
                     rock, combination of volatiles, containers, etc.).
                   - allows one to specify a range of REDOX-conditions.
      F5: Search - starts the SEARCH through the database.
                   - same as ESC, quits the current operation.
      F9: Ouit
   Esc - exit
                                                                         🗕 PgUp/PgDn 📙
26 25
                          On-Line Help Facility =
                                                                                   9
                         SELECTING PHASE ASSEMBLAGE
        Using the Function Keys in SPA-mode (1st step, DATA FILTERING)
      F1: Help
                  - same as in other MAIN MENU options.
                   - same as in other MAIN MENU options.
      F2: Notes
      F3: Window - go to another active box (window). F4: On/Off - enables/disables flag "\sqrt{}" in front of a string addressed
                     to Phase Assemblages or References.
      F5: Exclude - deletes strings marked to disable. F9: Phases - Go to 2nd step.
        Using the Function Keys in SPA-mode (2nd step, CHOOSING PHASES)
      F1,F2:
                   - same as in other MAIN MENU options.
      F3-F8:
                   - Don't use the keys! (see User's Manual for information).
      F9: Quit
                   - quits the operation and returns to Main Menu.
      NOTE: Move the highlighted field to a phase field and press "ENTER".
   Esc - exit
                                                                        PqUp/PqDn
```



11

F7: Save - saves current settings to a file.
F9: Quit - quits the operation and returns to Main Menu.

PgUp/PgDn

EXCHANGE INFOREX DATA

On-Line Help Facility =

Here you can:

- EXPORT all data linked with a reference to a separate file just press <F2> key, then select a reference and type a file name to be created. The file exported is now ready to import, if you are going to send it to another user of INFOREX.
- IMPORT a data set directly to the INFOREX database press <F3> key, select a file to be imported (".imp" extension!) and press "Enter".

NOTE: All strings in the imported files must be formatted exactly as described in the User's Manual: you will find there commentaries on all error messages.

One can UNDO the import operations IMMEDIATELY after a reference was imported - just press <F5> key as prompted at the bottom of EID-option screen (the undo procedure doesn't work if you have leaved the EID-option after import operations).

Esc - exit - PgUp/PgDn 26 25

= On-Line Help Facility =

12

PETROLOGICAL CALCULATIONS

We have included in INFOREX few special functions which could be useful in petrological studies. Those can be called from MAIN MENU by pressing the <F7> key followed by selection of one of 5 function:

- 1. MINERAL-MELT GEOTHERMOMETERS: development of expressions describing mineral-melt equilibria in terms of temperature, pressure, fO2, and liquid (glass) composition.
- 2. Fe3+/Fe2+ RATIO CALCULATION: calculations of ferric-ferrous species in silicate melts as a function of temperature, fO2, and composition.
- 3. WATER SOLUBILITY: calculations of water contents in experimental glasses obtained in water-saturated conditions.
- 4. SPECIAL TOOLS: recast of experimental glass compositions to 100% with an ability to create files for testing COMAGMAT and METEOMOD models.
- 5. TERNARY DIAGRAMS: calculations of projection parameters for plotting experimental glass compositions onto ternary diagrams.

Esc - exit PgUp/PgDn

26 25

On-Line Help Facility =

13

Petrological Calculations: MINERAL-MELT GEOTHERMOMETERS

This option allows one to create a system of equations describing 1 atm pressure mineral-melt equilibria for OL, PL, AUG, PIG, OPX and SPINels. Following the extraction of a set of experimental data using INFOREX (including tables of coexisting mineral and glass compositions - F5 key) the calibration of the equations is accomplished by multiple regression in the form of a modified Arrhenius relation:

ln Ki = A/T + BlgfO2 + C ln(Fe3+/Fe2+)L + DlnRL + E,

where Ki is an equilibrium constant of a mineral-forming reaction, distribution coefficient or exchange constant for component-i, T is the absolute temperature, fO2 is oxygen fugacity, RL is a melt structure-chemical parameter, such as Si/O, NBO/T, Al/Si, (Na+K)/Al, etc., while the parameters A, B, C, D and E are the regression constants.

After the regression parameters are calculated, the option calculates

the statistics on goodness of fit, comparing the calculated and experi-Esc - exit 🚽 PgUp/PgDn 📙 26 25 On-Line Help Facility = 14 Petrological Calculations: Fe3+/Fe2+ RATIO CALCULATION 1.Sack R.O., Carmichael I.S.E., Rivers M., Ghiorso M.S. (1980) Ferric-ferrous equilibria in natural silicate liquids at 1 Bar. CMP 75: 369-376. 2.Kilinc A., Carmichael I.S.E., Rivers M., Sack R.O. (1983) The Ferric-ferrous ratio of natural silicate liquids equilibrated in air. CMP 83: 136-140. 3.Kress V.C. & Carmichael I.S.E.(1988) Stoichiometry of the iron oxidation reaction in silicate melt. Amer. Miner. 73: 1267-1274. 4.Borisov A.A.&Shapkin A.I.(1990) A new empirical equation rating Fe3+/Fe2+ in magmas to their composition, oxygen fugacity, and temperature. Geochem. Intern. $27\,(1)$; 111-116. 5.Mysen B.O. (1991) Relations between structure, redox equilibria of iron, and properties of magmatic liquids. Physical Chemistry of Magmas. Advances in Physical Geochemistry. V. 9. New-York (Springer-Verlag): 41-98. 6.Kress V.C.&Carmichael I.S.E.(1991)The compressibility of silicate liquids containing Fe2O3 and the effect of composition, temperature, oxygen fugacity and pressure on their redox states. CMP 108: 82-92. 7.Nikolaev G.S., Borisov A.A., Ariskin A.A. (1996) Calculation of the Fe3+/ Fe2+ ratio in magmatic melts: testing and additional calibration of empirical equations for various petrochemical series. Geochem. Intern. 34(8): 641-649. desc - exit⊨ 26 25 On-Line Help Facility = 15 Petrological Calculations: WATER SOLUBILITY This option allows one to calculate water content in experimental glasses saturated with H2O and test experimental H2O contents at given conditions. MODELS INVOLVED IN THE CONSIDERATION: Burnham C.W. (1994) Development of the Burnham model for prediction of the H2O solubility in magmas. Volatiles in Magmas. Reviews in Mineralogy. V. 30: 123-129. IMPORTANT: ONLY for basalts and andesites, not for granites and rhyolites! Dixon J.E., Stolper E.M., Holloway J.R. An experimental study of water and carbon dioxide solubilities in mid-ocean ridge basaltic liquids. J. Petrol. in press. ALGORITHM: Holloway & Blank (1994) Reviews in Miner. V. 30: 187-230. Almeev R.R., Ariskin A.A. (1996) Mineral-melt equilibria in a hydrous basaltic system: computer modeling. Geochem. Intern. 34(7): 563-573. EMPIRIC:lnC(H2O) = 4.39+(48438*Si/O-14710)/T+0.59*lnP-21.45*Si/O+3.89*Al/Si =|Esc - exit|=

26 25

On-Line Help Facility =

16

Petrological Calculations: SPECIAL TOOLS

er cent including M20 as follow

- 1 Recast to 100 per cent including H2O as follows: $\Sigma = Sio2+..+H2O$, SiO2'=100*(SiO2/ Σ) etc Name of output file: RECAST_1.DAT
- 2 Recast to 100 per cent excluding H2O (dry residue) as follows: $\Sigma{=}SiO2{+}..Cr2O3,~SiO2{\,}'{=}100{\,}'(SiO2/\Sigma)$ etc Name of output file: RECAST 2.DAT
- 3 Create file for testing the mineral-melt equilibrium temperatures to be calculated with COMAGMAT program (Ariskin et al.,1993):
 Will be created from experimental glass compositions, including 12 major


```
element contents with the value of experimental TEMP in the end of row.
     Name of output file: ** THERM.DAT
   Esc - exit
                                                                           🗕 PgUp/PgDn 📙
26 25
                           = On-Line Help Facility =
                                                                                    17
                 Petrological Calculations: TERNARY DIAGRAMS
     The purpose of the function is to calculate parameters for projecting
   experimental glass compositions onto ternary diagrams routinely used in
   petrology. To use this function one should first create datafiles, like
   OL_LQ.DAT, PL_LQ.DAT, etc. by means of SPA-option.
   The list of ternary diagrams included is followed:
   1. AB-AN-QTZ
                    (Tormey et al., 1987) 6. OLIV-DI-PLAG (Tormey et al., 1987)
                     (Tormey et al., 1987)
                                              7. OLIV-CPX-QTZ (Walker et al.,1979)
   2. AB-AN-CPX
  3. PLAG-CPX-QTZ (Tormey et al.,1987) 8. OLIV-PLAG-CPX (Walker et al.,1979) 4. OLIV-PLAG-QTZ (Tormey et al.,1987) 9. OLIV-PLAG-QTZ (Walker et al.,1979)
  5. OLIV-DI-QTZ (Tormey et al.,1987) 10. AFM, FeO/MgO, ALK, MG#, CA#, AL#
      Details of the calculations one can find in papers given in the window
  called "FULL REFERENCE". The output files are saved in D:\INFOREX\SAVE\
   directory. These files contain parameters necessary to put experimental
  points onto ternary diagrams using any standard graphics package.
      The name of an output file includes first two letters to be the same as
  in initial datafile, followed by TRN plus the number of diagram in the list above (e.g., OL_TRN_3 if OL_\overline{\text{LQ}}.DAT was used an initial datafile).
  desc - exit⊨
26 25
                           On-Line Help Facility =
                                                                                    18
                              ADDITIONAL COMMENTS
  desc - exit⊨
                                                                                PgUp
```

INFOREX abbreviations (NOTES.TXT file)

```
21 38
                     IGNEOUS ROCK AND SYSTEM INDICATORS
   N (No volatiles)
                          V (with Volatiles)
                                                / F (Fluid saturated)
   - (undozed) / + (dozed with minor or major elements, crystals, etc.)
  TERRA()/LUNAR(L)/METEORitic(M)/PLANETary(P)-Martian, Venusian/SYNTHetic(S)
                                                     29 Eucrite
                    -> PER
                             15 Basalt
                                             -> BAS
                                                                     -> EUC
       Peridotite
                            16 AlkalineBAS -> ALB
                    -> PXT
                                                     30 Howardite
                                                                     -> HOW
     2 Pyroxenite
                            17 High-Al_BAS
18 Andesit.BAS
                    -> LHR
                                             -> HAB
                                                     31 Diogenite
                                                                     -> DGN
     3 Lherzolite
                    -> HAR
                                            -> ANB
                                                     32 Shergottite
     4 Harzburgite
                                                                     -> SNC
                    -> KTM
                            19 Andesite
                                             -> AND
                                                     33 WIncAllende
     5 Kimberlite
                                                                     -> WTA
                                             -> DAC
     6 Picrite
                    -> PTC
                             20 Dacite
                                                     34 Chondrite
                                                                     -> CHN
                                                     35 ANT series
     7 Komatiite
                    -> KOM
                            21 Rhyolite
                                             -> RHY
                                                                     -> ANT
                    -> MAR
                                            -> SYE
     8 MiAmRu-rock
                            22 Syenite
                                                     36 Green glass
                                                                    -> GGL
     9 Amphibolite -> AMP
                            23 Trachyte
                                            -> TRA
                                                    37 Pyrolite
                                                                     -> PYR
    10 Eclogite
                    -> ECL
                            24 Diorite
                                             -> DIO 38 Synthetic
                                                                     -> SYN
```

```
11 Melilitite
                      -> MEL
                              25 Granodior.
                                               -> GRD
                                                      39 Aug Minette
                                                                       -> MNT
    12 Nephelinite
                    -> NEP
                              26 Granite
                                               -> GRN
                                                       40 Tonalite
                                                                        -> TON
                     -> NOR
                              27 Gneiss
                                               -> GNE
                                                       41 Charnockite
                                                                       -> CHA
     13 Norite
                              28 Sedimentary
    14 Boninite
                      -> BON
                                               -> SED
                                                       42
   Esc - exit
                                                                          PgDn
7 66
      NAMES OF PHASE COMPOSITION FILES AND THEIR INDECIES IN CONDIT.EXP
                                                                               2
                         File
                                 Index
                                                              File
                                                                      Index
              Phase
                                            Ν
                                                   Phase
        Ν
                         LIO
                                             19
                                                 AlkFeldspar
                                                              ALFS
         1
            Glass (melt.)
                                  LO
                                                                       ΑF
         2
                                             20
            Plagioclase
                         PLAG
                                  PT.
                                                              APAT
                                                                       AΡ
                                                Apatite
            Olivine
                         OLIV
                                             21
                                                              ARM
         3
                                  OT.
                                                Armalcolite
                                                                       AR
         4
            Augite+Di
                         AUG
                                  ΑU
                                             22
                                                 Carbonate
                                                              CARB
                                                                       CB
                                             23
         5
            Pigeonite
                         PIG
                                  PG
                                                Epidote
                                                              EPID
                                                                       EΡ
            Clino-Px(*)
                                  CP
         6
                         CPX
                                             2.4
                                                Alumina-sil
                                                              ALST
                                                                       AS
            Ortho-Px
                         OPX
                                  OP
                                             25
                                                 Quartz
                                                              QUAR
                                                                       QU
         8
            ProthoPx(*)
                         PRPX
                                  PR
                                             26
                                                Perovskite
                                                              PERV
                                                                       ΡV
         9
            {\tt Ilm-Hem\_sol}
                         ILM
                                  IL
                                             27
                                                Ps-brookite
                                                              PSBR
                                                                       PB
        10
            MagnUlv_sol
                         MAGN
                                  ΜТ
                                             28
                                                Sulfide
                                                              SULF
                                                                       SF
        11
            Spinel
                          SPIN
                                  SP
                                             29
                                                Ca-perovsk
                                                              CAPV
                                                                       РС
        12
            Garnet
                         GARN
                                  GR
                                             30
                                                Mg-perovsk
                                                              MGPV
                                                                       {\tt Pm}
        13
            Nepheline
                         NEPH
                                  ΝP
                                             31
                                                Tourmaline
                                                              TOUR
                                                                       TU
        14
            Leucite
                          LEUC
                                  LC
                                             32
                                                 Corundum
                                                              CRND
                                                                       CO
        15
            Melilite
                         MEL
                                  ML
                                             33
                                                 Zircon
                                                              ZRCN
                                                                       ZR
        16
           Ortho-Amph
                         AMPH
                                  AM
                                             34
                                                Monticellit
                                                              MONT
                                                                       MO
        17
           Clino-Amph
                         HORN
                                  ΗВ
                                             35
                                                Rutile
                                                              RUT
                                                                       RU
           Mica (*)
                         MICA
                                                              SPHE
                                                                       SH
        18
                                  MΙ
                                             36
                                                Sphene
           see additional comments
                                                 see additional comments
   Esc
       - exit
                                                                     17 45
                = INDECIES OF CONTAINERS AND OXYGEN BUFFERS =
                                                                               3
                                Index
                                                                      Index
           Containers
                                               Oxygen buffers
       1 Pt loops
                                 PTL
                                             Iron-wustite
                                                                       ΙW
       2 Graphite_caps.le
                                 GRA
                                             Wustite-magnetite
                                                                       WM
       3 Fe_(\pmPt_sat)c.psules
                                 FEC
                                                                       TM
                                             Iron-magnetite
       4 Mo_capsule
                       psule
                                 MOC
                                            Magnetite-hematite
                                                                       MH
       5 Pt-Ag_capsule psule
                                 РТА
                                                                       QFM
                                             Ouartz-Favalite-Magn.
                                             Iron-quartz-fayalite
       6 Pt_capsule
                      psule
                                 PTC.
                                                                       IOF
       7 Ag-Pd_capsule psule
                                 AGP
                                            Nickel
                                                      -NiO
                                                                       NNO
      8 Gold_capsule psule
                                 AUC
                                             Cobalt
                                                       -CoO
                                                                       CCO
                                                       -MnO
       9 Iridium_loops
                                 IRL
                                            Mn304
                                                                       MMO
     10 Peridotite cappsule
                                 PER
                                             Graphite-CO-CO2-O2 *
                                                                       COC
                                             Graphite-CO-CO2-O2 **
      11 Al203 crucibles
                                 COR
                                                                       COC
                      psule
     12 Rh capsule
                                 REC
                                             Graphite-COH system
                                                                       GRS
     13 Pt-Rh capsule psule
                                 PRH
                                             Graphite-methane (ND)
                                                                       GCH
     14 Au-Pd_capsule_psule
                                 AUP
                                           * >5 kbar, ** <5 kbar(see below)
   Esc - exit⊨
                                                                      d PgUp/PgDn d
```

Buffer	lgf02 = a0	+ a1/TK + a2*	Deference	
Buller	a0	a1	a2	References
1 IW 2 WM 3 IM 4 MH 5 QFM 6 IQF 7 NNO 8 CCO 9 MMO 10 COC>5 11 COC<5 12 GRA	6.471 16.092 8.990 13.480 8.290 6.396 9.360 7.936 13.380 2.740 -0.044 4.620	-26834.7 -36951.3 -29260.0 -23847.6 -24441.9 -27517.5 -24930.0 -25070.0 -25680.0 -19559.0 -20586.0 -22324.0	0.055 0.083 0.061 0.019 0.092 0.050 0.046 0.055 * 0.081 0.130 -0.028 0.189 #	Myers, Eugster, 1983 Myers, Eugster, 1983 Huebner J.S., 1971 Myers, Eugster, 1983 Myers, Eugster, 1983 Myers, Eugster, 1983 Huebner, Sato, 1970 Myers, Gunter, 1979 Huebner, Sato, 1970 Woermann et al., 1977 French B.H., 1966 Ulmer, Luth, 1991



Appendix 2

List of INFOREX-4.0 (1997) experimental studies

- Stolper E. Experimental petrology of eucritic meteorites // Geochim. Cosmochim. Acta, 1977, V. 41, P. 587-611.
- 2 Grove T.L., Vaniman D.T. Experimental petrology of very low-Ti (VLT) basalts // In: "Mare Crisium: the view from Luna-24". Pergamon Press: 1978, P. 445-471.
- 3 Walker D., Kirkpatrick R.J., Longhi J., Hays J.F. Crystalization history of lunar picritic basalt sample 12002: phase equilibria and cooling-rate studies // Geol.Soc.Amer.Bull., 1976, V. 87, N 5, P. 646-656.
- 4 **Biggar G.M., O'Hara M.J., Peckett A., Humphries D.J.** Lunar lavas and the achondrites:petrogenesis of protohypersthene basalts in the maria lava lakes // In: Proc. Sec. Lunar Sci. Conf. The Mit Press: **1971**, V. 1, P. 617-643.

- 5 Donaldson C.H., Usselman T.M., Williams R.J., Lofgren G.E. Experimental modeling of the cooling history of Apollo 12 olivine basalts // In: Proc. Lunar Sci. Conf. 6th. Pergamon Press: 1975, V. 1, P. 843-869.
- 6 Walker D., Longhi J., Lasaga A.C., Stolper E.M., Grove T.L., Hays J.F. Slowly cooled mickrogabbros 15555 and 15065 // In: Proc. Lunar Sci. Conf. 8th. Pergamon Press: 1977, V. 2, P. 1521-1547.
- **Grove T.L.,Bence A.E.** Experimental study of pyroxene-liquid interaction in quartz-normative basalt 15597 // In: Proc. Lunar. Sci. Conf. 8th. Pergamon Press: **1977**, V. 2, P. 1549-1579.
- **Akella J., Boyd F.R.** Partitioning of Ti and Al between coexisting silicates, oxides and liquids // In: Proc. Lunar Sci. Conf. 4th. Pergamon Press: **1973**, V. 1, P. 1049-1959.
- *Huebner J.S., Lipin B.R., Wiggins L.B.* Partitioning of chromium between silicate crystals and melts // In: Proc. Lunar Sci. Conf. 7th. Pergamon Press: **1976**, V. 2, P. 1195-1220.
- *Grove T.L., Raudsepp M.* Effect of kinetics on the crystallization of quartz normative basalt 15597: an experimental study // In: Proc.Lunar Planet. Sci. Conf. 9th. Pergamon Press: **1978**, V. 1, P. 585-599.
- **Hess P.C., Rutherford M.J., Campbell H.W.** Ilmenite crystallization in nonmare basalt: genesis of KREEP and high-Ti mare basalt // In: Proc. Lunar Planet. Sci. Conf. 9th. Pergamon Press: **1978**, V. 1, P. 705-724.
- **Rhodes J.M., Lofgren G.E., Smith D.P.** One atmosphere melting experiments on ilmenite basalt 12008 // In: Proc. Lunar Planet. Sci. Conf. 10th. Pergamon Press: **1979**, V. 1, P. 407-422.
- 13 Grove T.L., Beaty D.W. Classification, experimental petrology and possible volcanic histories of the Apollo 11 high-k basalts // In: Proc.Lunar Planet.Sci.Conf.11th. Pergamon Press: 1980, V.1, P.149-177.
- **O'Hara M.J., Biggar G.M., Hill P.G., Jefferies B., Humphries D.J.** Plagioclase saturation in lunar high titanium basalt // Earth and Planet. Sci. Lett., **1974**, V. 21, N 3, P. 253-268.
- **Weill D.F., McKay G.A.** The partitioning of Mg, Fe, Sr, Ce, Sm, Eu and Yb in lunar igneous systems and a possible origin of kreep by equilibrium partial melting //ln: Proc.Lunar Sci.Conf. 6th. Pergamon Press: **1975**, V. 1, P. 1143-1158.
- *Biggar G.* Chemistry of protopyroxene, orthopyroxene and pigeonite, crystallized from liquids to chondrule composition // Bull. Miner., **1986**, V. 109, P. 529-541.
- *Irving A.I., Merrill R.B., Singleton D.E.* Experimental partitioning of rare earth elements and scandium among armalcolite, ilmenite, olivine and mare basalt liquid // In: Proc. Lunar Planet. Sci. Conf. 9th. Pergamon Press: **1978**, V. 1, P. 601-612.
- **Delano J.W.** Chemistry and liquidus phase relations of Apollo 15 red glass: implications for the deep lunar interior // In: Proc. Lunar Planet. Sci. Conf. 11th. Pergamon Press: **1980**, V. 1, P. 251-288.
- 19 Ringwood A.E. Limits on the bulk composition of the Moon // Icarus, 1976, V. 28, N 3, P. 325-349.
- *Longhi J., Pan V.* A reconnaisance study of phase boundaries in low-alkali basaltic liquids // J. Petrol., **1988**, V. 29, Part 1, P. 115-147.
- *Grover J.E., Lindsley D.H., Bence A.E.* Experimental phase relations of olivine vitrophyres from breccia 14321: the temperature and pressure-dependence of Fe/Mg partitioning for olivine and liquid in a highlands melt-rock // In: Proc. Lunar Planet. Sci. Conf. 11th. Pergamon Press: **1980**, V. 1, P. 179-196.
- 22 Akella J., Williams R.J., Mullins O. Solubility of Cr, Ti and Al in co-existing olivine, spinel and liquid at 1 atm // In: Proc. Lunar. Sci. Conf. 7th. Pergamon Press: 1976, V. 2, P. 1179-1194.
- 23 Longhi J., Walker D., Hays J.F. The distribution of Fe and Mg between olivine and lunar basaltic liquid // Geochim. Cosmochim. Acta, 1978, V. 42, N 10, P. 1545-1558.
- *Nielsen R.L., Davidson P.M., Grove T.L.* Pyroxene-melt equilibria: an updated model // Contribs Mineral. and Petrol., **1988**, V. 100, N 3 , P. 361-373.
- 25 Drake M.J. Plagioclase-melt equilibria // Geochim. Cosmochim. Acta, 1976, V. 40, N 5, P. 457-465.
- *Walker D., Shibata T., Delong S.E.* Abyssal tholeiites from the Oceanographer Fracture Zone II. Phase equlibria and mixing // Contribs Mineral. and Petrol., **1979**, V. 70, N 2, P. 111-125.
- *Grove T.L., Gerlach D.C., Sando T.W.* Origin of calc-alkaline series lavas at Medicine Lake volcano by fractionation, assimilation and mixing // Contribs Mineral. and Petrol., **1982**, V. 80, N 2, P. 160-182.
- *Grove T.L., Bryan W.B.* Fractionation of pyroxene-phyric morb at low pressure: an experimental study // Contribs Mineral. and Petrol., **1983**, V. 84, N 4, P. 293-309.
- **Mahood G.A., Baker D.R.** Experimental constraints on depths of fractionation of midly alkalic basalts and associated felsic rocks: Pantelleria, strait of Sicily // Contribs Mineral. and Petrol., **1986**, V. 93, N 2, P. 251-264.

- *Tormey D.R., Grove T.L., Bryan W.B.* Experimental petrology of normal MORB near the Kane Fracture Zone: 22-25 N, Mid-Atlantic Ridge // Contribs Mineral. and Petrol., **1987**, V. 96, N 2, P. 121-139.
- **Sack R.O., Walker D., Carmichael I.S.E.** Experimental petrology of alkalic lavas: constraints on cotectics of multiple saturation in natural basic liquids // Contribs Mineral. and Petrol., **1987**, V. 96, N 1, P. 1-23.
- **Brey G., Green D.H.** Systematic study of liquidus phase relations in olivine melilitite+H2O+CO2 at high pressures and petrogenesis of an olivine melilitite magma // Contribs Mineral. and Petrol., **1977**, V. 61, N 2, P.141-62.
- *Thompson R.N.* Primary basalts and magma genesis II. Snake river plain, Idaho, U.S.A. // Contribs Mineral. and Petrol., **1975**, V. 52, N 13, P. 213-232.
- *Takahashi E., Scarfe C.M.* Melting of peridotite to 14GPa and the genesis of komatiite // Nature, 1985, V.315, N 6020, P. 566-568.
- **Takahashi E.** Melting of a dry peridotite KLB-1 up to 14GPa:implications on the origin of peridotitic upper mantle // J. Geophys. Res., **1986**, V. 91B, N 9, P. 9367-9382.
- **Conrad W.K., Nicholls I.A., Wall V.J.** Water-saturated and undersaturated melting of metaluminous and peraluminous crustal compositions at 1 kb: evidence for the origin of silicic magmas in the Taupo volcanic Zone, New Zeland, and other occurrences // J. Petrol., **1988**, V. 29, Part 4, P. 765-803.
- *Elthon D., Scarfe C.M.* High-pressure phase equilibria of a high-magnesia basalt and the genesis of primary oceanic basalts // Amer. Miner., **1984**, V. 69, N 1/2, P. 1-15.
- **Stolper E.** A phase diagram for mid-ocean ridge basalts: preliminary results and implications for petrogenesis // Contribs Mineral. and Petrol., **1980**, V. 74, N 1, P. 13-27.
- *Takahashi E., Kushiro I.* Melting of a dry peridotite at high pressures and basalt magma genesis // Amer. Miner., **1983**, V. 68, N 9/10, P. 859-879.
- *Jaques A.L., Green D.H.* Anhydrous melting of peridotite at 0-15 kb pressure and genesis of tholeitic basalts // Contribs Mineral. and Petrol., **1980**, V. 73, N 3, P. 287-310.
- **Bender J.F., Hodges F.N., Bence A.E.** Petrogenesis of basalts from the Project Famous Area: experimental study from 0 to 15 kbars // Earth and Planet. Sci. Lett., **1978**, V. 41, N 3, P. 277-302.
- **Green D.H.** Experimental melting studies on a model upper mantle composition at high pressure under water-saturated and water-undersaturated conditions // Earth and Planet. Sci. Lett., **1973**, V. 19, N 1, P.37-53.
- **Thompson R.N.** Primary basalts and magma genesis I. Skye, north-west Scotland // Contribs Mineral. and Petrol., **1974**, V. 45, N 4, P.317-341.
- **Bickle M.J., Ford C.E., Nisbet E.G.** The petrogenesis of peridotitic komatiites:evidence from high pressure melting experiments // Earth and Planet. Sci. Lett., **1977**, V. 37, N 1, P. 97-106.
- *Green D.H., Hibberson W.O., Jaques A.L.* Petrogenesis of mid-ocean ridge basalts // In: The Earth. Its origin structer and evolution. London: Academic Press, **1979**, P. 265-299.
- *Kushiro I., Thompson R.N.* Origin of some abyssal tholeiites from the mid-Atlantic ridge // Ann. Rept Dir. Geophys.Lab. 1971/1972, Washington D.C., **1972**, P. 403-406.
- *Fujii T., Bougault H.* Melting relations of a magnesian abyssal tholeiite and the origin of MORBs // Earth and Planet. Sci. Lett., **1983**, V. 62, N 2, P. 283-295.
- *Duncan R.A., Green D.H.* The genesis of refractory melts in the formation of oceanic crust // Contribs Mineral. and Petrol., **1987**, V. 96, N 3, P. 326-342.
- *Falloon T.J., Green D.H.* Anhydrous partial melting of MORB pyrolite and other peridotite compositions at 10 kbar: implications for the origin of primitive MORB glasses // Miner. and Petrol., 1987, V. 37, N 3-4, P. 181-219.
- *Falloon T.J., Green D.H., Hatton C.J., Harris K.L.* Anhydrous partial melting of a fertile and depleted peridotite from 2 to 30 kb and application to basalt petrogenesis // J. Petrol., **1988**, V. 29, Part 6, P.1257-1282.
- **Sen G.** Composition of basaltic liquids generated from a partially depleted lherzolite at 9 kbar pressure // Nature, **1982**, V. 299, N 5881, P. 336-338.
- Gust D.A., Perfit M.R. Phase relations of a high-Mg basalt from the Aleutian island arc: implications for primary island arc basalts and high-Al basalts // Contribs Mineral. and Petrol., 1987, V. 97, N 1, P. 7-18.
- **Shimizu H., Sangen K., Masuda A.** Experimental study on rare-earth partitioning in olivine and clinopyroxene formed at 10 and 20 kb in basaltic systems // Geochem. J., **1982**, V. 16, N 3, P. 107-117.
- **Baker D.R., Eggler D.H.** Compositions of anhydrous and hydrous melts coexisting with plagioclase, augite, and olivine or low-Ca pyroxene from 1 atm to 8 kbar: application to the Aleutian volcanic centre of Atka // Amer.Miner., 1987, V. 72, N 1/2, P. 12-28.

- **Tatsumi Y.** Origin of high-magnesian andesites in the Setouchi volcanic belt, southwest Japan.II. Melting phase relations at high pressures // Earth and Planet. Sci. Lett., **1982**, V. 60, N 2, P. 305-317.
- Johnston A.D. Anhydrous P-T phase relations of near-primary high-alumina basalt from the south Sandwich islands. Implications for the origin of island arcs and tonalite-trondhjemite series rocks // Contribs Mineral. and Petrol., 1986, V. 92, N 3, P. 368-382.
- **Spulber S.D., Rutherford M.J.** The origin of rhyolite and plagiogranite in oceanic crust: an experimental study // J. Petrol., **1983**, V. 24, Part 1, P. 1-25.
- *Mysen B.O., Kushiro I.* Compositional variations of coexisting phases with degree of melting of peridotite in the upper mantle // Amer. Miner., **1977**, V. 62, N 9/10, P. 843-856.
- **Cawthorn R.G., Davies G.** Experimental data at 3 kbars pressure on parental magma to the Bushveld Complex // Contribs Mineral. and Petrol., **1983**, V.83, N 1/2, P. 128-135.
- **Johnston A.D., Wullie P.J.** The system tonalite-peridotite-H2O at 30 kbar with applications to hybridization in subduction zone magmatism // Contribs Mineral. and Petrol., **1989**, V. 102, P. 257-264.
- **Stolper E**. Crystallization sequences of Ca-Al-rich inclusions from Allende: an experimental study // Geochim. Cosmochim. Acta, **1982**, V. 46, N 11, P. 2159-2180.
- **Stolper E., Paque J.M.** Crystallization sequences of Ca-Al-rich inclusions from Allende: the effects of cooling rate and maximum temperature // Geochim. Cosmochim. Acta, **1986**, V. 50, N 8, P. 1785-1806.
- *Arndt N.T.* Partitioning of nickel between olivine and ultrabasic and basic komatiite liquids // Ann. Rept Dir. Geophys. Lab. 1976/1977, Washington D.C., **1977**, P. 553-557.
- **Roeder P.L.** Activity of iron and olivine solubility in basaltic liquids // Earth and Planet. Sci. Lett., 1974, V. 23, N 3, P. 397-410.
- *Lukanin O.A., Kadik A.A., Biggar G.M., Fedotov S.A.* Physical-chemical conditions of crystallization of 1975-1976 BTTI basalts // Volcanology and Seysmology, **1980**, N 3, P. 16-50 (in Russian).
- *Chen H.-K., Lindsley D.H.* Apollo 14 very low titanium glasses:melting experiments in iron-platinum alloy capsules // In: Proc. Lunar and Planet. Sci. Conf. 14th. J. Geophys. Res., **1983**, V. 88B, P. 335-342.
- **Grove T.L., Juster T.C.** Experimental investigations of low-Ca pyroxene stability and olivine-pyroxene-liquid equilibria at 1-atm in natural basaltic and andesitic liquids // Contribs Mineral. and Petrol., 1989, V.103, N 3, P. 287-305.
- *Gee L.L., Sack R.O.* Experimental petrology of melilite nephelinites // J. Petrol., **1988**, V. 29, Part 6, P.1233-1255.
- *Kinzler R.J.,Grove T.L.* Crystallization and differentiation of archean komatiite lavas from northeast Ontario: phase equilibrium and kinetic studies // Amer.Miner., **1985**, V. 70, N 1/2, P. 40-51.
- **Ussler III W., Glazner A.F.** Phase equilibria along a basalt-rhyolite mixing line: implications for the origin of calc-alcaline intermediate magmas// Contribs Mineral. and Petrol., **1989**, V. 101, N 2, P. 232-244.
- **Juster T.C., Grove T.L.** Experimental constraints on the generation of the FeTi basalts, andesites, and rhyodacites at the Galapagos Spreading Center, 85 W and 95 W // J.Geophys.Res., 1989, V. 94B, N 7, P.9251-9274.
- **Biggar G.M.** The composition of diopside solid solutions, and of liquids, in equilibrium with forsterite, plagioclase, and liquid in system Na2O-CaO-MgO-Al2O3-SiO2 and in remelted rocks from 1 bar to 12 kbar // Mineral. Magazine, **1984**, V. 48, Part 4, P. 481-494.
- **Barnes S.J.** The distribution of chromium among orhtopyroxene, spinel and silicate liquid at atmospheric pressure//Geochim. Cosmochim. Acta, **1986**, V. 50, N 9, P. 1889-1909.
- *Fisk M.R., Bence A.E.* Experimental crystallization of chrome spinel in FAMOUS Basalt 527-1-1 // Earth and Planet Sci. Lett., **1980**, V. 48, N 1, P. 111-123.
- **Bartels K.S., Grove T.L.** High-pressure experiments on magnesian eucrite compositions: constrains on magmatic processes in the eucrite parent body // In: Proc. Lunar Planet. Sci. Conf. 21th, **1991**, P. 351-365.
- *Murck B.W., Campbell I.H.* The effects of temperature, oxygen fugacity and melt composition on the behaviour of chromium in basic and ultrabasic melts // Geochim. Cosmochim. Acta, **1986**, V.50, N 9, P.1871-1887.
- **Longhi J., Pan V.** The parent magmas of the meteorites // In: Proc. Lunar Planet. Sci. Conf. 19th. Cambridge University Press: **1989**, P. 451-464.
- *Hill R., Roeder P.* The crystallization of spinel from basaltic liquid as a function of oxygen fugacity // J. Geology, **1974**, V. 82, N 6, P. 709-729.
- **Sekine T.** Liquidus temperature of plagioclase and pyroxene in andesitic melts at one atmosphere // Geochem. J., **1986**, V. 20, N 6, P. 287-296.

- *Thy P.* High and low pressure phase equilibria of a mildly alkalic lava from the 1965 Surtsey eruption: Experimental results // Lithos, **1991**, V. 26, P. 223-243.
- *Thy P., Lofgren G.E., Imsland P.* Melting relations and the evolution of the Jan Mayen magma system // J. Petrol., 1991, V. 32, Part 2, P. 303-332.
- *Fram M.S., Longhi J.* Phase equilibria of dikes associated with Proterozoic anorthosite complexes // Amer. Miner., **1992**, V. 77, N 5/6, P. 605-616.
- **Snyder D., Carmichael I.S.E., Wiebe R.A.** Experimental study of liquid evolution in an Fe-rich, layered mafic intrusion: Constraints of Fe-Ti oxide precipitation on the T-fO2 and T-p paths of tholeitic magmas // Contribs Mineral. and Petrol., **1992**, V. 113, N 1, P. 73-86.
- *Grove T.L., Kinzler R.J., Bryan W.B.* Fractionation of mid-ocean ridge basalt. In RIDGE Summer Institute Volume, AGU Special Publication, **1992** (in press).
- *Kinzler R.J., Grove T.L.* Primary magmas of mid-ocean ridge basalts,I: Experiments and methods // J. Geophys. Res., **1992**, V. 97B, N 5, P.6885-6906.
- *Draper D.S., Johnston A.D.* Experimental generation of arc-like high-alumina basalt from anhydrous, primitive olivine tholeiite: An experimental study from 1 atm to 20 kbar // Contribs Mineral. and Petrol., **1992**, V. 112, N 4, P. 501-519.
- *Meen J.K.* Elevation of potassium content of basaltic magmas by fractional crystallization: the effect of pressure // Contribs Mineral. and Petrol., **1990**, V. 104, N 3, p. 309-331.
- *Meen J.K.* Formation of shoshonites from calcalkaline basalt magmas: geochemical and experimental constraints from the type locality // Contribs Mineral. and Petrol., **1987**, V. 97, N 3, P. 333-351.
- *Fisk M.R., Upton B.G.J., Ford C.E., White W.M.* Geochemical and experimental study of the genesis of magmas of Reunion Island, Indian Ocean // J. Geophys. Res., **1988**, V. 93B, N 5, P. 4933-4950.
- **Shi P.,Libourel G.** The effect of FeO on the system CMAS at low pressure and implications for basalt crystallization processes // Contribs Mineral. and Petrol., **1991**, V. 108, N 1/2, P. 129-145.
- **Bartels K.S., Kinzler R.J., Grove T.L.** High pressure phase relations of primitive high-alumina basalts from Medicine Lake volcano, northern California // Contribs Mineral. and Petrol., **1991**, V. 108, N 3, P. 253-270.
- *Grove T.L., Kinzler R.J., Bryan W.B.* 2. Natural and experimental phase relations of lavas from Serocki volcano // Proceedings of the Ocean Drilling Program, Scientific Results, **1990**, V. 106/109, P. 9-17.
- *Roeder P.L., Reynolds I.* Crystallization of chromite and chromium solubility in basaltic melts // J. Petrol., **1991**, V. 32, Part 5, P. 909-934.
- *Ryerson F.J., Weed H.C., Piwinskii A.J.* Pheology of subliquidus magmas 1.Picritic compositions // J. Geophys.Res.,**1988**, V. 93B, N 4, P. 3421-3436.
- *Eggins S.M.* Petrogenesis of Hawaiian tholeiites: 1. Phase equilibria constraints // Contribs Mineral. and Petrol., **1992**, V. 110, N 2/3, P. 387-397.
- *Thy P., Lofgren G.E.* Experimental constraints on the low-pressure evolution of transitional and mildly alkalic basalts: multisaturated liquids and coexisting augites // Contribs Mineral. and Petrol., **1992**, V.112, N 2/3, P. 196-202.
- *Forsythe L.M., Fisk M.R.* Comparison of experimentally crystallized and natural spinels from leg 135 //ln: Proc.ODP Sci.Results. Leg 135: College Station, TX (submitted **1992**).
- **Gaetani G.A., Grove T.L., Bryan W.B.** Experimental phase relations of basaltic andesite from hole 839B under hydrous and anhydrous conditions // In: Proc. ODP Sci. Results. Leg 135: College Station, TX (submitted **1992**).
- *Walker D., Jurewicz S., Watson E.B.* Adcumulus dunite growth in a laboratory thermal gradient // Contribs Mineral. and Petrol., **1988**, V. 99, N 3, P. 306-319.
- *Kilinc A., Mehmet Z.C.* Experimental study of low pressure mineral-melt equilibria in alkaline lavas (unpublished data, **1989**).
- 101 Forsythe L., Nielsen R.L. Experimental study of magnetite-melt equilibria in andesite-basaltic systems (unpublished data, 1992).
- **Nielsen R.L.** Experimental study of high temperature plagioclase-melt equilibria in high magnesian tholeiitic systems (unpublished data,**1992**).
- *Hirose K., Kushiro I.* Partial melting of dry peridotites at high pressures: determination of compositions of melts segregated from peridotite using aggregates of diamond // Earth and Planet. Sci. Lett., **1993**, V.114, N 4, P. 477-489.
- *Houghton D.R., Roeder P.L., Scinner B.J.* Solubility of sulfur in mafic magmas // Econ. Geol., **1974**, V. 69, N 4, P. 451-467.
- **Shima H., Naldrett A.J.** Solubility of sulfur in an ultramafic melt the relevance of the system Fe-S-O // Econ.Geol., **1975**, V.70, P. 960-967.

- **Buchanan D.L., Nolan J.** Solubility of sulfur and sulfide immiscibility in synthetic tholeitic melts and their relevance to Bushveld-complex rocks // Can. Mineral., **1979**, V. 17, P. 483-494.
- **Buchanan D.L., Nolan J., Wilkinson N., de Villiers J.P.R.** An experimental investigation of sulfur solubility as a function of temperature in synthetic silicate melts // Geol. Soc. South Afr., Spec. Publ., **1983**, V. 7, P. 383-391.
- *Rutherford M.J.* The May 18, 1980, Eruption of Mount St. Helens 1.Melt composition and experimental phase equilibria // J.Geophys.Res., **1985**, V. 90, N B4, P.2929-2947.
- *Olafsson M.* Partial melting of peridotite in the presence of small amounts of volatiles, with special reference to the low-velocity zone // M.S.Thesis, Pennsylvania State University, **1980**, 59 pp.
- **Johnston A.D., Wyllie P.J.** Constraints on the origin of Archean trondhjemites based on phase relationships of Nuk gneiss with H2O at 15 kbar // Contribs Mineral. and Petrol., **1988**, V. 100, N 1, P. 35-46.
- *Luhr J.F.* Experimental phase relations of water- and sulfur-saturated arc magmas and the 1982 eruptions of El Chicon Volcano // J. Petrol., **1990**, Part 5, P. 1071-1114. & *Housh T.B., Luhr J.F.* Plagioclase-melt equilibria in hydrous systems // Amer.Miner., **1991**, V. 76, P. 472-492.
- **Beard J.S.,Lofgren G.E.** Dehydration melting and water-saturated melting of basaltic and andesitic greenstones and amphibolites at 1, 3, and 6-9 kb // J. Petrol., **1991**, V. 32, Part 2, P. 365-401 (plus unpublished data).
- **Patino D.A.E., Johnston A.D.** Phase equilibria and melt productivity in the pelitic system: implications for the origin of peraluminous granitoids and aluminous granulites // Contribs Mineral. and Petrol., **1991**, V. 107, P. 202-218.
- *Rushmer T.* Partial melting of two amphibolites:contrasting experimental results under fluid-absent conditions // Contribs Mineral. and Petrol., **1991**, V. 107, P. 41-59.
- *Holtz F., Pichavant M., Barbey P., Johannes W.* Effects of H2O on liquidus phase relations in the haplogranite system at 2 and 5 kbar // Amer. Miner., **1992**, V. 77, N 11/12, P. 1223-1241.
- **Gerke T.L., Kilinc A.I.** Enrichment of SiO2 in rhyolites by fractional crystallization: an experimental study of peraluminous granitic rocks from the St. Francois Mountains, Missouri, USA // Lithos, **1993**, V. 29, N 3/4, P. 273-283.
- **Sisson T.W., Grove T.L.** Temperatures and H2O contents of low-MgO high-alumina basalts // Contribs Mineral. and Petrol., **1993**, V. 113, N 2, P. 167-184.
- **Sisson T.W., Grove T.L.** Experimental investigations of the role of H2O in calc-alkaline differentiation and subduction zone magmatism // Contribs. Mineral. and Petrol., **1993**, V. 113, N 2, P. 143-166.
- **Sweeney R.J., Thompson A.B., Ulmer P.** Phase relations of a natural MARID composition and implications for MARID genesis, lithospheric melting and mantle metasomatism // Contribs Mineral. and Petrol., **1993**, V. 115, P. 225-241.
- **Wolf M.B., Wyllie P.J.** Dehydration-melting of amphibolite at 10 kbar: the effects of temperature and time // Contribs Mineral. and Petrol., **1994**, V. 115, P. 369-383.
- *Tatsumi Y.* Melting experiments on a high-magnesian andesite // Earth Planet. Sci. Lett., **1981**, V. 54, P.37-53.
- *Tatsumi Y.* Generation of Arc Basalt Magmas and Thermal Structure of the Mantle Wedge in Subduction Zones // J. Geophys. Res., **1983**, V. 88, N B7, P. 5815-5825.
- *Kushiro I.* Partial melting of mantle wedge and evolution of island arc crust // J. Geophys. Res., **1990**, V.95, N B10, P. 15929-15939.
- 124 Edgar A.D., Arima M., Baldwin D.K., Bell D.R., Shee S.R., Skinner E.M.W., Walker E.C. High-pressure high-temperature melting experiments on a SiO2-poor aphanitic kimberlite from the Wesselton mine, Kimberley, South Africa // Amer. Miner., 1988, V. 73, P. 524-533.
- *Edgar A.D., Charbonneau H.E.* Melting experiments on a SiO2-poor, CaO-rich aphanitic kimberlite from 5-10 GPa and their bearing on sources of kimberlite magmas // Amer. Miner., **1993**, V. 78, N 1/2, P. 132-142.
- *Umino S., Kushiro I.* Experimental studies on boninite petrogenesis // Boninites (ed. A.J.Crawford), London, Unwin Hyman, **1989**, P. 89-111.
- 127 Van der Laan S.R., Flower M.F.J., Koster Van Groos A.F. Experimental evidence for the origin of boninites: near-liquidus phase relations to 7.5 kbar // Boninites (ed. A.J.Crawford), London, Unwin Hyman, 1989, P. 112-147.
- **Mysen B.O., Boettcher A.L.** Melting of a hydrous mantle: I. Phase relations of natural peridotite at high pressures and temperatures with controlled activities of water, carbon dioxide, and hydrogen // J. Petrol.,
 - , V. 16, Part 3, P. 520-548. + 129 (CONTINUED)

- **Mysen B.O., Boettcher A.L.** II. Geochemistry of crystals and liquids formed by anatexis of mantle peridotite at high pressures and temperatures as a function of controlled activities of water, hydrogen, and carbon dioxide // J. Petrol., **1975**, V.16, Part 3, P. 549-593. + 128 (BEGINNING)
- *Mysen B.O., Boettcher A.L.* Melting of a hydrous mantle: III. Phase relations of garnet websterite+H2O at high pressures and temperatures // J. Petrol., **1976**, V. 17, Part 1, P. 1-14.
- 131 Longhi J., Wooden J.L., Coppinger K.D. The petrology of High-Mg dikes from the Beartooth Mountains, Montana: a search for the parent magma of the Stillwater Complex // J.Geophys.Res., 1983, V. B88 (Supplement), P. 53-69 (+ unpublished data).
- **Baker D.R., Eggler D.H.** Fractionation paths of Atka (Aleutians) high-alumina basalts: constraints from phase relations // J. Volcanol. Geotherm. Res., **1983**, V. 18, P. 387-404.
- *Helz R.T.* Phase relations of basalts in their melting ranges at PH2O=5 kb as a function of oxygen fugacity. Part I. Mafic phases // J.Petrol., **1973**, V. 14, Part 2, P. 249-302. & Part II. Melt compositions // J. Petrol., **1976**, V. 17, Part 2, P. 139-193.
- *Holloway J.R., BurnhamC.W.* Melting relations of basalt with equilibrium water pressure less than total pressure // J. Petrol., **1972**, V. 13, Part 1, P. 1-29.
- *Allen J.C., Boettcher A.L., Marland G.* Amphiboles in andesite and basalt: I. Stability as a function of P-T-fO2 // Amer. Miner., **1975**, V. 60, N 11/12, P. 1069-1085.
- *Allen J.C., Boettcher A.L.* Amphiboles in andesite and basalt: II. Stability as a function of P-T-fH2O-fO2 // Amer. Miner., 1975, V.63, N 11/12, P. 1074-1087.
- **Nehru C.E., Wyllie P.J.** Compositions of glasses from St. Paul's peridotite partially melted at 20 kilobars // J. Geology, **1975**, V. 83, N 4, P. 455-471.
- **Sekine T., Wyllie P.J.** The system granite-peridotite-H2O at 30 kbar, with applications to hybridization in subduction zone magmatism // Contribs Mineral. and Petrol., **1982**, V. 81, N 3, P. 190-202.
- *Carroll M.R., Wyllie P.J.* The system tonalite-H2O at 15 kbar and the genesis of calc-alkaline magmas // Amer. Miner., 1990, V. 75, P. 345-357.
- *Adam J.* The geochemistry and experimental petrology of sodic alkaline basalts from Oatlands, Tasmania // J. Petrol., **1990**, V. 31, Part 6, P. 1201-1223.
- *Burnham C.W., Jahns R.H.* A method for determining the solubility of water in silicate melts // Amer. J. Sci., **1962**, V. 260, N 10, P. 721-745.
- *Hamilton D.L., Burnham C.W., Osborn E.F.* The solubility of water and effects of oxygen fugacity and water content on cristallization in mafic magmas // J. Petrol., 1964, V. 5, N 2, P. 21-39.
- *Khitarov N.I., Kadik A.A., Lebedev E.B.* Solubility of water in basalt melt // Geokhimiya, **1968**, N 7, P.763-772 (in Russian).
- **Bohlen S.R., Boettcher A.L., Wall V.J.** The system albite-H2O-CO2:a model for melting and activities at high pressures // Amer. Mineral., **1982**, V. 67, N 5/6, P. 451-462.
- *Mironov A.G., Epel'baum M.B., Chekhmir A.S.* Experimental determination of the relative water solubility in the granite and basalt melts at 900-1100°C and 2 kbar with tritium autoradiographic method // Geokhimiya, **1993**, N 4, P. 487-498 (in Russian).
- **Nesbitt R.W., Hamilton D.L.** Cristallisation of an alkali-olivine basalt under controlled pO2, pH2O conditions // Phys. Earth Planet. Interiors, **1970**, V. 3, P. 309-315.
- *Fudali R.F.* Oxyden fugacities of basaltic and andesitic magma // Geochim. Cosmochim. Acta, **1965**, V.29, N 9, P. 1063-1075.
- **Shibata K.** The Oxyden Partial Pressure of the Magma from Mihara Volcano, O-sima, Japan // Bull. Chem. Soc. Japan, **1967**, V. 40, N 4, P. 830-834.
- *Thornber C.R., Roeder P.L., Foster J.R.* The effect of composition on the ferric-ferrous ratio in basaltic liquids at atmospheric pressure // Geochim. et Cosmochim. Acta, **1980**, V. 44, N 3, P. 525-532.
- **Sack R.O., Carmichael I.S.E., Rivers M., Ghiorso M.S.** Ferric-Ferrous Equilibria in Natural Silicate Liquids at 1 Bar // Contribs Mineral. and Petrol., **1980**, V. 75, N 4, P. 369-376.
- *Kilinc A., Carmichael I.S.E., Rivers M.L., Sack R.O.* The Ferric-Ferrous Ratio of Natural Silicate Liquids Equilibrated in Air // Contribs Mineral. and Petrol., **1983**, V. 83, N 1/2, P. 136-140.
- *Kress V.C., Carmichael I.S.E.* Stoichiometry of the iron oxidation reaction in silicate melts // Amer. Miner., 1988, V. 73, N 11/12, P. 1267-1274.
- *Kress V.C., Carmichael I.S.E.* The compressibility of silicate liquids containing Fe2O3 and the effect of composition, temperature, oxygen fugacity and pressure on their redox states // Contribs Mineral. and Petrol., **1991**, V. 108, N 1/2, P. 82-92.
- *Fujii T., Scarfe C.M.* Composition of liquids coexisting with spinel lherzolite at 10 kbar an the genesis of MORBs // Contribs Mineral. and Petrol., **1985**, V. 90, N 1, 18-28.

- 155 Wei K., Tronnes R.G., Scarfe C.M. Phase Relations of Aluminum-Undepleted and Aluminum-Depleted Komatiites at Pressures of 4-12 GPa // J.Geophys. Res., 1990, V. 95, N B10, P.15817-15827.
- *Canil D.* Orthopyroxene stability along the peridotite solidus and the origin of cratonic lithosphere beneath southern Africa // Earth Planet. Sci. Lett., **1992**, V. 111, N 1, P.83-95.
- *Kadik A.A., Lebedev E.B, Khitarov N.I.* The water in magmatic melts // Nauka Publishing House, Moscow, **1971**, 267 P. (in Russian)
- **Baker M.B., Stolper E.M.** Determining the composition of high-pressure mantle melts using the diamond aggregates // Geochim. Cosmochim. Acta, **1994**, V. 58, N 13, P. 2811-2827.
- *Falloon T.J., Green D.H., Jaques A.L.* Refractory magmas in back-arc basin settings: experimental constraints on a Lau Basin exemple // Miner. Magazine (Extended abstracts to V.M. Goldschmidt conference, Edinburg), **1994**, V. 58A, P. 263-264.
- *Herzberg C., Gasparik T., Sawamoto H.* Origin of Mantle Peridotite: Constraints From Melting Experiments to 16.5 GPa // J.Geophys. Res., **1990**, V. 95, N B10, P. 15779-15803.
- *Herzberg C., Gasparik T.* Garnet and pyroxenes in the mantle: a test of the majorite fractionation hypothesis // J.Geophys. Res., **1991**, V. 96, N B10, P. 16263-16274.
- *Herzberg C.* Depth and degree of melting of komatiites // J.Geophys. Res., **1992**, V. 97, N B4, P. 4521-4540.
- **Stolper E., McSween Jr. H.Y.** Petrology and origin of the shergottite meteorites // Geochim. et Cosmochim. Acta, **1979**, V. 43, N ?, P.1475-1498.
- *Jurewicz A.J.G., Mittlefehldt D.W., Jones J.H.* Partial melting of the Allende (CV3) meteorite: Implications for origins of basaltic meteorites // Science, **1991**, V. 252, P. 695-698.
- 165 Jurewicz A.J.G., Mittlefehldt D.W., Jones J.H. Experimental partial melting of the Allende (CV) and Murchison (CM) chondrites and the origin of asteroidal basalts // Geochim. Cosmochim. Acta, 1993, V. 57, N 9, P. 2123-2139.
- *Jurewicz A.J.G., Mittlefehldt D.W., Jones J.H.* Experimental partial melting of the St. Severin (LL) and Lost City (H) chondrites // Geochim. Cosmochim. Acta, **1995**, V. 59, N 2, P. 391-408.
- 167 Walker D., Grove T. Ureilite smelting // Meteoritics, 1993, V. 28, N5, P. 629-636.
- *Tronnes R.G., Canil D., Wei K*. Element partitioning between silicate minerals and coexisting melts at pressures of 1-27 GPa, and implications for mantle evolution // Earth and Planet. Sci. Lett., **1992**, V. 111, N 2/4, P. 241-255.
- *Kennedy A.K., Grove T.L., Johnson R.W.* Experimental and major element constraints on the evolution of lavas from Lihir Island, Papua New Guinea // Contrib. Mineral. and Petrol., **1990**, V. 104, N6, P.722-734.
- *Kuehner S.M.* Experimental determination of the emplacement pressure of Proterozoic mafic dikes within portions of the East Antarctic Shield // Lithos, **1992**, V. 28, N 1, P. 21-42.
- *Pan V., Longhi J.* The system Mg2SiO4-Ca2SiO4-CaAl2O4-NaAlSiO4-SiO2: one atmosphere liquidus equilibria of analogs of alkaline mafic lavas // Contrib. Miner. and Petrol., **1990**, V. 105, P.569-584.
- **Soulard H., Provost A., Boivin P.** CaO-MgO-Al2O3-SiO2-Na2O (CMASN) at 1 bar from low to high Na2O contents: Topology of an analogue for alkaline basic rocks // Chemical Geology, **1992**, V. 96, P. 459-477.
- **Shi P.** Low-pressure phase relationships in the system Na2O-CaO-FeO-MgO-Al2O3-SiO2 at 1100°C, with implications for the differentiation of basaltic magmas // J. Petrol., **1993**, V. 4, N 4, P. 743-762.
- *Thy P., Lofgren G.E.* Experimental constraints on the low-pressure evolution of transitional and mildly alkalic basalts: the effect of Fe-Ti oxide minerals and the origin of basaltic andesites // Contrib. Mineral. and Petrol., **1994**, V. 116, N 3, P. 340-351.
- **Bertka C.M., Holloway J.R.** Anhydrous partial melting of an iron-rich mantle I: subsolidus phase assemblages and partial melting phase relations at 10 to 30 kbar. + II: primary melt compositions at 15 kbar // Contrib. Mineral. and Petrol., 1994, V.115, P.313-322 + P.323-330.
- *Toplis M.J., Libourel G., Carroll M.R.* The role of phosphorus in crystallization processes of basalt: an experimental study // Geochem. Cosmochim. Acta, **1994**, V. 58, N 2, P. 797-810.
- **Auwera J.V., Longhi J.** Experimental study of a jotunite (hypersthene monzodiorite): constraints on the parent magma composition and crystallization conditions (P, T, fO2) of the Bjerkreim-Sokndal layered intrusion (Norway)//Contrib.Miner. and Petrol., **1994**, V.118, N1, P.60-78.
- **Baker M.B., Grove T.L., Price R.** Primitive basalts and andesites from the Mt.Shasta region, N.California: products of varying melt fraction and water content // Contrib. Mineral. and Petrol., **1994**, V. 118, N 2, P. 111-129.
- *Nielsen R.L., Forsythe L.M., Gallahan W.E., Fisk M.R.* Major- and trace-element magnetite-melt equilibria // Chemical Geology, **1994**, V. 117, P. 167-191.

- *Panjasawatwong Y.,Danyushevsky L.V.,Crawford A.J.,Harris K.L.* An experimental study of the effects of melt composition on plagioclase-melt equilibria at 5 and 10 kbar: implications for the origin of magmatic high-An plagioclase // Contrib.Mineral.Petrol., **1995**, V. 118, N 4, P. 420-432.
- *Longhi J.* Liquidus equilibria of some primary lunar and terrestrial melts in the garnet stability field // Geochim. et Cosmochim. Acta, **1995**, V. 59, N 11, P. 2375-2386.
- **Baker M.B., Hirschmann M.M., Ghiorso M.S., Stolper E.M.** Compositions of near-solidus peridotite melts from experiments and thermodynamic calculations // Nature, **1995**, V. 375, P. 308-311.
- **Kushiro I.** Partial melting of a fertile mantle peridotite at high pressures: an experimental study using aggregates of diamond // AGU Monograph 95. Earth Processes: Reading the Isotopic Code, **1996**, P. 109-122.
- *Takahashi E.* Melting of a Yamato L3 chondorite (Y-74191) up to 30 kb // Proc. of 8th Symp. on Antarctic Meteorites, Mem. Nat. Inst. Polar Res., Spec. Issue N 30, **1983**, P. 168-180.
- *Kato T., Ringwood A.E., Irifune T.* Experimental determination of element partitioning between silicate perovskites, garnets and liquids: constraints on early differentiation of the mantle // Earth Planet. Sci. Lett., **1988**, V. 89, N 1, P. 123-145.
- **Yasuda A., Fujii T., Kurita K.** Melting phase relations of an anhydrous mid-ocean ridge basalt from 3to 20 GPa: implications for the behavior of subducted oceanic crust in the mantle // J.Geoph.Res., **1994**, V. 99, N B5, P. 9401-9414.
- **Zhang J., Herzberg C.** Melting experiments on anhydrous peridotite KLB-1 from 5.0 to 22.5 GPa // J.Geophys.Res., **1994**, V. 99, N B9, P. 17729-17742.
- *Ohtani E., Nagata Y., Suzuki A., Kato T.* Melting relations of peridotite and the density crossover in planetary mantles // Chem. Geol., **1995**, V. 120, N 3-4, P. 207-221.
- **Sack R.O., Carmichael I.S.E.** Fe2+=Mg2+ and TiAl2=MgSi2 exchange reactions between clinopyroxenes and silicate melt // Contribs Mineral. and Petrol., **1984**, V. 85, N ??, P. 103-115.
- **Sack R.O., Ghiorso M.S.** Thermodynamics of multicomponent pyroxenes: III. Calibration of Fe2+(Mg)-1, TiAl2(MgSi2)-1, TiFe3+(MgSi2)-1, ... and Ca(Mg)2-1 exchange reactions between pyroxenes and silicate melts // Contribs Mineral. and Petrol., **1994**, V. 118, N 3, P. 271-296.
- 191 Wagner T.P., Donnelly-Nolan J.M., Grove T.L. Evidence of hydrous differentiation and crystal accumulation in the low-MgO, high-Al2O3 Lake basalt from Medicine Lake volcano, California // Contribs Mineral. and Petrol., 1995, V. 121, N ??, P. 201-216.
- *Moore G., Righter K., Carmichael I.S.E.* The effect of dissolved water on the oxidation state of iron in natural silicate liquids // Contribs Mineral. and Petrol., **1995**, V. 120, N 1, P. 170-179.
- **Toplis M.J., Carroll M.R.** An experimental study of the influence of oxygen fugacity on Fe-Ti oxide stability, phase relations, and mineral-melt equilibria in ferro-basaltic systems // J. Petrol., **1995**, V. 36, N 5, P. 1137-1170.
- *Thy P.* Low-pressure experimental constraints on the evolution of komatiites // J. Petrol., **1995**, V. 36, N 6, P.1529-1548.
- *Girnis A.V.,Brey G.P.,Ryabchikov I.D.* Origin of group 1A kimberlites: fluid-saturated melting experiments at 45-55 kbar // Earth Planet.Sci. Lett., **1995**, V. 134, N 1, P. 283-296.
- *Hirose K., Kawamoto T.* Hydrous partial melting of lherzolite at 1 Gpa: The effect of H2O on the genesis of basaltic magmas // Earth Planet. Sci. Lett., **1995**, V.133, N ??, P.463-473.
- **Thy P.** Experimental constraints on the evolution of transitional and mildly alkalic basalts: crystallization of spinel // Lithos, **1995**, V. 36, N ??, P.103-114.
- *Klingenberg B.M.E.T., Kushiro I.* Melting of a chromite-bearing harzburgite and generation of boninitic melts at low pressures under controlled oxygen fugacity // Lithos, **1996**, V. 37, N ??, P. 1-13.
- *Toplis M.J., Dingwell D.B., Libourel G.* The effect of phosphorus on the iron redox ratio, viscosity, and density of an evolved ferrobasalt // Contribs Mineral. and Petrol., **1994**, V. 117, N ??, P. 293-304.
- *Dixon J.E., Stolper E.M., Holloway J.R.* An experimental study of water and carbon dioxide solubilities in mid-ocean ridge basaltic liquids. Part 1: Calibration and solubility models // J. Petrol., 1995, V. 36, N 6, P. 1607-1631.
- *Kadik A.A., Maksimov A.P. Ivanov B.V.* Physical-chemical conditions of crystallization and the genesis of andesits (an example of the Klyuchevskoy group volcanoes) // Moscow, Nauka publ., 1986. 158 pp. (in Russian).
- *Biggar G.M.* Phase equilibrium studies of the chilled margins of some layered intrusions // Contribs Mineral. and Petrol., **1974**, V. 46, N 2, P. 159-167.
- *Hoover J.D.* The chilled marginal gabbro and other contact rocks of the Scaergaard intrusion // J. Petrol., **1989**, V. 30, N 2, P. 441-476.
- *Putirka K., Johnson M., Kinzler R., Longhi J., Walker D.* Thermobarometry of mafic igneous rocks based on clinopyroxene-liquid equilibria, 0-30 kbar // Contribs Mineral. and Petrol., **1996**, V. 123, N ??, P. 92-108.

- *Righter K., Carmichael I.S.E.* Phase equilibria of phlogopite lamprophyres from western Mexico: biotite-liquid equilibria and P-T estimates for biotite-bearing igneous rocks // Contribs Mineral. and Petrol., **1996**, V. 123, N 1, P. 1-21.
- *Merzbacher C., Eggler D.H.* A magmatic geohygrometer: Application to Mount St.Helens and other dacitic magmas // Geology, **1984**, V. 12, N 10, P. 587-590.
- *Huang W.-L., Wyllie P. J.* Phase relationships of gabbro-tonalite-granite-water at 15 kbar with applications to differentiation and anatexis // Amer. Miner., 1986, V. 71, N 3/4, P. 301-316.
- *Clemens J.D., Holloway J.R., White A.J.R.* Origin of an A-type granite: Experimental constraints // Amer. Miner., **1986**, V. 71, N 3/4, P. 317-324.
- **Skjerlie K.P., Douce A.E.P., Johnston A.D.** Fluid absent melting of a layered crustal protolith:implications for the generation of anatectic granites//Contribs Mineral. and Petrol., **1993**, V. 114, N 3, P. 365-378.
- **Skjerlie K.P., Douce A.E.P.** Anatexis of interlayered amphibolite and pelite at 10 kbar: effect of diffusion of major components on phase relations and melt fraction // Contribs Mineral. and Petrol., **1995**, V. 122, N ??, P. 62-78.
- **Beard J.S., Lofgren G.E., Sinha A.K., Tollo R.P.** Partial melting of apatite-bearing charnockite,granulite, and diorite: Melt compositions, restite mineralogy, and petrologic implications // J. Geophys. Res., **1994**, V. 99B, N 11, P. 21591-21603.
- *Vielzeuf D., Montel J.M.* Partial melting of metagreywackes. Part 1. Fluid-absent experiments and phase relationships // Contribs Mineral. and Petrol., **1994**, V. 117, N ??, P. 375-393.
- **Sen C., Dunn T.** Dehydration melting of a basaltic composition amphibolite at 1.5 and 2.0 GPa: implications for the origin of adakites // Contribs Mineral. and Petrol., **1994**, V. 117, N 4, P. 394-409.
- **Scaillet B., Pichavant M., Roux J.** Experimental crystallization of leucogranite magmas // J. Petrol., **1995**, V. 36, N 3, P. 663-705.
- **Johnston A.D., Wyllie P.J.** Interaction of granitic and basic magmas: experimental observations on contamination processes at 10 kbar with H2O // Contribs. Mineral. and Petrol., **1988**, V. 98, N 3, P. 352-362.
- *Vielzeuf D., Holloway J.R.* Experimental determination of the fluid-absent melting relations in the pelitic system (consequences for crustal differentiation) // Contribs Mineral. and Petrol., **1988**, V. 98, N 3, P. 257-276.
- **Tsuchiyama A.** Experimental study of olivine-melt reaction and its petrological implications // J. Volcanol. Geotherm. Res., **1986**, V. 29, N 1-4, P. 245-264.
- **Wolf M.B., Wyllie P.J.** Dehydration-melting of solid amphibolite at 10 kbar: Textural development, liquid interconnectivity and applications to the segregation of magmas // Mineral. and Petrol., **1991**, V. 44, N ??, P. 151-179.
- *Bradford S.W., Wyllie P.J.* Solidus reactions in synthetic Iherzolite-H2O-CO2 from 20-30 kbar, with applications to melting and metasomatism // J. Volcanol. Geotherm. Res., **1992**, V. 50, N ??, P. 117-130.
- *Van der Laan S.R., Wyllie P.J.* Constraints on archean trondhjemite genesis from hydrous crystallization experiments on Nuk gneiss at 10-17 kbar // J. Geology, **1992**, V. 100, N ??, P. 57-68.
- *Wolf M.B., Wyllie P.J.* Garnet growth during amphibolite anatexis: Implications of a garnetiferous restite // J. Geology, **1993**, V. 101, N ??, P. 357-373.
- *Rapp R.P., Watson E.B., Miller C.F.* Partial melting of amphibolite/eclogite and the origin of Archean trondhjemites and tonalites // Precambrian Research, **1991**, V. 51, N 1-4, P. 1-25.
- *Holloway J.R., Burnham C.W.* Melting relations of basalt with equilibrium water pressure less than total pressure // J. Petrol., **1972**, V. 13, Part 1, P. 1-29.
- *Allen J.C., Boettcher A.L.* The stability of amphibole in andesite and basalt at high pressures // Amer. Miner., 1983, V. 68, P. 307-314.
- *Hacker B.R.* Amphibolite-facies-to-granulite-facies reactions in experimentally deformed, unpowdered amphibolite // Amer. Miner., **1990**, V. 75, N ??, P. 1349-1361.
- *Blank J.G., Stolper E.M., Carroll M.R.* Solubilities of carbon dioxide and water in rhyolitic melt at 850°C and 750 bars // Earth Planet. Sci. Lett., **1993**, V. 119, N ??, P. 27-36.
- *Rutherford M.J., Devine J.D.* The May 18, 1980, Eruption of Mount St. Helens: 3. Stability and chemistry of amphibole in the magma chamber // J. Geophys. Res., **1988**, V. 93B, N 10, P. 11949-11959.
- *Johnson M.C., Rutherford M.J.* Experimentally determined conditions in the Fish Canyon Tuff, Colorado, Magma Chamber // J. Petrol., **1989**, V. 30, Part 3, P. 711-737.
- *Johnson M.C., Rutherford M.J., Hess P.C.* Chassigny pertogenesis: melt compositions, intensive parameters, and water contents of Martian (?) magmas // Geochim. Cosmochim. Acta, 1991, V. 55, N 1, P. 349-366.

- 230 *Rutherford M.J., Hill P.M.* Magma ascent rates from amphibole breakdown: An experimental study applied to the 1980-1986 Mount St. Helens eruptions // J. Geophys. Res., **1993**, V. 98B, N 11, P. 19667-19685.
- 231 **Baker L.L., Rutherford M.J.** The effect of dissolved water on the oxidation state of silicic melts // Geochim. Cosmochim. Acta, **1996**, V. 60, N 12, P. 2179-2187.
- 232 **Beard J.S., Abitz R.J., Lofgren G.E.** Experimental melting of crustal xenoliths from Kilbourne Hole, New Mexico and implications for the contamination and genesis of magmas // Contribs Mineral. and Petrol.. **1993**. V. 115. N 1. P. 88-102.
- Petrol., **1993**, V. 115, N 1, P. 88-102. 233 *Ono S., Yasuda A.* Compositional change of majoritic garnet in a MORB composition from 7 to 17 GPa and 1400 to 1600°C // Phys. Earth and Planet. Interiors, **1996**, V. 96, N ??, P. 171-179.
- 234 *Cawthorn R.G., Biggar G.M.* Crystallization of titaniferous chromite, magnesian ilmenite and armalcolite in tholeiitic suites in the Karroo Igneous Province // Contribs Mineral. and Petrol., **1993**, V. 114, N ??, P. 221-235.
- 235 *Cawthorn R.G.* Re-evaluation of magma compositions and processes in the uppermost Critical Zone of the Bushveld Complex // Mineral. Magazine, **1996**, V. 60, N ??, P. 131-148.

List of INFOREX-4.2 (2001) experimental studies added to the main archive

- 236 Yang H.-J., Kinzler R.J., Grove T.L. Experiments and models of anhydrous, basaltic olivine-plagioclase-augite saturated melts from 0.001 to 10 kbar // Contribs Mineral. and Petrol., 1996, V. 124, N 1, P. 1-18.
- 237 Kawamoto T. Experimental constraints on differentiation and H2O abundance of calc-alkaline magmas // Earth and Planet. Sci. Lett., 1996, V. 144, N ??, P. 577-589.
- 238 Khodorevskaya L.I., Zharikov V.A. Experimental simulation of amphibolit and ultrabasic rock interaction in subduction zones // Petrology, 1997 V. 5, N 1, P. 2-7.
- 239 Patino Douce A.E., Beard J.S. Effects of P, f(O2) and Mg/Fe ratio on dehydration melting of model metagreywackes // J.Petrol., 1996, V.37, N 5, P. 999-1024.
- 240 Patino Douce A.E. Effects of pressure and H2O content on the compositions of primary crustal melts // Transactions of the Royal Society of Edinburgh: Earth Sciences, 1996, V. 87, N ??, P. 11-21.
- 241 Rapp R.P. Amphibole-out phase boundary in partially melted metabasalt, its control over liquid fraction and composition, and source permeability // J. Geophys. Res., 1995, V. 100, N 8B, P. 15601-15610.
- Grove T.L., Donnelly-Nolan J.M., Housh T. Magmatic processes that generated the rhyolite of Glass Mountain, Medicine Lake volcano, N.California // Contribs Mineral. and Petrol., 1997, V. 127, N ??, P. 205-223.
- 243 Wagner T.P., Grove T.L. Experimental constraints on the origin of luna high-Ti ultramafic glasses // Geochim.Cosmochim.Acta, 1997, V. 61,N 6, P. 1315-1327.
- 244 Gudfinnsson G.H., Presnall D.C. Melting relations of model lherzolite in the system CaO-MgO-Al2O3-SiO2 at 2.4-3.4 GPa and the generation of komatiites // J. Geophys. Res., 1996, V. 101, N B12, P. 27701-27709.
- 245 Skjerlie K.P., Johnston A.D. Vapour-absent melting from 10 to 20 kbar of crustal rocks that contain multiple hydrous phases: Implications for anatexis in the deep to very deep continental crust and active continental margins // J. Petrol., 1996, V. 37, N 3, P. 661-691.
- 246 Singh J., Johannes W. Dehydration melting of tonalites. Part II. Composition of melts and solids // Contribs Mineral. and Petrol., 1996, V. 125, N 1, P. 26-44.
- 247 Kinzler R.J. Melting of mantle peridotite at pressures approaching the spinel to garnet transition: Application to mid-ocean ridge basalt petrogenesis // J. Geophys. Res. 1997, V. 102, N. B.1, P. 853-874
- petrogenesis // J. Geophys. Res., 1997, V. 102, N B1, P. 853-874.
 248 Edgar A.D., Charbonneau H.E., Mitchell R.H. Phase relations of an
 armalcolite-phlogopite lamproite from Smoky Butte, Montana:
 Applications to lamproite genesis // J. Petrol., 1992, V. 33, Part 3,
 P. 505-520.
- 249 Edgar A.D., Vukadinovic D. Potassium-rich clinopyroxene in the mantle:

- An experimental investigation of a K-rich lamproite up to 60 kbar //
- Geochim. Cosmochim. Acta, 1993, V. 57, N 23/24, P. 5063-5072. 250 Edgar A.D., Mitchell R.H. Ultra high pressure-temperature melting experiments on an SiO2-rich lamproite from Smoky Butte, Montana: Derivation of siliceous lamproite magmas from enriched sources deep in the continental mantle // J. Petrol., 1997, V. 38, N 4, P. 457-477.
- Patino Douce A.E. Experimental generation of hybrid silicic melts by reaction of high-Al basalt with metamorphic rocks // J. Geophys. Res., 1995, V. 100, N B8, P. 15623-15639.
 Rapp R.P., Watson E.B. Dehydration melting of metabasalt at 8-32 kbar:
- 252 Implications for continental growth and crust-mantle recycling // J. Petrol., 1995, V. 36, N 4, P. 891-931.
- Hirose K. Melting experiment on lherzolite KLB-1 under hydrous 253 conditions and generation of high-magnesian andesitic melts // Geology, 1997, V. 25, N 1, P. 42-44.
- 254 Walter M.J. Melting of garnet peridotite and the origin of komatiite and depleted lithosphere // J. Petrol., 1998, V. 39, N 1, P. 29-60.
- 255 Dunn T., Sen C. Mineral/matrix partition coefficients for orthopyroxene, plagioclase, and olivine in basaltic to andesitic systems: a combined analytical and experimental study // Geochim. Cosmochim. Acta, 1994, V.58, N2, P.717-733.
- Khodorevskaya L.I., Zharikov V.A. Experimental study of amphibolite melting concerning to genesis of tonalite-trondhemite magmatic series // In: Experimental and theoretical modelling of the processes of mineral formation (in Russian). Moscow, Nauka, 1998: p. 11-31.
- Moore G., Carmichael I.S.E. The hydrous phase equilibria (to 3 kbar) of an andesite and basaltic andesite from western Mexico: constraints on water content and conditions of phenocryst growth // Contribs Mineral. and Petrol., 1998, V. 130, N??, P. 304-319.
- Falloon T.J., Green D.H., Jacques A.L., Hawkins J.W. Refractory 258 magmas in Back-Arc basin settings - experimental constraints on the petrogenesis of a Lau Basin example // J. Petrol., 1999, V. 40, N 2, P. 255-277.
- Konzett J., Sweeney R.J., Thompson A.B., Ulmer P. Potassium amphibole 259 stability in the upper mantle: an experimental study in a peralkaline KNCMASH system to 8.5 GPa // J. Petrol., 1997, V. 38, N 5, P.537-568.
- 260 Konzett J. Phase relations and chemistry of Ti-rich K-richteritebearing mantle assemblages: an experimental study to 8.0 GPa in a Ti-KNCMASH system // Contribs Mineral. and Petrol., 1997, V. 128, N 4, P. 385-404.
- 261 Nimis P., Ulmer P. Clinopyroxene geobarometry of magmatic rocks Part 1: An expanded structural geobarometer for anhydrous and hydrous, basic and ultrabasic systems // Contribs Mineral. and Petrol., 1998, V. 133, N 1/2, P. 122-135.
- Falloon T.J., Green D.H., O'Neill H.St.C., Hibberson W.O. Experimental 262 tests of low degree peridotite partial melt compositions: implications for the nature of anhydrous near-solidus peridotite melts at 1 GPa // Earth Planet. Sci. Lett., 1997, V. 152, N ??, P. 149-162.

 Blatter D.L., Carmichael I.S.E. Plagioclase-free andesites from
- Zitacuaro (Michoacan), Mexico: petrology and experimental constraints // Contribs Mineral. and Petrol., 1998, V.132, N??, P.121-138.
- Kawamoto T., Holloway J.R. Melting temperature and partial melt chemistry of H2O-saturated mantle peridotite to 11 Gigapascals // Science, 1997, V.276, N??, P.240-243.
- Foden J.D., Green D.H. Possible role of amphibole in the origin of andesite: some experimental and natural evidence // Contribs Mineral. and Petrol., 1992, V.109, P.479-493.
- 266 Winther K.T. An experimentally based model for the origin of tonalitic and trondhjemitic melts // Chem. Geol., 1996, V.127, N.1, P.43-59.
- Draper D.S., Green T.H. P-T phase relations of silicic, alkaline, 267 aluminous mantle-xenolith glasses under anhydrous and C-O-H fluid-saturated conditions // J. Petrol., 1997, V.38, N9, P.1187-1224.
- Herzberg C., Zhang J.Z. Melting experiments on komatiite analog 268 compositions at 5 Gpa // Amer. Miner., 1997, V. 82, N3/4, P. 354-367.
- Springer W., Seck H.A. Partial fusion of basic granulites at 5 to

- 15 kbar: implications for the origin of TTG magmas // Contribs
- Mineral. and Petrol., 1997, V. 127, N 1/2, P. 30-45.

 Metrich N., Rutherford M.J. Low pressure crystallization paths of H2O-saturated basaltic-hawaiitic melts from Mt Etna: implications for open-system degassing of basaltic volcanoes // Geochim. Cosmochim. Acta, 1998, V. 62, N 7, P. 1195-1205.
- Pineau F., Shilobreeva S., Kadik A., Javoy M. Water solubility and D/H fractionation in the system basaltic andesite-H2O at $1250 \,\mathrm{mC}$ and
- between 0.5 and 3 kbars // Chem. Geol., 1998, V.147, N 1/2, P.173-184. Robinson J.A.C., Wood B.J., Blundy J.D. The beginning of melting of fertile and depleted peridotite at 1.5 GPa // Earth and Planet. Sci. 272 Lett., 1998, V. 155, N ??, P. 97-111.
- Kogiso T., Hirose K., Takahashi E. Melting experiments on homogeneous 273 mixtures of peridotite and basalt: application to the genesis of ocean island basalts // Earth and Planet. Sci. Lett., 1998, V. 162, N 1, P. 45-61.
- 274 Takahahshi E., Nakajima K., Wright T.L. Origin of the Columbia River basalts: melting model of a heterogeneous plume head // Earth and Planet. Sci. Lett., 1998, V. 162, N 1/4, P. 63-80. Kaszuba J.P., Wendlandt R.F. Effect of carbon dioxide on dehydration
- melting reactions and melt compositions in the lower crust and the origin of alkaline rocks // J. Petrol., 2000, V. 41, N 3, P. 363-386.
- Becker A., Holtz F., Johannes W. Liquidus temperatures and phase compositions in the system Qz-Ab-Or at 5 kbar and very low water activities // Contribs Mineral. and Petrol., 1998, V.130, P.213-224.
- 277 Wagner T.P., Grove T.L. Melt/harzburgite reaction in the petrogenesis of tholeiitic magma from Kilauea volcano, Hawaii // Contribs Mineral. and Petrol., 1998, V.131, N1, P.1-12.
- Putirka K. Garnet + liquid equilibrium // Contribs Mineral. and 278 Petrol., 1998, V.131, N??, P.273-288.
- Gaetani G.A., Grove T.L. The influence of water on melting of mantle peridotite // Contribs Mineral. and Petrol., 1998, V. 131, N ??, P.323-346.
- Thy P., Lesher C.E., Fram M.S. Low pressure experimental constraints on the evolution of basaltic lavas from site 917, southeast Greenland continental margin // Proc. Ocean Drilling Program Scientific Results, 1998, V. 152, P. 359-372.
- 281 Thy P., Lesher C.E., Mayfield J.D. Low-pressure melting studies of basalt and basaltic andesite from the southeast Greenland continental margin // Proc. Ocean Drilling Program Scientific Results, 1999, V. 163, P. 95-112.
- 282 McCoy T.J., Lofgren G.E. Crystallization of the Zagami shergottite: an experimental study // Earth and Planet. Sci. Lett., 1999, V. 173, P.397-411.
- 283 Sato H., Nakada S., Fujii T., Nakamura M., Suzuki-Kamata K. Groundmass pargasite in the 1991-1995 dacite of Unzen volcano: phase stability experiments and volcanological implications // J. Volcanol. Geotherm. Res., 1999, V. 89, N 1/4, P. 197-212.
- 284 Rapp R.P., Shimizu N., Norman M.D., Applegate G.S. Reaction between slab-derived melts and peridotite in the mantle wedge: experimental
- constraints at 3.8 GPa // Chem. Geol., 1999, V. 160, N ?, P.335-356. Martel C., Pichavant M., Holtz F., Scaillet B. Effects of fO2 and H20 285 on andesite phase relations between 2 and 4 kbar // J. Geophys. Res., 1999, V. 104B, N 12, P. 29453-29470.
- 286 Gardien V., Thompson A.B., Ulmer P. Melting of biotite+plagioclase+ quartz gneisses: the role of H2O in the stability of amphibole // J. Petrol., 2000, V. 41, N 5, P. 651-666.
- Cottrell E., Gardner J.E., Rutherford M.J. Petrologic and experimental 287 evidence for the movement and heating of the pre-eruptive Minoan rhyodacite (Santorini, Greece) // Contribs. Mineral. and Petrol., 1999 V. 135, N ??, P. 315-331.
- 288 Longhi J., Auwera J.V., Fram M.S., Duchesne J.-C. Some phase equilibrium constraints on the origin of proterozoic (massif)

INFOREX - 4.2 USER

- anorthosites and related rocks // J. Petrol., 1999, V. 40, N 2, P. 339-362.
- 289 Scaillet B., Evans B.W. The 15 June 1991 eruption of Mount Pinatubo. I. Phase equilibria and pre-eruption P-T-f02-fH2O conditions of the dacite magma // J. Petrol., 1999, V. 40, N 3, P. 381-411.
- dacite magma // J. Petrol., 1999, V. 40, N 3, P. 381-411.

 290 Dall'agnol R., Scaillet B., Pichavant M. An experimental study of a lower proterozoic A-type granite from the Eastern Amazonian craton, Brazil // J. Petrol., 1999, V. 40, N 11, P. 1673-1698.